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**Monitoring a patient's vital parameters and functions during treatment or outpatient surgery is key for all professionals.**

Today, non-invasive systems and devices for quick, safe monitoring are available, ensuring patient control and protection. The basic instrument normally used by all dentists and doctors is the pulse oximeter - an essential, accurate and cost-effective device.

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**The pulse oximeter, oximeter or saturation meter is a medical device to measure the amount of bound hemoglobin in the blood in a non-invasive way. The device detects the percentage of bound hemoglobin.**

Normally, hemoglobin binds with oxygen, so it is possible to measure the amount of oxygen in the blood. The device usually consists of a probe which takes the reading and a unit which calculates and displays the result. **The probe** of a standard pulse oximeter comprises a "clip", normally attached to the patient's finger. This clip is connected to the **calculation unit** which displays the reading via a **monitor**. The probe consists of diodes which emit bands of light in the red and infrared spectra and a photocell where light is conveyed after the bands have penetrated the skin and circulation of the patient. The oxygenated hemoglobin typically absorbs light at these wave lengths. With a known amount of initial and final light, the device can calculate oxygen saturation in the patient (abbreviated to **SpO<sub>2</sub>**). The probe is normally used in a surface circulation area, such as the finger or earlobe, because if circulation is too "deep", the bands of light cannot penetrate and cross the area and the reading cannot be taken. Oxygen saturation, pulse rate and intensity can be displayed at each heart beat.

The device is **non-invasive** (i.e. no actions or analyses are necessary that penetrate the patient's tissues) and it can detect **hypoxia** in relation to cyanosis, early on, allowing for a diagnosis of oxygen desaturation to be made prior to the onset of serious complications. The device **can be freely used**. It is normally employed by healthcare professionals (doctors and paramedics) as well as by other operators.

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**Vital Test is Tecno-Gaz new and last generation pulse oxymeter which is distinguishable for its solidity, precision and advanced features.**

**Vital Test** is small, easy transportable, equipped with shoulder bag and easy to be used both by surgeons and assistants. The item is **simple and professional**, it can be used in intensive care units, surgeries, dentist's surgeries, emergency rooms, since it guarantees immediate, accurate and non-invasive monitoring. **Vital Test** has pleasant and ergonomic shape, it is easy to grip and to be transported, equipped with wide display and with simple and inductive controls. **Vital Test** allows SpO<sub>2</sub> measurement, heart rate (pulse rate), waveform visualization and also last 10 minutes trend value of these two parameters. Supplied with rechargeable battery so that it can be used during transportation too.

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**LCD wide display** allows easy reading even at long distance, and it is also equipped with a control useful to activate or to remove retro-lighting in order to make it readable with all light conditions. Despite visualizing basic parameters it also visualizes battery level and alarms.

### Settings are the following:

- Wave visualization or value decimal reading
- Adults and children reading selection
- On/off sound and adjustable volume
- SpO<sub>2</sub> min/max alarms
- Pr min/max alarms
- Patient's identification
- Data reading
- Volume pulsation

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

The basic monitor indicates all parameters and functions: measured parameters included the shape of wave, date, hour, charge of the battery. The reading is simple, clear, direct and thanks to the backlight the reading can be facilitated in every situation.

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Pressing the upper needle placed on the board of the pulse oxymeter you can display the trend of SPO2

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
Pressing the lower needle placed on the board of the pulse oxymeter you can display the trend of the pulsations during the last 10 minutes.

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Set up of the graphics, definition of the patient adult/ child, alarm activation, alarm volume.

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
Set up alarms min/max indicating saturation and cardiac beat.

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Identification of the patient, date of data storage, date of transfer, volume of the pulsations.

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**Vital Test** can be supplied with proper support **code 11255**, planned to be fixed to dentist's unit, wall or trolley in order to allow a continuous vision even when hands are busy.

**An excellent solution both for dentists and outpatient surgeons.**

Interface program between computer and **Vital Test** can be supplied on demand in order to record and print all collected parameters.

**Vital Test** is compatible with **Nellcor** and **BCI D89** sensors.

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All rights reserved. Variations can be done without notice. Tecno-Gaz S.p.A. is to be considered not responsible for damages caused by the lack or the wrongness of the information here mentioned.