europa B pro
europa B pro 20
europa B pro 24

<i>Istruzioni per l'uso</i>	<i>ITALIANO</i>
<i>Instructions for use</i>	<i>ENGLISH</i>
<i>Mode d'emploi</i>	<i>FRANCAIS</i>
<i>Instrucciones de uso</i>	<i>ESPAÑOL</i>
<i>Bedienungsanleitung</i>	<i>DEUTSCH</i>

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INDEX

- 1. INTRODUCTION**
- 2. USE AND DESTINATION OF USE OF THE AUTOCLAVE**
- 3. SAFETY**
 - 3.1 SAFTY MARKING
 - 3.2 SAFETY DEVICES
 - 3.3 SAFETY NOTES
 - 3.4 DISPOSAL
- 4. TECHNICAL DATA**
- 5. OUTPUTS AND INDICATORS LIST**
- 6. UNPACKAGING**
- 7. ACCESSORIES**
- 8. INSTALLATION**
- 9. PROGRAMMING DISPLAY**
- 10. USE INSTRUCTION**
 - 10.1 TURNING ON THE AUTOCLAVE AND BAROMETRIC ALIGNMENT
 - 10.2 CLEAN WATER TANK MANUAL LOAD
 - 10.3 CHARACTERISTICS OF THE WATER TO BE USED
 - 10.4 LOADING OF THE MATERIALS INTO THE AUTOCLAVE
 - 10.5 STARTING THE STERILIZATION CYCLE
 - 10.6 END OF THE STERILISATION CYCLE
 - 10.7 UNLOADING OF THE STERILIZED MATERIALS
 - 10.8 DISCHARGING USED WATER
 - 10.9 INTERRUPTION OF A STERILIZATION CYCLE
- 11. STERILIZATION CYCLES**
 - 11.1 CYCLES DESCRIPTION
 - 11.2 CYCLE DIAGRAM
 - 11.3 READING OF THE CYCLE REPORT
- 12. SOFTWARE VISUALIZATION CYCLES**
- 13. AUTOCLAVE CHECK TEST**
 - 13.1 CHEMICAL INTEGRATORS
 - 13.2 BIOLOGICAL INTEGRATORS
 - 13.3 BOWIE & DICK TEST
 - 13.4 HELIX TEST
 - 13.5 VACUUM TEST
- 14. OSMOSI SYSTEM INSTALLATION**
- 15. MAINTENANCE**
- 16. ALARM MESSAGES**

17. SOLUTION TO OPERATING PROBLEMS

- 17.1 THE AUTOCLAVE DOES NOT DRY CORRECTLY
- 17.2 THE AUTOCLAVE CHAMBER HAS TURNED WHITE
- 17.3 THE AUTOCLAVE CHAMBER HAS GREEN STAINS
- 17.4 THE STERILIZATION CYCLE INTERRUPTS
- 17.5 THE AUTOCLAVE DOES NOT RECEIVE CONTROLS
- 17.6 STAINS ON INSTRUMENTS

18. PROCEDURE FOR SERVICE AND ASSISTANCE

A. CONSUMABLE SPARE PARTS SUMMARY

Dear Client,

Thank you for having chosen our autoclave, we know how to exchange your fidelity, with maximum attention and service definitely corresponding to your expectations.

Before using this autoclave, we invite you to read with maximum attention the user's manual and then keep it in a place accessible to all operators in charge of STERILISATION.

Sterilise means adopting a specific working methodology and adhering to precise operational protocols:

DISINFECTION obligatory phase, to ensure operator safety, to be done with immersion in chemical liquids or thermo-disinfection

CLEANSING the most important phase that ensures the removal of all types chemical and organic residues. The most suitable instruments are ultrasound baths;

DRYING essential phase, which prevents corrosion of the instruments and interference of the STERILISATION cycle;

ENVELOPING essential phase for sterility maintenance over time;

STERILISATION final step steam STERILISATION.

The autoclave is the key point of this methodology.

We remind you that failure to carry out all the phases of the STERILISATION process may invalidate the final result.

For installation, maintenance and assistance ask exclusively for a technician authorized. We invite you to use and ask for exclusively original spare parts.

2 USE AND DESTINATION OF USE OF THE AUTOCLAVE

The autoclave is able to sterilize the three types of load provided for by the standard EN13060:2009, especially:

	AUTOCLAVE 18L	AUTOCLAVE 20L	AUTOCLAVE 24L
<u>METAL OR SOLID MATERIALS</u> Instruments with no cavities and no obstacles to the penetration of steam	<i>max kg. 5</i>	<i>max kg. 6</i>	<i>max kg. 7</i>
<u>POROUS OBJECTS</u> Simple or composite materials that can absorb fluids (fabrics, gowns, surgical gauzes, dressings, etc)	<i>max kg. 1,5</i>	<i>max kg. 1,75</i>	<i>max kg. 2</i>
<u>HOLLOW OBJECTS</u> Materials or devices with cavities, obstructions, etc. These are subdivided into two types, classified according to the length and diameter of the cavity. Approximately: TYPE B: cannulas, tubes or devices with large passages. TYPE A: turbines, hand pieces and devices with blind or small holes	<i>max kg. 5</i>	<i>max kg. 6</i>	<i>max kg. 7</i>

**Only for European Countries*


The autoclave shall only have to be used for the sterilization of tools and materials being compatible with the steam sterilization system. Always make sure that the loads that need to undergo sterilization can stand the temperatures scheduled for the selected cycle.

3

SAFETY

3.1 Safety marking

 <p>HAZARDOUS VOLTAGE</p>	 <p>ATTENZIONE ATTENTION ATTENTION ACHTUNG</p> <p>ALTA TEMPERATURA HIGH TEMPERATURES TEMPERATURES ELEVÉES HOHE TEMPERATUR</p> <p>HIGH TEMPERATURE</p>
 <p>ATTENZIONE TOGLIERE TENSIONE PRIMA DI RIMUOVERE IL COPERCHIO</p>  <p>WARNING DISCONNECT THE MAINS SUPPLY BEFORE REMOVING THIS COVER</p> <p>DISCONNECT THE POWER BEFORE REMOVING THE LID</p>	 <p>EARTH CONNECTION</p>

3.2 Safety devices

The following safety devices are installed:

-) Safety valve set at 2.4 bar 0/+10%
-) Electromagnetic lock to prevent the door from opening while the cycle is running
-) Resistance over temperature thermostats

3.3 Safety notes

- The manufacturer is liable for the marketed product in accordance with current regulations. The manufacturer's liability will expire when operations are carried out on the device, or a part of it, by unskilled personnel or using non-original spare parts.
- There should be no potential risk of explosion and/or fire in the room where the autoclave is installed.
- The autoclave should be installed in a special well-ventilated room.

3.4 Disposal

This product is subject to Directive 2002/96/EC of the European Parliament and the Council of the European Union on Waste of Electrical and Electronic Equipment (WEEE) and, in jurisdictions adopting that Directive, is marked as being put on the market after August 13, 2005, and should not be disposed of as unsorted public waste. Please utilise your local WEEE collection facilities in the disposition and otherwise observe all applicable requirements.

4
TECHNICAL DATA

		18L	20L	24L
MECHANICAL DATA	Working temperature	+5°C ÷ +40°C		
	Maximum altitude	2.000 m		
	Max relative humidity at 30°C	80%		
	Max relative humidity at 40°C	50%		
	Dimensions of space occupied (L x H x P) mm	528 x 400 x 638	504 x 400 x 730	
	Space occupied with open door	300 mm		
	Weight (empty tanks)	60kg.	62kg.	64kg.
	Weight (full tanks + room full)	70 kg.	72kg.	74kg.
	Weight of area of support	2058 N/m ²		
	Potential sound level	< 70 db A		
ELECTRICAL DATA	Power voltage	230 V a.c. +/-10 % single phase		
	MAX power	1,5 kW	2,2 kW	
	Frequency	50 / 60 Hz		
	Power cord	2 + 1 x 1mm ²		
	Fuses	5x20 10A		
	Heat transmitted	3.6 E ⁶ J / hour		
CHAMBER	MAX working pressure	2.4 bar (relative)		
	MAX empty	- 0.9 bar (relative)		
	MAX temperature	138 °C		
	Material	Inox AISI 304		
	Size (mm)	Ø 245 x 318	Ø 245 x 430	Ø 245 x 500
CLEAN WATER TANK	Volume	4,5 l		
	Usable cycles	4	3	2
	Material	polyethylene		
USED WATER TANK	Volume	4,5 l		
	Usable cycles	4	3	2
	Material	polyethylene		
	MAX temperature used water	50°C		
BACTERIO LOGICAL FILTER	Diameter	56 mm		
	Filtering capacity	0.3 µm		

05

OUTPUTS AND INDICATORS LIST

OUTPUTS AND INDICATORS LIST (Fig.A)	00	Spacer
	01	Tap for emptying used water tank
	02	Tap for emptying clean water tank / Tap filling up trough osmosis
	03	Exceeding full for used water tank
	04	Pipe fitting for emptying used water tank
	05	Feeding socie with fuses
	06	Socket for connection of osmosis system
	07	Bacterial filter
	08	Pipe fitting for manual willing up of clean water
	09	Pipe fitting for discharging used water
	10	Display
	11	SD CARD Slot
	12	RS232 serial port
	13	Main switch
	14	121°C Cycle (valido per materiale imbustato e non i mbustato)
	15	134°C Cycle (valido per materiale imbustato e non i mbustato)
	16	Helix / Bowie&Dick Test
	17	<i>MAX level clean water tank</i>
	18	<i>MIN level clean water tank</i>
	19	<i>MAX level used water tank.</i>
	20	134°C Prion Cycle
	21	134°C Flash Cycle
	22	Vacuum Test
	B-M1	<i>Multifunction push button 1</i>
	B-M2	<i>Multifunction push button 2</i>
	B-M3	<i>Multifunction push button 3</i>
	PUMP WATER	<i>Push button for filling up of water</i>
	SELECT CYCLE	<i>Push button for selection of cycle</i>

06

UNPACKAGING

The autoclave is shipped in a suitable package to be transported and moved easily and to protect its contents.

The packaging must not be subject to impact, must be handled with care and avoid dropping it or roll it.

In case autonomous handling means are not available handle the packaged autoclave always in two persons.

Autoclave is placed on wooden pallet and packed with corrugated cardboard internal and external application.

For remove the autoclave to its pack please remove first the corrugated cardboard. For lift the autoclave please use the belts.

The autoclave shall have to be handled by at least two people and by using belts only.

Do not lift the autoclave from the inferior part of the door or command panel, this incorrect operation can create problems of a mechanical nature.

ATTENTION: Always conserve original packaging.

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
- **USER'S MANUAL:** to be read with attention and kept in a place available to all operators assigned to sterilization.
- **CERTIFICATION:** which must be conserved.
- **INSTALLATION REPORT - TESTING AND GUARANTEE CONDITIONS:** must be completed upon installation of the machine following the instructions indicated on the form.
- **QUICK USE GUIDE:** should be kept in the vicinity of the machine.


07

ACCESSORIES


TRAY HOLDER


	18L	20L	24L
Material	Aluminium anodized		
Size (L x H x P) (mm)	192 x 165 x 280	192 x 165 x 370	192 x 200 x 470
Picture			
	Fig.1		
Envelope standard	1		
Code	DPLA048	SNDA048	DXLA348



DUAL TRAY HOLDER (optional)


	18L	20L	24L
Material	Aluminium anodized		
Size (L x H x P) (mm)	193 x 200 x 280 / (200 x 193 x 280)	-	-
Picture			
	Fig.2		
Envelope standard	1	-	-
Code	DPLA048 + DPLA058	-	-


TRAY

	18L	20L	24L
Material	Aluminium anodized		
Size (L x H x P) (mm)	184 x 286 x 17	370 x 185 x 17	184 x 17 x 286 184 x 17 x 140
Picture			
	Fig.3		
Envelope standard	4		4 + 4
Code	DANA049	DHLA049	DANA049 + DXLA348

TRAY (optional)			
	18L	20L	24L
Material	Aluminium anodized		
Size (L x H x P) (mm)	-	-	185 x 17 x 460
Picture	-	-	
	Fig.4		
Envelope standard	-	-	4
Code	-	-	DXLA349


TRAY EXTRACTION AND DOOR ADJUSTMENT WRENCH		
Use for extract the trays and for door adjustment (<i>par. 15</i>)		
Picture	 Fig.5	 Fig.6
Envelope standard	1	
Code	DANA008	

CHAMBER AND DOOR GASKET CLEARING SPONGE	
Use to clean sterilization chamber and door gasket (<i>par. 15</i>)	
Picture	 Fig.7
Envelope standar	1
Code	CPMG004

WATER FILLING PIPE WITH FILTER AND PIPE FITTING	
Use for manual water loading (<i>par. 10.2</i>)	
Picture	 Fig.8
Envelope standard	1
Code	DANA099 + DXBA711 + CPRG117


WATER DISCHARGE PIPE

Use to drain the water used from the faucet on the front of the autoclave (*Fig.A-pos.09*) - (*par. 10.8*)

Picture	 Fig.9
Envelope standard	1
Code	DANA130


BLACK PLASTIC SPACER PIN

Put the spacer in the autoclave's back panel (*FIG.1-pos.00*) It's necessary for guarantee a good ventilation if you place the autoclave near a wall.

Picture	 Fig.10
Envelope standard	1
Code	CPAP014

FITTING FOR REAR TAPS DISCHARGE


Tighten the loading faucet (*Fig. A-pos.02*) to use the osmosi system, or to empty the loading tank, and tighten the drain faucet (*Fig. A-pos.04*) to empty the drain tank.

Picture	 Fig.11
Envelope standard	2
Code	CPRG096

PIPE FOR DISCHARGE UTILITIES


1 - Clean water overflow hose: connect one end of the hose to the clean water overflow hose fitting (*Fig. A-pos.01*), and the other in a water recovery container.

2 - Used water overflow hose: connect one end of the hose to the fitting (*Fig. A-pos.03*), and the other to a used water recovery container.

Picture	 Fig.12
Envelope standard	2
Code	SXBA799

POWER SUPPLY CORD

Take the power supply cord provided (pos.B) and insert the female plug (FIG.A-pos.05) in the socket of the back panel of the autoclave then insert the male plug (pos.A) in the electric plug of the system.

Picture	 <p style="text-align: center;">Fig.13</p>
Envelope standard	1
Codee	CECG006

MEMORY CARD


To be used to memorize the autoclave cycles (for the best compatibility with the device, it is advisable to always use the original memory card).

WARNING: The card contains the software reading log cycle - Carry out the rescue and installed on the PC before the commissioning of the autoclave (see par.12)

Picture	 <p style="text-align: center;">Fig.14</p>
Envelope standard	1
Code	CEGS001

08

INSTALLATION

	18L	20L	24L
A	640 mm	755 mm	
B	504 mm		
C	400 mm		
D Max. hatch opening	382,5 mm		
E	426 mm		
F	455 mm		
G	402 mm	435 mm	
H	min. 50 mm		
I	72 mm	79 mm	
L	164 mm	245 mm	
M	525 mm	564 mm	

- 1 – Install the autoclave in environments suitable to sterilization.
- 2 – The premises have to be suitably lit and ventilated, in compliance with regulations in force.
- 3 – Install the autoclave away from heat sources and water splashes.
- 4 – Position the autoclave on a surface, suitable to support the weight (80 Kgs.) and with adequate dimensions.
- 5 – Place the autoclave at a height allowing easy intervention by the operator for inspection and clearing of the whole sterilization chamber.
- 6 - Open the door of the autoclave and remove all packages containing the single accessories from inside the sterilization chamber.
- 7 – Leave inside the sterilization chamber only the tray carrier with the trays. All other accessories should be positioned in a separate space available to operators.
- 8 – Don't put anything on the autoclave.
- 9 – Do not lean on the door.
- 10 – Leave a space of at least 5 cm in the rear of the autoclave, using the spacer (*FIG.A-pos.00 / FIG.10*) and in the sides of the unit to ensure the ventilation required.
- 11 - Make connections of pipes in the back (Chapter 7)
- 12 – Always make sure the electrical system to which the autoclave is to be connected is in conformity with the legislation in force and sized to suit the specifications of the said device.
- 13 – Connect the power supply cord to the socket on the rear panel of the autoclave (*FIG.A-pos.5*)
- 14 - Connect the electrical plug to the system ensuring that it is adequate to the supply of the machine.

NOTE:	Do not use extensions, reducers or adapters for connection as this could cause micro interruptions with consequent generation of alarm signals.
--------------	---

- 15 - Turn the autoclave on by pressing the mains switch (*Fig. A-pos.13*) and open the hatch of the autoclave itself. Wait a few seconds, there will be two acoustic signals to inform you of the acquisition of the parameters for the automatic barometric alignment, at the same time the text DOOR OPEN will appear on the display.

NOTE:	Never select a command before hearing the two sound signals, the autoclave will not accept the chosen programming.
--------------	--

9

PROGRAMMING DISPLAY

From the first page press **Setup** to enter the autoclave setup menu

LANGUAGE

Press the central button to change the language of menus and vocal messages

Press the arrow to shift to the next item

DATE AND TIME

Press the central button to enter the date/time setup.

00:00

00/00/1900

When the box flashes press the arrows to select the desired value.

Press once more the central button to move among boxes and the arrows to select a value.

Go on until the last value is selected.

Press the central button for the last time to display the final selection page.

00:00
00/00/1900

Press to confirm selected values and go back to the setup menu

Press to restart the procedure


Press to cancel selected values and go back to the setup menu

Press the arrow to move to the next item

EXTERNAL PRINTER


After the printer installation (*optional accessory*), is possible to print the traceability adhesive labels to be stuck on the wrap before starting sterilization is entered.

When the box flashes press the arrows to select the number of labels to be printed.

Press the central button to move to the next selection

When the box flashes press the arrows to select the number of months before expiry.

Press the central button for the last time to display the final selection page.

Press the arrow to move to the next item

OSMOSIS CONTROL

- OSMOSIS SYSTEM : Activate / Disconnect with the button  the supply system with osmosis device (*optional*). When the system is activated the charge pump is disconnected.
- CYCLES FROM THE CHANGE OF THE FILTER: Cycles from the last change of filters are visualized.
- RESET THE COUNTER: It allows you to reset the counter when you replace the filters.

SERVICE MODE

Access to the service menu (after password entry). This mode is reserved exclusively for the setup by an **authorised personal**. The manufacturer shall not be held responsible for any tampering with or injury to unauthorised staff.

10**USE INSTRUCTION**

After installing the autoclave proceed with preparation and use.

10.1 Turn on the autoclave and barometric alignment

Press the main switch (*Fig.A-pos.13*). After the display of the logo, the autoclave checks the memory and the connections. Once checks are over the autoclave shifts to the operating setup.

Open the door and wait for some seconds until an acoustic signal informs that the automatic barometric alignment values have been acquired; at the same time the display shows the message DOOR OPEN.

THE AUTOCLAVE IS READY FOR USE**10.2 Clean water tank filling**

Connect the hose supplied (*Fig. 8*) to the front fitting of the autoclave (*Fig. A-pos.08*).

Put the other end of the hose with the filter inside the demineralised or distilled water container.

At this point, press the **PUMP WATER** button to operate the water loading pump and keep it pressed until the countdown appears.

The pump loads the clean water tank inside the autoclave. If the maximum level is not achieved within 180 seconds the pump stops automatically and it will be necessary to press button **B-PUMP**.

The pump stops automatically when the maximum level is achieved.

10.3 | Characteristics of the water to be used


TABLE SHOWING THE QUALITY LEVELS LAID BY THE UNI EN 13060: 2009

CEN STANDARD UNI EN 13060 : 2009			
Evaporation residue	≤	10	mg/l
Silicon oxide	≤	1	mg/l
Iron	≤	0.2	mg/l
Cadmium	≤	0.005	mg/l
Lead	≤	0.05	mg/l
Remains of heavy metals apart from iron, cadmium, lead	≤	0.1	mg/l
Chloride (Cl')	≤	2	mg/l
Phosphate (P20s)	≤	0.5	mg/l
Conductivity (at 20°C)	≤	15	μs/cm
Ph value (acidity level)	5 ÷ 7,5		
Appearance	Transparent, clear, without deposits		
Hardness (and alkaline earth ions)	≤	0.02	mmol/l

10.4 | Loading of the materials in the autoclave

Arranging the materials to be sterilized on the provided trays, as follows:

- Do not superpose the materials
- Arrange the wrapped materials with the paper side facing upwards
- Never bring the materials into contact with the sterilization chamber or the autoclave door
- Put scissors and dental forceps with blades open

After loading the materials, close the door of the autoclave. The display will show the icon and the message DOOR CLOSED.

10.5 Starting the sterilization cycle

After the steps listed above, choose the most suitable STERILISATION program for the prepared load, by pressing the **SELECT CYCLE** button.

After choosing the program, start the cycle by pressing the **Start** button. The door will lock automatically, and the cycle will begin.

During the cycle the display shows all parameters and information related to the cycle in progress. In this setup the display shows: the type of cycle, the cycle status, the time remaining before the end of cycle (for the Vacuum Test it identifies the whole cycle, while for all the other cycles it identifies the sterilization phase plus that of drying), the number of cycles done by the machine and the button **Info** allowing access to the list of working parameters.

10.6 End of cycle

An acoustic signal will inform the operators about the finished STERILISATION cycle and the display will show the END CYCLE icon and message.

Unlock the door by pressing the **Unlock** button displayed on the display by one of the three multi-function buttons. Should there be pressure present inside the chamber, the button will not trigger the unlocking. Wait for the complete depressurization of the chamber and repeat the operation. With the door unlocked, pull the door handle and open.

10.7 Unloading the sterilized materials

Wear personal protective equipment in compliance with regulations on safety and hygiene at work. Extract the trays by using the spanner provided (FIG.5/6), let rest the instruments and store them in environments which are not exposed to contamination.

10.8 Discharge used water

When the used water level LED (Fig. A-pos.19) lights up, the used water collection tank must be emptied.

If nothing is done, the operation of the autoclave is blocked.

Take the hose supplied (Fig.9), and insert it into the used water drain fitting located on the front of the autoclave (Fig.A-pos.09). Put the other end of the hose in a container and unscrew the ring nut, turning it anticlockwise, the water will fall into the container, emptying the tank.

IMPORTANT:

A - The hose located in the receiving container, must never touch or be immersed in the discharged water, otherwise sucking up will occur.



B – Always wait until the discharge water has been completely drained. The used water maximum level LED turns off when some water still remains into the tank, consequently do not consider it as a reference for this operation.

When drained, screw in the ring nut and remove the hose.

10.9 Interruption of a sterilization cycle

A sterilization cycle can be voluntarily interrupted by pressing the button **Start** for at least 2 seconds.

The autoclave will send out a beep and the display will show the alarm message A001.

For reset the alarm, press at the same time, the multifunction push buttons on the Reset bar up to the cancelling of the same.

11

STERILIZATION CYCLES

11.1 Descriptions cycles

The autoclave has three series of cycles:

A – OPERATION CYCLES

All operation cycles have a system of fractioned vacuum, which therefore permits sterilization of materials that are hollow, porous, solid, both free and packaged. Temperatures can be selected from 121°C – 134°C. Normally the cycles of 121°C are used for thermoplastics or sensitive materials, while the 134°C cycles are used for all other materials. In all cases always follow the indications given by the manufacturer of the instruments or devices to be sterilized. There are also two special cycles:

- Prion cycle, studied for the Crutzfeldt-Jacobs disease (mad cow disease)
- Flash cycle created for rapid sterilization of unpacked instruments and devices.

B – NIGHT CYCLES

The autoclave is equipped with a special economiser device.

This device allows you to perform the STERILISATION cycles in the absence of the operator.

At the end of the cycle, if the door is not opened, the autoclave is stabilised, it then turns off automatically, only the mains switch will remain on (Fig. 1-pos.13)

When the operators arrive, simply press any button to turn the autoclave back on and read the result of the cycle on the display. In addition, the printer will normally have prepared the written report on the cycle.

C – TEST CYCLES

The available test cycles are the following:

- Helix test
- Bowie&Dick test
- Vacuum test

11.2 Cycle diagram


TABLE OF TESTS TYPE AS RULE EN13060:2009	
<i>Type test</i>	<i>Operative cycle</i>
<i>Dynamic sterilizer chambre pressure</i>	X
<i>Air Leakage</i>	X
<i>Empty chamber</i>	X
<i>Solid load</i>	X
<i>Small porous items</i>	X
<i>Small porous loads</i>	X
<i>Full porous load</i>	X
<i>Hollow load B</i>	X
<i>Hollow load A</i>	X
<i>Multiple wrapping</i>	X
<i>Dryness, solid load</i>	X
<i>Dryness, porous load</i>	X



12

SOFTWARE VISUALIZATION CYCLES


Insert the memory card on the front of autoclave into a memory card reader connected to your PC. Open “My computer” (or similar) and double-click on external disk.

Open *index.htm*

Make sure that the ActiveX control is active to provide access to resources.

Press Apri file log and browse through the folders until you reach the one called LOG containing the save file of cycles made by the machine.

After selecting the desired cycle, in the mask are loaded the contents of all files present on the memory card containing the generic data cycles (date, cycle number, cycle type, outcome).

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Fine

Risorse del computer 100%

Selecting a cycle from a list, you will get a second mask containing all the details of the cycle.

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
Pres Anteprima di stampa, will be displayed a page in text format showing autoclave data and autoclave data cycle in a structure suitable to be printed.

Riga	Data	Time left	T1	T2	PI	V rete	Cycle step
1	17/12/2009 14:37:48	0.00	73.74	75.34	-0.011	234.7	START-UP
2	17/12/2009 14:37:56	0.00	73.97	75.47	-0.010	232.8	START-UP
3	17/12/2009 14:37:56	960.00	74.47	75.86	0.000	234.6	VACUUM
4	17/12/2009 14:38:05	960.00	74.55	76.15	-0.151	226.5	VACUUM
5	17/12/2009 14:38:15	960.00	74.76	76.85	-0.210	227.4	VACUUM
6	17/12/2009 14:38:26	960.00	75.70	77.56	-0.342	223.5	VACUUM
7	17/12/2009 14:38:36	960.00	76.09	77.74	-0.400	225.8	VACUUM
8	17/12/2009 14:38:46	960.00	75.57	77.09	-0.483	225.1	VACUUM
9	17/12/2009 14:38:56	960.00	74.37	75.94	-0.542	226.4	VACUUM
10	17/12/2009 14:39:05	960.00	72.87	74.00	-0.608	225.8	VACUUM
11	17/12/2009 14:39:15	960.00	72.25	73.06	-0.638	226.7	VACUUM
12	17/12/2009 14:39:26	960.00	69.47	70.00	-0.690	225.9	VACUUM
13	17/12/2009 14:39:36	960.00	66.98	67.40	-0.718	226.7	VACUUM
14	17/12/2009 14:39:46	960.00	65.91	66.30	-0.728	226.7	VACUUM
15	17/12/2009 14:39:56	960.00	64.92	65.31	-0.739	224.0	VACUUM
16	17/12/2009 14:40:06	960.00	64.24	64.60	-0.744	227.8	VACUUM
17	17/12/2009 14:40:17	960.00	63.95	64.14	-0.760	224.9	VACUUM
18	17/12/2009 14:40:27	960.00	62.86	62.44	-0.777	225.4	VACUUM
19	17/12/2009 14:40:37	960.00	62.63	61.97	-0.776	226.8	VACUUM
20	17/12/2009 14:40:47	960.00	62.29	60.93	-0.793	226.7	VACUUM

Selecting: File->Stampa you can select the desired printer and start the process of printing of the document.

Before removing the memory card you should disconnect the device, to avoid errors or file system corruption on the card itself.

We recommed to make a copy of the entire contents of the memory card on your computer at least one a week and arrange for appropriate storage backups.

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A	Indication of date and time when the machine makes a reading of parameters
B	Time remaining at the end of cycle (in seconds – is the sum of the time of sterilization and the drying)
C	Reading the two temperature probes
D	Reading of chamber pressure
E	Reading of voltage
F	Cycle information, serial machine and firmware version
G	Operational phase of the cycle

13

AUTOClave VALIDATION TESTS

13.1 Chemical integrators

These tests (spare code 200/S and 215-S) are based on the properties of indicator dyes that are able to modify their own aspect if they are adequately exposed to heat and pressure according to given times, since they interact with all the cycle parameters (pressure, temperature, time). If the test is successfully completed, the correct stability is validated.

USE	<p>The chemical integrators (spare code 200/S o 215-S) are to be put into the sterilization chamber before starting the cycle, even with the load already in. The integrators can either be unwrapped and positioned into the trays or they can be wrapped into the sterilization bags together with the materials to be sterilised. In any case it is advisable to wrap the integrators where there are many operators. As a matter of fact by doing so the sterilization of every single tool will be certified.</p>	
RESPONSE		<p>If the colour change remains within the white area, the test <u>failed</u></p>
		<p>If the colour change enters the SAFE window, the test <u>succeeded</u></p>
PERIODICITY	<p>There is no preset testing interval, apart from areas to which specific regional regulations apply. In any case all the chemical integrators are the cheapest and quickest available tests, hence in order to constantly monitor the autoclave performance, it is advisable to use them on all of the cycles or at least once a day.</p>	

13.2 Biological indicators

The test (spare code 262-S) is used to prove the autoclave capability as to the complete destruction of all microorganisms. It is represented by a standardized spore preparation, having biological characteristics, high resistance to heat, which is to be deemed an absolutely safe control tool (ATCC 7953).

It consists of spore ampoules. It is non pathogenic, non toxic and non pyrogenic.

USE	<p>The ampoules shall have to be put into the sterilization chamber even when loaded.</p> <p>Carry out the cycle and then take out the ampoule by taking care since it is hot and pressurized.</p> <p>Leave it cool down for about 10 minutes and then activate it by using the supplied crusher by taking care so as to keep the ampoule always upright.</p> <p>The chemical indicator placed on the label of the ampoule will have changed its colour from blue into black. Subsequently incubate the ampoule into a biological incubator at 57°C for 48 hours, together with an unprocessed test ampoule, that has anyway be activated through the relevant crusher.</p>
------------	---

RESPONSE	<p>After 48 hours take the processed ampoule out of the incubator and assess the test result. If the ampoule has changed its colour to yellow this means that the autoclave has failed the test and there is bacterial growth. If the ampoule remains purple, this means that there is no bacterial growth and hence the test succeeded. The test ampoule of course will always change to yellow since it has not been processed and will only be used as term of comparison.</p> <p>At the end of the test dispose of the ampoule in the solid urban waste. It is anyway advisable to perform another sterilization cycle at 121° C on the ampoule.</p>
PERIODICITY	<p>There is no preset testing interval, apart from areas to which specific regional regulations apply. In any case the biological tests are the more concrete ones. Hence in order to constantly monitor the autoclave performance, it is advisable to perform the test at least once every 90 days.</p>

13.3 Bowie & Dick test



It is a physical test (spare code 268/S) enabling to verify steam penetration capability into porous loads. For the test a “test pack” is used which is standardized and complying with the technical standards being currently in force.

USE	<p>Perform the test with empty chamber. The Bowie-Dick test pack shall have to be positioned on the autoclave central tray. Perform the cycle, indicated in the display of the autoclave, and at the end remove the pack, open it and check the sheet colour change through the chemical indicator placed inside.</p>	
RESPONSE		<p>The assessment is simple and rapid. If the colour change turns out to be uniform (as shown in the figure) the test perfectly succeeded. On the contrary the test failed, hence the autoclave cannot properly sterilize porous bodies.</p>
PERIODICITY	<p>There is no preset testing interval, apart from areas to which specific regional regulations apply. Hence in order to constantly monitor the autoclave performance, it is advisable to perform the test at least once every 30 days.</p>	

13.4 Helix test

It is a physical test (spare code 267-S) enabling to verify steam penetration capability into hollow loads. For the test a standardized system is used which complies with the technical standards in force.

ENGLISH

USE	Put the test device into the empty chamber. Fit the suitable indicator strip into the test device end cap, then position the test device on the central tray of the autoclave. Perform the suitable cycle, indicated in the autoclave display, then remove the test device, open the cap and check the indicator strip colour change.	
RESPONSE	 ART.NO. 134 °C - 3,5 MIN. / 121 °C - 15 MIN. 	BEFORE THE TEST
	 ART.NO. 134 °C - 3,5 MIN. / 121 °C - 15 MIN. 	TEST FAILED
	 ART.NO. 134 °C - 3,5 MIN. / 121 °C - 15 MIN. 	TEST SUCCEEDED
PERIODICITY	There is no preset testing interval, apart from areas to which specific regional regulations apply. Hence in order to constantly monitor the autoclave performance, it is advisable to perform the test at least once every 30 days.	

13.5 Vacuum test

It is a chamber leak test or a vacuum leak test. This test aims at verifying that during the cycle no air infiltration takes place through the chamber seals (gaskets, valves, etc...)

USE	Perform the test with empty chamber. Select the suitable cycle indicated on the autoclave display. The cycle is enabled. The autoclave automatically performs the cycle by complying with precise technical procedures. Then the final result is printed out through the printer connection.
RESPONSE	The test result is immediate and is printed by the printer which prints out all the values relating to the performed cycle as well as the final assessment.
PERIODICITY	There is no preset testing interval, apart from areas to which specific regional regulations apply. Hence in order to constantly monitor the autoclave performance, it is advisable to perform the test at least once per week.

**WE RECOMMEND TO USE ONLY TEST PROPOSED
BY TECNO - GAZ S.p.A**

14

INSTALLATION OSMOSIS SYSTEM

Setting osmosis system

The autoclave is also designed for loading demineralised water through an external osmosi demineralisation system (*optional accessory*).

The operator before installing the system, must schedule the autoclave following the instructions below:

Turn the autoclave on by pressing the mains switch (*Fig. A-pos.13*).

When the autoclave is brought up on the operating screen, press **Setup** and enter the osmosi management menu.

Activate by setting the OSMOSI SYSTEM item to **ON**.

Exit the menu and return to the operating screen.

NOTE

Connected by a demineraliser, if the maximum water level has not been reached, the operation of the autoclave will be blocked.

WARNING:

The number on the osmosi management screen indicates how many STERILISATION cycles have been performed since changing filters.

When you reach the maximum number of cycles, a message on the display will inform the user. It's a good idea to remember to reset the cycles counter when changing the filter on the system.

Connection of the demineralizer

Here as follows the detail of the specific connection to the autoclave of the water supply hose and electrical plug connection:

Turn off the autoclave, if ON (*Fig.A-pos.13*)

- Close the cock upstream the demineralizer
- Install the demineralizer as indicated in its manual;
- Wrap the male thread of the hose coupler with teflon or another component ensuring water tightness;
- Screw the hose coupler (*Fig.11*) on the female thread of the clean water discharge (*FIG.A-pos.2*)
- Insert the hose from the demineralizer into the hose coupler and screw it to the autoclave;
- Insert the feeding coupler of the demineralizer into the port (*FIG.1-pos.6*) on the rear panel of the autoclave;
- Open the valve upstream the demineralizer;
- Be sure there is no water leakage ;

ENGLISH

- Switch on the autoclave;
- Carry out one or several sterilization cycles to check for the correct connection and mainly to be sure of the absence of any leak.

At the end of every day always close the cock upstream the demineralizer

Connect the demineralizers exclusively to suitable autoclaves

NOTE

For the connection of the demineralizers to the autoclaves, see also the demineralizer manual.

15

MAINTENANCE

Correct maintenance of the autoclave assures correct functioning of it and a secure saving in terms of time and costs for assistance and maintenance. The following operations are compulsory feasible by operators.

Cleaning of chamber

Every 20 cycles or one a week


Clean periodically the chamber, remove eventual deposits or debris, thus avoiding the introduction in the discharging circuit of material which can cause obstructions. For correct cleaning use only demineralized water and the abrasive sponge provided (*non abrasive side* – Fig. 7).

To execute absolutely with cold chamber to avoid burns – Never use solvents, detergents, chemical solutions, descaling agents or other similar products.

Cleaning of chamber filter

Every 20 cycles or once a week

Pull up the filter (*spare code DXBA091*), paying attention not to damage it, wash with demineralized water and dry with cloth. Then replace in the seat, making sure that protrudes about 15 mm.

Cleaning of tray and tray holder

Every 20 cycles or once a week

Clean with sponge supplied (*non abrasive side*) and demineralized water.

Change of the bacteriological filter

Every 200 cycles or when it assumes a dark colour


Substitute the bacteriological filter (*Fig.A–pos.08*) turning it counterclockwise to unscrew and clockwise to screw. Use only original filter (spare code DAVA101).

Cleaning of the door gasket

Every 20 cycles or once a week

Periodically remove eventual residuals which are deposited on the circumference of the seal (*spare code DANA038*) using water and the sponge provided (*non abrasive side*), or else a moistened cloth.

To regulate the door lock pressure, it is necessary to adjust the regulator of the same door using the double function wrench (extraction tray + wrench) (*Fig.5*) provided (spare code *DANA008*). Turn anti clockwise (*pos.A*) to increase lock pressure. Turn clockwise (*pos.B*) to decrease the lock pressure.

The following steps are recommended by the manufacturer. These operations must be performed only by trained technicians authorized by the manufacturer.

<i>Replacement of the water loading filter</i>	<i>every 1 year / every 1200 cycles*</i>
<i>Control of the water filter of flowmeter inlet</i>	<i>every 1 year / every 1200 cycles*</i>
<i>Cleaning Y brass filter downstream of the radiator</i>	<i>every 1 year / every 1200 cycles*</i>
<i>EV4 (water loading solenoid valve) cleaning</i>	<i>every 1 year / every 1200 cycles*</i>
<i>Vacuum pump valves replacement</i>	<i>every 1 year / every 1200 cycles*</i>
<i>Cleaning of the radiator</i>	<i>every 1 year / every 1200 cycles*</i>
<i>Fans cleaning</i>	<i>every 1 year / every 1200 cycles*</i>
<i>Control integrity of the pneumatic circuit</i>	<i>every 1 year / every 1200 cycles*</i>
<i>Control integrity of the electric circuit</i>	<i>every 1 year / every 1200 cycles*</i>
<i>Control thermal insulation of the chamber</i>	<i>every 1 year / every 1200 cycles*</i>
<i>Control tightening screws</i>	<i>every 1 year / every 1200 cycles*</i>
<i>Control temperature probes</i>	<i>every 1 year / every 1200 cycles*</i>
<i>Safety valve check</i>	<i>Follow the guidelines of the manufacturer (see instructions)</i>
<i>Grease the moving part of the lever</i>	<i>every 1 year / every 1200 cycles*</i>
<i>Level probes check</i>	<i>every 1 year / every 1200 cycles*</i>
<i>State of tank check</i>	<i>every 1 year / every 1200 cycles*</i>

* valid the condition of the reached first

Regional/national regulations fix the mandatory of autoclave validation that checks sterilization parameters. Periodic validation must always be performed by trained technicians using appropriate measurement instruments. In the case of non-mandatory, validation is a voluntary operation which aims to certify that the autoclave is fully functional and fit for sterilization of loads for which it was designed and built.

16

MESSAGES OF ERRORS

Messages of error are emphasized through an Alfa-numerical code, consisting in a letter and three numbers.

The "A" code is relative to alarms, irregularities of the autoclave, after having carried out the remedy, if the problem persists it is necessary to request telephonic assistance.

If an alarm message is displayed (code "A") the cycle is to be considered aborted: it will be necessary to repeat all the preparation and sterilization operations.

To reset alarms and errors hold down at the same time the buttons under the bar indicated with **Reset**.

ERROR	CAUSE	SOLUTION
A 001	<i>Cycle interrupted by the user</i>	<i>Reset and restart the system</i>
A 101	<i>Vacuum not achieved in 10 min.</i>	<i>Reset and restart the system</i>
A 111	<i>Vacuum not maintained on the first phase of VACUUM TEST</i>	<i>Reset and restart the system</i>
A 121	<i>Vacuum not maintained on the second phase of VACUUM TEST</i>	<i>Reset and restart the system</i>
A132	<i>Error in the flow-meter operation</i>	<i>Reset and restart the system</i>
A133	<i>Pressure above the allowable limit during the water recall for the warmup phases</i>	<i>Reset and restart the system</i>
A 200	<i>Error in the control of the EV operation</i>	<i>Reset and restart the system</i>
A 400	<i>Error in the locking operation</i>	<i>Reset and restart the system</i>
A 401	<i>Error in the locking operation</i>	<i>Reset and restart the system</i>
A 403	<i>Error in the locking operation</i>	<i>Reset and restart the system</i>
A 405	<i>Error in the locking operation</i>	<i>Reset and restart the system</i>
A 551	<i>Pressure out of limit</i>	<i>Reset and restart the system</i>

A 637	<i>Error in the access to the memory card</i>	<i>Reset and restart the system</i>
A 651	<i>Reading of the T1 probe during sterilization above the maximum limit</i>	<i>Reset and restart the system</i>
A 653	<i>Reading of the T2 probe during sterilization above the maximum limit</i>	<i>Reset and restart the system</i>
A 661	<i>Error in the probe reading</i>	<i>Reset and restart the system</i>
A 662	<i>Error in the probe reading</i>	<i>Reset and restart the system</i>
A 701	<i>Error for failed pressure achievement during the first two phases of the warm-up</i>	<i>Reset and restart the system</i>
A 711	<i>Error for failed pressure achievement during the third phase of the warm-up</i>	<i>Reset and restart the system</i>
A 751	<i>Reading of the T1 probe during sterilization below the minimum limit</i>	<i>Reset and restart the system</i>
A 753	<i>Reading of the T2 probe during sterilization below the minimum limit</i>	<i>Reset and restart the system</i>
A 781	<i>Temperature of the 121°C cycle out of the maximum limit</i>	<i>Reset and restart the system</i>
A 782	<i>Temperature of the 134°C cycle out of the maximum limit</i>	<i>Reset and restart the system</i>
A 801	<i>Errors of out of maximum time during the first phases of discharge</i>	<i>Reset and restart the system</i>
A 811	<i>Errors of out of maximum time during the last phase of discharge</i>	<i>Reset and restart the system</i>
A 901	<i>Cycle interrupted due to power supply failure</i>	<i>Reset and restart the system Check the power supply system of the machine and of the room.</i>

In the event of one of the alarms persisting after some time, consult technical assistance.

17**SOLUTION TO OPERATING PROBLEMS**

In most cases alarms or errors are caused by lack of attention or lack of familiarity with some technical and operational aspects. You will find below the list of some anomalies with associated solutions.

17.1 The autoclave does not dry materials correctly

- Replace the bacteriological filter with a new original one
- Non original trays, without holes or with different holes made of different material have been used. It is advisable to use original trays only.
- Instruments have not been arranged correctly. Carefully comply with instructions of par.8.4

17.2 The autoclave chamber has turned white

- Change immediately the type of water used, use exclusively demineralized or distilled water, as indicated in the previous chapters and then clean the chamber .
- The whitish colour may be a consequence of the evaporation of organic materials from the instruments. Instruments should undergo a more suitable and deeper cleaning .
- Check the demineralizer installed

17.3 The autoclave chamber has bluish-green stains

- Instruments have not been cleaned correctly after the cleaning, rinse more carefully, if the stains remain apply to the phone technical service .

17.4 The sterilization cycle interrupts without any patent reason

- Check whether the autoclave is connected to the power mains through extensions, reducers or adapters, if so remove such accessories and connect the autoclave directly to the power socket.

17.5 The autoclave does not receive controls

- If the autoclave is carrying out the automatic barometric alignment, wait for the double acoustic signal after the door opening , then set functions.
- The demineralized water tank is empty, the Led indicating the minimum level is ON, fill the tank with pure water.
- The used water tank is full, the Led indicating the maximum level is ON, discharge used water .

17.6 Stains on instruments

- Instruments get yellow due to the deposit of residual chemical fluid which has fixed on the instruments due to heating. The rinsing was not suitable.
- The sterilization chamber has yellow stains. Some instruments with traces of chemical fluid have been loaded into the chamber, such fluid has fixed due to heating. The rinsing was not adequate.
- Instruments have whitish stains, the rinsing has been done with very limy water and the instruments have not been carefully dried. For the final rinsing we suggest the use of demineralized water and a careful drying of the instruments
- Instruments have blackened since the material they are made of contain a high percentage of carbon.

18 PROCEDURE FOR SERVICE AND ASSISTANCE

In the event of failure, revision, validation, contact the **Mediline ITALY s.r.l** technical helpline.

PHONE	+39 0522 94.29.96
FAX	+39 0522 94.47.98
@	service@tecnogaz.com

Assistance will assess the return at headquarters or with the intervention of a technician, and having viewed the machine in order to draw up a cost estimate, which will be forwarded to the distributor customer who will forward it to the final customer, for acknowledgement and signing.
















After receiving prior written acceptance of the cost estimate, the autoclave will be serviced and reshipped according to times indicated on the cost estimate.

In case the autoclave must be shipped for repairs, controls, reactivation, revisions, validations follow the obligatory indications below:

- Use the original packaging; if this is no longer in your possession, use adequate packaging. The merchandise travels at risk to the sender.
- Ship the autoclave only (do not include any component contained in the accessories kit).
- Carefully clean the sterilization chamber and autoclave in general before shipping. In case it arrives dirty or with residual the autoclave will be returned without being repaired, or it will be put through a cleansing action and disinfection.
- Always empty the clean water tank through the attachment located on the back of the autoclave (*FIG. 1-pos.02*)
- Always empty the used water tank through the attachment located on the back of the autoclave (*FIG. 1-pos.04*)
- Indicate by letter and insert in the package a document which indicates precisely the irregularity or service desired.
- Ship at your expense, otherwise you will be billed for shipping

All non original packaging which arrives will be disposed of.

Autoclave will be returned with new and original packing to insure maximum protection for your autoclave during shipping. Cost of packaging will be charged to client.

A RIEPILOGO RICAMBI CONSUMABILI				
18L		TRAY HOLDER		TRAY
		DPLA048		DANA049
20L		TRAY HOLDER		TRAY
		SNDA048		DHLA049
24L		TRAY HOLDER		TRAY
		DXLA348		DANA049 + DXLA348
		TRAY EXTRACTION		POWER SUPPLY CORD
		DANA008		CECG006
		SPACER		BACTERIAL FILTER
		CPAP014		DAVA101
		FILTER CHAMBER		DOOR GASKET
		DXBA091		DANA038
		PIPE FOR CHARGING WATER		PIPE FOR DISCHARGING USED WATER
		DANA099 + DXBA711 + CPRG117		DANA130
		PIPE FITTING FOR DISCHARGING USED WATER		PIPES FOR DISCHARGE UTILITIES
		CPRG096		SXBA799
		SPONGE		
		CPMG004		
CHEMICAL INTEGRATOR		200/S o 215-S	BOWIE & DICK TEST	268/S
BIOLOGICAL INDICATOR		262-S	HELIX TEST	267-S

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TECNO-GAZ S.p.A.
MEDILINE ITALIA S.r.l.