**MANUALE UTENTE
USER MANUAL
ANWENDERHANDBUCH**

NYBRA

Art. 4600/S

ITALIANO

ENGLISH

DEUTCH

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
- Sensore

R:

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Declaration Of Conformity	3
Warranty condition	4
Loss of warranty	5
Delivery terms	6
General warnings and safety rules	7
Useful contacts and addresses	7
Symbols	7
Technical Features	9
Instructions For Packaging, Transportation, Storage And Lifting	10
Machine Use	10
Operator's Safety	11
Operating Principle	12
Machine Description	13
List Of Mechanic Spare Parts	14
Installation Instructions	17
List Of Supplied Accessories	17
Control List	17
Materials Used For Sandblasting	18
List Of Tool Spare Parts	19
Tool Description And Operation	20
Instructions For The Machine Use	22
Emissions Of Toxic Substances During Operation	24

Routine Maintenance	25
Extraordinary Maintenance And Repair	27
Instructions For Machine Demolition	27
List Of Electric System Components	28
List Of Components Of The Pneumatic System	28
Summary Of Instructions To Comply With	29
Drawings And Diagrams	30

DECLARATION OF CONFORMITY

(according to ISO/IEC 17050-1)

Complies with the following pertinent regulations in force: Machine Directive 2006/42/CEE, Electromagnetic Compatibility Directive 2004/108/CEE, Low Tension Directive 2006/95/CEE and their amendments and updating.

The original declaration of conformity is provided in attached to the manual.

WARRANTY CONDITION

- 1) **DURATION:** the product is covered by a 12-month (twelve) guarantee.
- 2) **EXCLUSION:** the guarantee excludes:
 - a) faults originating from non-performance of scheduled maintenance and due to negligence or improper use of the product by the User;
 - b) periodic checks and maintenance;
 - c) repair or replacement of parts subject to wear that are fragile or have an unpredictable life, unless they prove faulty at the time of their delivery;
 - d) faults caused by personnel or originating from technical personnel transfer and transportation;
 - e) faults or damage originating from improper use or operating errors;
 - f) faults or damage originating from contamination in the water and air infeed systems, extraordinary chemical agents or electrical events;
 - g) faults or damage originating from the use of detergents, disinfecting agents, sterilizing fluids or processes that are not described explicitly in the instruction manual for operation and maintenance;
 - h) natural color change of plastic components.
- 3) **INSTALLATION AND TESTING REPORT:** an essential prerequisite for guarantee validity is the return of the installation and testing report after all its sections have been filled in by the User and the Fitter. The report must be returned within 15 days after installation, otherwise the guarantee becomes null and void.
- 4) **LIMITS:** the guarantee entitles its holder to repairs or replacement of faulty components free of charge. The right to replace the entire machine is excluded. As for components either applied or added to TECNO-GAZ products or by third parties, having their own guarantee certificate, the conditions, limits and exclusions indicated in the relevant certificates apply.

Under no circumstances shall carrying out of one or more than one repair during the guarantee term change the guarantee expiry date.
- 5) **CLAIMS:** in the event of the Purchaser making a claim relating to the application of the guarantee and the quality or conditions of the delivered equipment, the said Purchaser shall in no way suspend and/or delay payments.
- 6) **DISPUTES:** if a dispute arises on the application and interpretation of this Certificate of guarantee, the Court of PARMA (Italy) shall be the court of jurisdiction, regardless of the location in which the purchase agreement was stipulated.
- 7) **DEROGATIONS:** any special derogations from the guarantee terms described in this section shall not imply the granting to the Purchaser of any right and shall be considered as granted for a specific case only.

8) MISCELLANEOUS: for any other issue, which is not covered by this Certificate of Guarantee, the rules described in the Civil Code and existing relevant Laws in the Republic of Italy apply.

LOSS OF WARRANTY

The guarantee becomes null and void in the event that:

- a) the equipment presents damage caused by a fall, exposure to flames, liquid spillage, natural events, bad weather conditions or by any other cause that is not due to manufacturing defects;
- b) installation has not been performed in line with TECNO-GAZ instructions and by authorized personnel;
- c) the equipment is repaired, modified or tampered with by the Purchaser or by non-authorized third parties;
- d) when asking for an intervention under guarantee, the product serial number has been removed, deleted, counterfeited, etc.;
- e) the installation and testing report is not returned suitably filled in and signed within 15 days from installation;
- f) the Purchaser suspends and/or delays payment of amounts due, for any reason whatsoever, relating to the purchase and/or maintenance of the equipment;
- g) scheduled periodic maintenance or any other instruction envisaged in the manual for operation and maintenance is not carried out.

DELIVERY TERMS

1. All goods delivered must be packaged in their original packaging. If the goods are delivered in a non-genuine packaging, a charge will be applied to restore the original packaging.
2. All goods must be delivered CARRIAGE PAID.
3. The client shall ask a written authorization to TECNO-GAZ sales dept. before returning the goods by means of the "Form for Client's Returned Goods Management" (CM-M-002) in compliance with the requirements of TECNO-GAZ Quality Management System (ref. CM-P-003, CM-I-009).
4. The document accompanying the returned goods shall show:
 - the TECNO-GAZ product ode/item no.;
 - the number of the TECNO-GAZ purchase invoice or shipment document;
 - the product serial number and batch number;
 - a description of the goods returned;
 - the number of the Form for Client's Returned Goods Management" (CM-M-002) or a copy of the same attached.
5. Returned goods are considered as accepted only after they have been examined. If they prove non-conforming, they may be sent back to the consignor at TECNO-GAZ full discretion.
6. Charges will be applied to restore and test the returned product.
7. Parts replaced under guarantee must be returned, carriage paid, to TECNO-GAZ.
8. If the parts are not returned, their cost will automatically be debited to the recipient.
9. TECNO-GAZ does not accept goods returns from end users.
10. Goods sent to TECNO-GAZ for repair are managed as returned goods.

GENERAL WARNINGS AND SAFETY RULES

-  Make sure that the equipment is supplied with the correct voltage as shown in the plate.
-  Do not remove the plate.
-  Make sure that the machine is equipped with a grounding system.
-  Use a dry cloth to clean the machine.
-  Before any operation is performed, disconnect the power cable from the socket.
-  Use original spare parts only.

Non-compliance with the rules above relieves TECNO-GAZ of all responsibility.

USEFUL CONTACTS AND ADDRESSES

For technical service please apply to the dealer that invoiced the device, or directly to:

TECNO-GAZ SpA

Str. Cavalli n°4 – Sala Baganza – Parma – ITALIA.

Tel. +39 521 8380 - Fax +39 521 833391

e-mail : info@tecnogaz.com

<http://www.tecnogaz.com>

SYMBOLS

Pursuant to DIRECTIVE 2002/96/EC, this symbol shows that the product must not be disposed of as urban waste at the end of its operating life.

General warning and information for the user.

Pay attention to the instructions provided alongside this symbol.

Conformity with the Directive 2004/108/CE and Directive 2006/95/CE and further amendments.

PE

External protection.

Protection GROUNDING SYSTEM.

(Operating) GROUNDING SYSTEM.

Alternate current.

Warning: high voltage. Risk of electrocution!

Read the instruction manual carefully.

POS.

Position (number identifying a component in the figures or diagrams).

Fg.

Figure

Ω

Ohm (Electric Resistance unit).

s

Seconds (time unit).

W

Watt (Power unit).

Hz

Hertz (Frequency unit).

mm

Millimeters (Length unit).

A

Ampere (Electric Current unit).

V

Volt (Electric Voltage unit).

Pa

Pascal (Atmospheric Pressure unit).

°C

Degrees Centigrade (Temperature unit).

kg

Kilograms (Weight unit).

TECHNICAL FEATURES

MECHANIC:		ELECTRIC:	
Height	345 mm	Supply Voltage	230 V±10%
Width	415 mm	Frequency	50 Hz
Depth	530 mm	Phases	1
Weight	29 Kg ≅ 284.49 N	Fuse Current	5 A
Material	Iron	Maximum Power	1 KW
		Protection	PE
		Class	1
		Power Cord 2 M	2X1.5 mm²+GND
		Connector (Plug)	16 A + GND
SOUND POWER LEVEL:		PNEUMATIC:	
Standard Operation Noise (Iso/Dis 3746))	83.1 dB(A)	Feeding	Dry non lubricated compressed air
Maximum Noise (Iso/Dis 3746)	95.0 dB(A)	Feeding Pressure:	8 Bar
		Max Consumption 12mc/H	200 NL/Min

INSTRUCTIONS FOR PACKAGING, STORAGE, TRANSPORTATION, AND LIFTING

The sandblasting machine has a weight of 29 kg and a compact size, hence it needs no special packaging.

It is bagged into a cellophane bag and then fitted into a foam polystyrene protection to be housed into a corrugated board box sealed with adhesive tape.

During transportation the package has to be protected from bumps or shocks and humidity.

ATTENTION:

Do not disperse packaging materials into the environment; recycle reusable products and give the remaining ones to a company specialising in waste disposal.

Packed machines have to be stored in dry places at a temperature of -10 to $+40^{\circ}\text{C}$.

Machines cannot be stacked.

CAUTION:

The only seizing point for the manual lifting and transportation of the machine is the lower plane.

MACHINE USE

Is a small compressed air-fed sandblasting machine to be used for the surface treatment of objects.

Such objects have to be of correct size so as to be housed in the closed cabins and allow a safe operation (closed door).

During processing, they must not react with materials used for sandblasting, that is to say they do not have to send out gas substances (that cannot be captured by filters).

The correct operation of the machine depends on the absence of humidity inside the cabin, consequently no wet object or fluid has to be introduced to avoid sands wetting.

OPERATOR'S SAFETY

During the sandblasting machine operation the operator may run two types of risk:

1. If handling the tool when the door is open, he could aim the jet to himself and send sand into his eyes.
2. If the sandblasting treatment is carried out with the door open, the dust suction does not occur and, consequently, dusts pollute the workplace.

To eliminate such risks, the sandblasting machine is equipped with two safety systems inhibiting the machine use in case one of the door is open or the dust aspirator is not working.

The first safety device is linked to the door lock pict. 15 and includes a pneumatic sensor and a sensor indicator located on the locking lever.

If the door is not perfectly locked, the sensor indicator does not close the small continuous leak nozzle of the sensor hence the latter does not allow the air feeding to the tools. To lock the door approach it to the seal, press the lock handle downward and rotate it clockwise by 90° - see pict. 13-14 – then release the handle.

Once the operation is completed, be sure the handle cannot rotate in the opposite direction unless it is pressed downward again.

The sensor indicator is a wear part, periodical checks are suggested.

The second safety device is a check of the aspirator start and of the doors locking, only in these conditions sandblasting is possible.

The aspirator conveys the dusts deriving from the sandblasting processing to a set of filters. If the maintenance is correct, the air getting into the environment will be clean.

OPERATING PRINCIPLE

The operation of the mini sandblasting machine is characterised by two actions:

The action of the sand on the object and the action of the aspirator on dusts produced during the processing.

Two independent circuits allow the accomplishment of the above described operations.

Sand contained in the tank **pos 133 pict 1** is sucked by an ejector housed in the tool **pos 215-217 pict 1**.

Through the venturi effect, the air conveys sand into the tool nozzle and gives it a speed.

Coming out the nozzle **pos 202 pict 3**, sand is thrown to the object.

The sand speed, weight, and toughness abrade the object and thus remove part of the material and fragment the sand.

Heavy particles fall into the container. Lighter dusts are removed from the cabin through the suction mouth positioned next to the lamp.

The aspirator keeps the cabin atmosphere clean and conveys air and dusts in a filtering bag **pos 154 pict 1**.

Dusts are filtered while the air passes through the bag and is conveyed to a second filter, **pos 143 pict 1**, that further cleans the air before releasing it into the atmosphere.

MACHINE DESCRIPTION

Is a mini sandblasting machine using two types of sand.

The cabin is equipped with two tools with built-in control button **pos 215-217**, a compressed air gun to clean objects **pos 55**, a lighting lamp **pos 8-9**, a two racks for objects **pos 152-153**, two sand tanks **pos 133** (for two different types of sand) and a sealed transparent door **pos 136 pict1**, for an easy visual inspection.

Gloves are of the super-resistant type, (any type of glove can be fitted in) and are fixed to the cabin by means of a screw ring nut that allows an easy replacement.

The processing dusts are sucked by the turbine **Pos 7** and held by a filtering bag **pos 154**.

Before reaching the environment the air is cleaned by the final filter **pos 143**, this filter is made in special paper.

Upon the cabin there is the control panel **pos 117** housing suction and lighting switches, the pressure reducers used to adjust the tools pressure and the manometers displaying the pressure values.

The cabin cleaning and the contaminated sands removal is carried out by suction.

The machine is equipped with a suction nozzle hose **pos 129** to be fitted instead of the plug **pos 119** to carry out the above described operation.

The handle **pos 5**, located on the control panel is used to adjust the turbine power from zero to the maximum value.

For dusts suction during processing the turbine power is kept at a low value (enough to guarantee a good visibility inside the cabin) while for the cleaning phase the turbine power has to be the maximum one.

LIST OF MECHANIC SPARE PARTS

Pos.	Q.ty	Description	Cod.
93	2	Door	SNYA093
94	1	Control board	SNYA094
115	1	Complete furniture	SNYA115
116	1	Grid under the furniture (casing)	SNYA116
118	1	Air filters box	SNYA118
119	1	Whole gloves valve	SSKA175
124	1	Air filter plug	SSKA141
126	1	Filter holder	SSKA117
127	2	Suction pipe ends	SSKA170
128	1	Suction pipe	CP01016
129	1	Suction nozzle	CM61002
130	2	Gloves counter nut	SSKA124
131	2	Gloves ring nut	SSKA144
132	2	Glove bushing	SSKA143
133	1	Double bloc sand tank	SNYA133
134	4	Sand plungers	SSKA164
135	2	Sand tank rubber plugs	SNYA135
136	1	Door glass	SNYA136
138	1	Filtering bag suction manifold	SNYA138
139	1	Glass protection film	SNYA139
140	1	Door lock	SSKA239

141	1	Door lock protection	SSKA232
142	1	Door seal	SNYA142
144	1	Final filter plug handle	CM30002
145	1	Final filter plug	SNYA145
146	1	Final filter bushing	SNYA146
147	3	Turbine spacer	SNYA147
148	1	Control board sticker	SNYA148
150	2	Glass hinge plate (LH)	SSKA150
151	2	Gloves	CM61001
152	1	RH grid	SNYA152
153	1	LH grid	SNYA153
154	1	Suction filtering bag	SSKA253
155	2	Suction plug seal	CM50001
156	1	Filtering bag seal	CM50003
159	1	Air filter spring	SNYA159
160	2	Chromium plated knob	CM30005
161	1	Lower grid cut sponge	SNYA161
162	1	Final filter plug seal	SNYA162
163	1	Final filter plug bushing seal	SNYA163
164	1	Final filter seal	SNYA164
165	4	Sand plunger seal	CM50027
166	1	Glass hinge (RH)	SSKA243
167	1	Glass hinge (LH)	SSKA244
215	1	Whole big handpiece Ø 3.5	SNYA215
217	1	Whole big tool Ø 1.5	SNYA217

230	1	Final filter	SNYA230
330	1	Final filter pin	SNYA330

INSTALLATION ISTRUCTIONS

Sandblasting machine has to be positioned on a stand at a suitable height for the operator.
Can be located in any position provided the access to the cabin is possible.

Be sure the stand supporting the machine has a capacity suitable to the machine weight (at least 50 kg) and is steady.

Has to be connected to:

- a monophase electric socket with 220 V+/-10% and 50 Hz frequency. Such socket has to be equipped with ground wire, magnetothermal switch with differential (overload cutout) and has to distribute a maximum power of 1 Kw.
- to a DRY and NON LUBRICATED compressed air source with a minimum pressure of 7 BARS and a distribution capacity of at least 200 nl/min.

LIST OF SUPPLIED ACCESSORIES

The sandblasting machine comes with the following accessories:

N° 1 Suction hose with nozzle **Pos.128-129 Pict.1**

N° 2 Filtering bags **Pos.154 Pict.1**

N° 3 Glass protective films **Pos.139 Pict.1**

N° 1 Pritt roller (glue)

N° 1 CE Use and Maintenance manual

CONTROL LIST

The control panel includes:

Pos. 1 Green luminous switch for dusts suction

Pos. 2 Green luminous switch for cabin lighting.

Pos. 5 Suction force adjusting knob

Pos. 51 Pressure reducer to adjust the small tool pressure.

Pos. 52 Pressure reducer to adjust the big tool pressure.

MATERIALS USED FOR SANDBLASTING

Abrasives used for sandblasting are numerous and each has different features to satisfy all requirements of the dental sector.

The most commonly used are:

CORUNDUMS used to remove investment and for metal surface treatment.

ALUMINIUM OXIDES used for pickling and frames preparation.

CRYSTAL MICROBALLS for glazing.

The above-described abrasives are available in different granulometries, the selection of the abrasive depends on the grade of processing to be obtained.

Coarse grain sands incise remarkably the object and create on the surface a high roughness while fine grain sands, although removing material, smooth surface.

The abrasion grade depends on the features of the used sand with regard to the processed type of metal and on the speed at which the sand strikes the object.

The machine gives sands the maximum speed when tools are fed at a 4-5 Bars (BIOXIDE and BALLS tools) and 6-6.5 Bars (CORUNDUM tool) pressure.

The sandblasting machine can be used with all abrasives available on trade provided their size does not exceed 500 my; for very fine sands, lower than 25 my, their compatibility with the venturi suction system has to be checked for.

LIST OF TOOL SPARE PARTS

Pos.	Q.ty	Description	Cod.
176	2	Big handle	SSKA176
178	2	Big tool body	SSKA178
180	1	Big tool tip	SSKA180
181	1	Big tool tip	SSKA181
183	1	Air nozzle diam. 2	SSKA183
185	1	Air nozzle diam. 0.8	SSKA185
188	2	Control lever	SSKA188
189	2	Button rod	SSKA189
190	2	Button plug	SSKA190
191	2	Hose coupler M 6	SSKA191
192	2	Hose coupler M 5	SSKA192
193	2	Washer diam.4x9	CM20002
194	2	O-ring 2012	CM50036
195	4	O-ring 2010	CM50050
196	2	O-ring 2056	CM50038
198	2	O-ring D.	CM50068
199	2	O-ring 2018	CM50048
200	2	Spring for big tool	SSKA200
202	1	Mix nozzle diam.7x3.5	SSKA202
204	1	Mix nozzle diam 4x1.5	SSKA204
206	2	Screw TCEI M3X8	CM15037
207	2	Hose elastollan neutral	CP01008
209	2	Hose elastollan light blue	CP01005
WHOLE TOOLS			
215	1	Big tool diam. 3.5	SNYA215
217	1	Big tool diam. 1.5	SNYA217

ON DEMAND BIG TOOL DIAM. 2.5 SNYA216 (kit 4528/S)

TOOLS DESCRIPTION AND OPERATION

The tool is comprised of:

Handle **pos 176**.

Tool body **pos 178**.

Control lever **pos 188**.

Air nozzle **pos 183-4-5**.

A tip with mix nozzle **pos 180-1**.

The tool body houses an ejector that sucks the sand from the container, conveys it to the tip and ejects it through the hard metal mix nozzle, **pos 202-4**.

Each mix nozzle corresponds to an air nozzle **pos 183-5**; if the diam. 3.5 **Pos 180** tool tip is replaced with a diam. 1.5 **pos 181** tip, the diam. 2 **pos 183** nozzle too has to be replaced with a diam. 0.8 **pos 185** nozzle.

ATTENTION:

When replacing the air nozzles **pos 183-5**, be sure the O-ring **pos198** is correctly housed with a perfect tightness.

The tool tip is screwed to the body; before screwing, remove dust, if any, from threads and seal **pos 196**, using a compressed air gun, then cover with talc to reduce the seal-screw friction.

If the stick is pressed, **pos 188**, the air flows out from the nozzles thus creating the vacuum inside the tip and, consequently, sucking the sand through the white tube **pos 207**.

The maximum vacuum value is 0.7 Bar corresponding to a tool feeding of 4-5 bars (diam. 0.8 nozzles) **pos 185** and 6-6.5 Bars (diam. 2 nozzles) **pos 183**.

An increase or a decrease in the feed pressure to the above mentioned values causes a drop in the ejector yield and consequently vacuum and dust suction drop too.

A further element affecting the tool yield is the position of the tip with regard to the air nozzle; indeed, if the tip is completely screwed to the body it prevents the sand from reaching the venturi.

In the above-described condition only air flows from the tip.

After unscrewing the tip, the sand flows through the venturi and is ejected.

The highest yield point is reached by unscrewing the tip of one or two turns to the “completely screwed” position pict 4.

When the sand jet is non constant but intermittent, it means the tip position is not correct.

IN CASE OF OBSTRUCTION OF THE TOOL DUCTS:

Before unloading the container and ducts we suggest to try to eliminate the obstruction by capping the tool tip by a finger toe.

In this situation, by pressing the control lever, ducts are pressurised and sand is forced to go back to the tank **pict 4**.

In case the above described intervention does not solve the problem, remove the tool tip and check for the presence of foreign bodies.

ATTENTION:

These operations have to be performed as if a processing phase, that is the door closed and the hands fitted into the gloves.

If there is no foreign body, the obstruction must be located near the sand tank.

In this case, the intervention with open door is necessary and it is hence necessary to stop the pneumatic feeding, unload the container and identify the obstruction cause.

INSTRUCTIONS FOR THE MACHINE USE

After machine installation remove the plug **pos.124**, (rotate it to the right or to the left while pulling it upward) and be sure that the filtering bag **pos 154** is correctly housed **pos 126**, otherwise follow "routine maintenance".

The cabin of the sandblasting machine is equipped with a big tool with a 3.5 mm nozzle to process coarse sands up to 500 my and a tool with a 1.5 mm nozzle to process sands up to 100 my.

To load sands open the door on the control panel left side **pos 114** pour not more than 1 Kg of sand into the tanks **pos.133**.

The tank marked by the transparent sticker:

- "CORUNDUM": sands up to 500 my,
- "BIOXIDE": sands up to 100 my.

Put an object into the cabin, close the door, turn on the internal lamp and the aspirator.

Adjust suction through the knob **pos 5**. The correct level of suction is indicated by gloves: when they are pushed inside the cabin without deformation of their volume (they do not inflate) they indicate the suction is correct.

Pressure adjustment:

- big tool Ø 3.5 to **6-6.5 Bars**, through pressure reducer **pos.51**.
- big tool Ø 1.5 to **4-5 Bars**, through pressure reducer **pos.52**.

Fit your hands into the gloves, hold the tool by the right hand and support the object by the left hand, press the tool control lever **pos.188 pict3** and aim the jet to the surface to be processed.

The strongest abrasive force of the sand is near the nozzle, the more you move away from it the slower the speed and the bigger the jet width. see **pict5**.

When the processing is over, before turning off the aspirator, unlock the door and slightly move the glass away from the seal.

In few seconds the aspirator can remove dust completely from the cabin and at door opening no dust will disperse in the environment.

Before every sandblasting intervention, it is important to check the filtering bag

efficiency.

It is impossible to decide an amount of working hours for filter replacement since it depends on the different working conditions.

Indeed, some types of processing generate a lot of dust, for instance the investment removal. The use of very fine sands, or recycled sands, wear the filter.

The operator decides for filter replacement according to carried out processes.

In any case, we suggest to replace the bag quite often since the machine operation and the operators' health depend on its correct action.

The cabin cleaning is carried out by means of a sucking tube **pos 128**, as follows:

Remove the glove protection valve **pos 119** and fit the tube, turn on the turbine and take it to the highest speed by means of the regulator **pos 5**, and with the nozzle **pos 129** suck sand from the cabin.

AFTER CLEANING THE CABIN, THE BAG HAS TO BE REPLACED AND THE ACCESSORIES CLEANED TOO.

The door is covered by a transparent film (replaceable) to protect it against abrasion.

EMISSIONS OF TOXIC SUBSTANCES DURING OPERATION

Materials used for sandblasting, are not harmful for the skin but can be very dangerous in case of inhalation

Sandblasting generates very fine dusts, the aspirator has to convey the cabin air, saturated with such particles, to a container (filter) made of porous material (paper). The paper filter pores have a smaller size than that of dust particles, the latter are hence held while the air passes through.

Is equipped with a second with characteristics specially conceived for health protection; the second filter cleans the air from the filtering bag once more before air can reach the environment.

The above described measures to protect people's health are effective if instructions are complied with, routine maintenance and extraordinary maintenance.

If filters are performing correctly, dust emissions during the machine use are negligible.

TECNO GAZ is not responsible for the bad functioning of the cleaning circuit in case it is due to the non respect of the relevant instructions or to the unskilfulness of the machine user.

The wrong functioning of the air cleaning system due to defects of the original components has to be notified to the local retailer in due time.

ATTENTION:

DO NOT USE THE MACHINE WITH DOORS OPEN.

DO NOT USE THE MACHINE WITHOUT GLOVES.

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Routine maintenance concerns the following components: **(PICT1)**:

- 1) FILTERING BAG. **pos 154**
- 2) FINAL FILTER **pos 143**
- 3) GLOVES **pos 151**
- 4) GLASS PROTECTION **pos 139**

1) FILTERING BAG **(PICT 6-7)**

The filtering bag has to be replaced frequently; if bags are not replaced frequently you might cause the obstruction of the machine ducts and put at risk the operators' health.

It is not possible to quantify the exact duration of a filter since it depends on many factors. operator has to decide on the filter replacement according to the type of processing: consider that when the filter is close to saturation **visibility inside the cabin is reduced. In this case, do not increase the turbine revolutions but replace the filter.**

ATTENTION:

After every cabin cleaning the filtering bag has to be replaced.

The saturated filter (obstructed) does not stop dust anymore and that means a very heavy work for the final filter.

The filtering bag cannot be cleaned nor reused: it is "ONE WAY".

To replace the filtering bag remove the plug **pos 124 pict 1** (rotate to the right or to the left while pulling it upward); hanging from the plug is the filter holder.

Take the assembly in a dry place where dust cannot be harmful and wear a protective nose- and-mouth mask.

Take the plug by one hand and the holder by the other **pict 6**.

Push the holder to the plug until reaching the end-of-stroke, the rotate the holder counterclockwise to the end-of-stroke, and separate slowly the holder from the plug.

When the filtering bag neck has come out of the holder, seize it with the fingers of your left hand and extract it completely (be sure the sand contained in the bag does not tear it).

Dispose of the full bag in a suitable dump and brush the holder and the plug.

To re-assemble the filter, open a new bag, fit the plug lower side and wrap the bag around the plug **pict 7**.

Fold the exceeding part of the bag clockwise on the plug.

Fit the bag into the holder without curling it up; pins along the plug diameter have to pierce the filter paper and be fitted into the slots on the upper edge of the holder.

Fix the holder on the plug, rotate it clockwise until reaching the end-of-stroke and then lock it by pulling it to the plug.

Before fitting the filter into the machine, clean the induction manifold with an aspirator **pos 138 pict 1**.

To facilitate the filter introduction cover seals with talc.

2) FINAL FILTER (**PICT 8**)

The final filter holds the dust escaped from previous treatments and protects people's health.

The final filter has to be serviced periodically; approximately every ten replacements of the filtering bag the final filter condition has to be checked.

The final filter can be cleaned and reused; through the gun blow the air from the outside to the inside.

ATTENTION: THE GUN AIR PRESSURE HAS TO BE LOW (LOWER THAN 1 BAR) AND THE JET HAS TO BE KEPT AWAY FROM THE PAPER NOT TO TEAR IT.

The performance of the above-mentioned filter is limited and depends on the performance of previous filters, the filter can be re-used only for a reasonable amount of times.

It is not possible to assess the filter condition visually, hence it is good rule to replace it to avoid environment contamination.

When intervening on this filter, the filter box **pos 118 PICT1** and the turbine duct have to be cleaned, using an aspirator.

To dismount the FINAL filter, open the door on the right hand of the control panel marked by 'air filter', unscrew the filter plug knob **pos 144 pict 8** and extract the cartridge **pos 143 pict 8**.

When the filter is fitted into the box, respect the seal position pos 164 pict 8.

3) GLOVES (**PICT 9**)

To replace gloves, unscrew by 1/4 of turn ccw the glove locking ring nut **pos 131** positioned inside the cabin, then remove and extract the glove.

To mount gloves, fit them into the bushing **pos 130** and fit the ring nut into the joint and rotate cw by 1/4 of turn. Thanks to the special system, any type of gloves can be fitted.

4) GLASS PROTECTION (**PICT 10**)

The transparent film beneath the glass is fixed with glue (supplied with the machine) on the external perimeter.

To remove the film lift an edge and pull it.

Clean the glass with alcohol, apply the glue to the new film, apply on the hinges holes too and be sure the film perfectly adheres to the glass.

EXTRAORDINARY MAINTENANCE AND REPAIR

Extraordinary maintenance: periodical interventions to check and clean the machine internal parts.

Such interventions are performed on-site by authorised staff or at TECNO GAZ premises.

Repair interventions too can be performed on-site by authorised staff.

Apply to the nearest retailer to ask for technical service.

INSTRUCTIONS FOR MACHINE DEMOLITION

The sandblasting machine is made of iron material and inert electromechanic components.

At the machine demolition no special instruction has to be complied with since all materials are reusable.

Do not dispose of the machine in unattended areas, apply to authorised dumps.

LIST OF ELECTRIC SYSTEM COMPONENTS

Pos.	Q.ty	Description	Cod.
1	1	Bipolar switch, luminous, green, 250 V 10 A	CEQ0001
2	1	Bipolar switch, luminous, green, 250 V 10 A	CEQ0001
6	1	Turbine regulation control unit 220 V 50 HZ with potentiometer	SSKA004
7	1	Motor with turbine 220V 50HZ with brushes 0.8 KW	CPG0001
8	1	Lamp E14 15Watt 230 V	CEE0005
9	1	Colourless lamp holder	CEE0004
11	1	Double fuses holder for min. 5x20 fuses	CEF0001
12	2	Fuses min. 5x20 5A	CEF0009
30	1	Passband filter	CEZ0001
83	1	Solenoid valve	SNYA083

LIST OF COMPONENTS OF THE PNEUMATIC SYSTEM


Pos.	Q.ty	Description	Cod.
51	2	Pressure reducer 1/8"	CPF0001
52			
53	2	Manometer diam. 53 1/8" flange and bracket	CPP0001
54			
55	1	Air gun	CP02001
83	1	Solenoid valve	SNYA083
98	1	Base	SSKA145
99	1	Electropilot	CPY0027
100	1	3-way valve	CPY0028
110	1	Safety sensor	SSKA110

SUMMARY OF INSTRUCTIONS TO COMPLY WITH














- 1** Has to be fed with DRY AND NON LUBRIFIED compressed air.
- 2** When sand is loaded in the tanks, comply with indications concerning the type and the amount of sand.
- 3** Do not introduce wet objects into the cabin.
- 4** DO NOT USE THE MACHINE IF DOORS ARE OPEN
- 5** Always use the aspirator at the set speed.
- 6** Be sure the filtering bag is always present.
- 7** **Replace the filtering bag at frequent intervals.**
- 8** **Always replace the filtering bag at each cabin cleaning.**
- 9** At each replacing of the filtering bag clean the plug, the holder, and seals.
- 10** Cover seals with talc.
- 11** Be careful not to tear the filtering bag.
- 12** Periodically check the final filter, clean and replace it.
- 13** Do not forget to apply the transparent film to the glass.
- 14** The feeding pressure of the big tool diam. 3.5 is 6-6.5 bars
- 15** The feeding pressure of the big tool diam. 1.5 is 4-5 bars
- 16** DO NOT USE THE MACHINE WITHOUT GLOVES.
- 17** IN CASE OF ANY INTERVENTION INSIDE THE MACHINE, UNPLUG THE MACHINE FROM POWER SUPPLY.
- 18** **In case of interventions on tools with doors open, disconnect the pneumatic feeding and discharge the machine system. (DETACH THE TUBE)**

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TECNO-GAZ S.p.A.