
# Europa B xp 24

Istruzioni per l'uso / Instructions for use / Mode d'emploi / Instrucciones de uso /  
Bedienungsanleitung

## Gentile Cliente

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## Dear Customer

The EUROPA B xp 24 is designed for dental instruments sterilization.

Before starting up the unit **read carefully the following instructions** to familiarise yourself with correct operations.

**The autoclaves must be used single and exclusively for the uses previewed from the constructor.**

**In order to avoid injuries to persons and damage to the equipment, pay special attention to Safety Notes.**

**CE** 0434

**MEDILINE  
ITALIA S.r.l.**

**Conformément aux dispositions de la Directive 93/42/CEE**

**COMPANY  
WITH QUALITY SYSTEM  
CERTIFIED BY DNV  
UNI EN ISO 9001  
= UNI CEI EN ISO 13485 =**

# List of Contents

- 0.-** Use and function foreseen for autoclave
- 1.-** Operating Controls and Displays
- 2.-** Accessories
- 3.-** Sterilization Cycles
  - 3.1.- Overview – Cycles description
    - 3.1.1.- Sterilization Cycles
    - 3.1.2.- Test Cycles
  - 3.2.- Sterilization Cycles Table
  - 3.3.- Test Cycles Table
- 4.-** Operations
  - 4.1.- Loading the unit
    - 4.1.1.- Disinfecting
    - 4.1.2.- Cleaning
    - 4.1.3.- Wrapping
    - 4.1.4.- Loading
  - 4.2.- Switching ON
  - 4.3.- Terminating a Program
  - 4.4.- End of Program – Unloading
  - 4.5.- Switching OFF
  - 4.6.- Programming
    - 4.6.1.- Language Selection
    - 4.6.2.- Hour Selection
    - 4.6.3.- Date Selection
    - 4.6.4.- Atmospheric pressure visualization
    - 4.6.5.- Water Feeding Selection
- 5.-** Error messages
  - 5.1.- Overview
  - 5.2.- List of Error Messages

## **6.- Maintenance and Care**

- 6.1.- Waste water discharge
- 6.2.- Replacing the Sterile Filter
- 6.3.- Cleaning the Chamber
- 6.4.- Regulation door
- 6.5.- Replacing Paper in Printer
- 6.6.- Printer Status
- 6.7.- Unpackaging and Installing
- 6.8.- Water feeding

## **7.- Safety**

- 7.1.- Safety Labels
- 7.2.- Safety Devices
- 7.3.- Safety Notes
- 7.4.- Waste Disposal

## **8.- Technical Data**

### **Annexes**

Water Quality Table

Cycle Diagrams

Declaration of Conformity

Guarantee

## 0 Use and function foreseen for autoclave

Autoclave mod. EUROPA B xp 24 is able to sterilize the 3 kind of charge indicates in the European Norm EN 13060, specifically :

<b>FERROUS MATERIALS OR SOLIDS</b>	Instruments or devices without cavity or obstacles for steam penetration
<b>POROUS PARTS</b>	Simple or composed materials they can absorb the fluid ( textures, overalls, gauzes, medications, etc...)
<b>HOLLOW PARTS</b>	Materials or devices with cavity, obstructions, etc... European Norm divides these kind of parts in 2 types, it depends by dimensions of cavity. Indicatively <b>TYPE B</b> : tubes with large hollow <b>TYPE A</b> : turbines, devices and strick with small holes or small dimensions

**Autoclave must be used, only and exclusively, for the sterilization of instruments and materials compatible with steam sterilization systems.**


This equipment can be used in the dental, medical, aesthetic fields and, generally, in all the fields where the sterilization of the instruments and materials is made

Always verify that the load to be sterilized can support the cycle temperature chosen.

# 1 Operating Controls and Displays


## 1.1.- Closed Door Front View

- |   |                       |    |                     |
|---|-----------------------|----|---------------------|
| 1 | Water level max       | 9  | Prion Cycle         |
| 2 | Water level min       | 10 | 134°C Wrapped Cycle |
| 3 | Waste water level max | 11 | 121°C Flash Cycle   |
| 4 | Helix Test            | 12 | Cycle Select        |
| 6 | 121°C Unwrapped Cycle | 13 | Water Inlet         |
| 6 | 121°C Wrapped Cycle   | 14 | ON - OFF            |
| 7 | 134°C Unwrapped Cycle | 15 | Open Door           |
| 8 | Vacuum Test           | 16 | Mains switch        |

## 1.2.- Open Door Front View


- 1 Manual filling for emergency water
- 2 Microbiological Filter
- 3 Waste water out

### 1.3 – Back View

- 1 Discharge of used water
- 2 Exceeding used water – condensate discharge
- 3 Discharge clean water – connection deionizer
- 4 Exceeding clean water

- 5 Plug for Electrical power ONLY AND EXCLUSIVELY FOR DEIONIZER
- 6 Safety valve
- 7 Fuse 12A
- 8 Fuse 12A
- 9 Fuse 500mA
- 10 Connection RS232
- 11 Principal electrical alimentation

## 2 Accessories

### 2.1.- Trays holder

### 2.2.- Tray

## 3 Sterilization Cycles

### 3.1.- Overview – Cycles description

The unit sterilizes and dries dental instruments according to the following cycles. The instruments may be wrapped and then positioned on the trays into the sterilization chamber.

Sterilization cycles are classified in two categories, one of normally cycles used for sterilizing instruments, and the other for testing the autoclave.

#### 3.1.1.- Sterilization cycles

- **121° C unwrapped** is related to a solid load without wrapping
- **121° C wrapped** is related to a wrapped solid load
- **134° C flash** is related to a cycle with fast drying
- **134° C unwrapped** is related to a solid load without wrapping
- **134° C wrapped** is related to a wrapped solid load
- **134° C Prion** is a cycle apt to sterilize CJD ( Creutzfeldt -Jacob Disease) contaminations.

#### 3.1.2.- Test Cycles

- **134° C Helix** is used for testing steam penetration into a hollow load. Helix Indicator System may be used.
- **134° C Bowie & Dick** is used for testing steam penetration into a porous load. Bowie & Dick packaging may be used.
- **134° C Vacuum** is used for testing vacuum efficiency of the autoclave

### 3.2.- STERILIZATION CYCLES TABLE

<b>CYCLES</b>	TIME EXPOSED (Minutes)	DRYING TIME (Minutes)	CHARGE* (Kg)	CHARGE** (Kg)	MINIMUM WORK PRESSURE ( bar, differential pressure )	VACUUM PHASES N°
<b>121° unwrap.</b>	18	12	1 + 0.5	0.50	1.04	3
<b>121° C wrapped</b>	18	13	0.75 + 0.35	0.50	1.04	3
<b>134° C flash</b>	4	4	0.75 + 0.35	0.50	2.02	3
<b>134° C unwrap.</b>	5	12	1 + 0.5	0.50	2.02	3
<b>134° C wrapped</b>	5	13	0.75 + 0.35	0.50	2.02	3
<b>134° C Prion</b>	20	12	0.75 + 0.35	0.50	2.02	3

\* The values indicate are the maximum load for every tray, for ferrous materials or solid or hollow part

\*\* The values indicate are the maximum load for every tray, for porous parts  
Is not possible to add the maximum loads indicates in the table over.

Warmig up time and vacuum fractionation can vary, because of autoclave's conditions, between 25 and 30 minutes. You have to an this time to sterilization's times showed in table 3.2

### 3.3.- TEST CYCLES TABLE

<b>CYCLES</b>	<b>STERIL. TIME ( Minutes)</b>	<b>DRYING TIME (Minutes)</b>	<b>MAX LOAD per Tray (Kg)</b>	<b>MINIMUM WORK PRESSURE ( bar, differential pressure )</b>
<b>134° C Helix/ Bowie&amp;Dick</b>	3,5	4		2.02
<b>134° C Vacuum</b>				-0,9

## 4 Operation

### 4.1.- Loading the Unit

A correct loading is very important to achieve a perfect sterilization, and the following operations are suggested.

#### 4.1.1.- Disinfecting

A disinfecting tank may be used; fill the tank at the half volume level with disinfecting solution, than insert trays and shut.

Manipulate contaminated instruments only using pliers.

A 30 minutes time is normally used to achieve a good disinfecting.

#### 4.1.2.- Cleaning

An Ultrasonic Cleaning is recommended, because it is very useful to remove organic pollution and up to 90% of microbiological contamination.

Be sure that all instruments are fully below the liquid level.

Do not pile up the instruments to avoid damage due to rubbing.

After cleaning a drying process is normally used.

#### 4.1.3.- Wrapping

Wrapping is obtained using a Thermosealing Unit, following its Operation Manual.

Be sure that a complete sealing is reached in all items.

#### 4.1.4.- Loading

It is very important that steam may enter into all instruments in order to obtain there sterilization.

Place the instruments over the trays and avoid to overlap its.

The maximum load for not porous parts is:

- 0.75 Kg ( packed cycles ) or 1.00 Kg ( unpacked cycles )  
for the trays with size 184 x 17 x 286 mm

- 0.35 Kg ( packed cycles ) or 0.50 Kg (unpacked cycles )  
for the trays with size 184 x 17 x 140 mm

If you trays don't match the size above then the maximum load is:

- 0.75 Kg + 0.35 Kg = 1.10 Kg ( packed cycles ); 1.00 Kg + 0.50 Kg ( unpacked cycles )  
for the trays with size 184 x 17 x 460 mm

## 4.2.- Switching ON

- Switch ON the pushbutton MAINS on the right side of the autoclave.

- Click the pushbutton labelled 


You may choose the sterilization program among the ones labelled on the frontal panel. One clicking makes a step to change one program into the successive.

- The words OPEN THE DOOR are displayed. Open the door.
- 5" after the door is opened, two "beep" warn that atmospheric pressure compensation is acquired . DOOR OPEN is displayed.
- Shut the door. DOOR CLOSED is displayed.

- Click the pushbutton labelled 

LOCKED DOOR is displayed and the chosen cycle is started.



## 4.3.- Terminating a Program

A running program may be interrupted clicking the pushbutton labelled  and keeping it pushed 2 seconds.

Wait till displayed pressure is 0 bar.

An error message A001 ( cycle interrupted) is displayed.

To reset error keep clicked together the below noted pushbuttons for 10 seconds:


### **WARNING !**

To avoid injuries take care not to touch hot instruments !

#### 4.4.- End of Program – Unloading

At the end of a program an acoustic signal warns the operator.  
If the cycle has been successfully completed END OF CYCLE is displayed.

Click  and the door is unlocked and may be open.

#### **WARNING !**

To avoid injuries take care not to closely approach to the open door, because some hot steam may be present.

Use protection gloves to manipulate the instruments.

#### 4.5.- Switching OFF


Switch off the MAINS switch.

If you do not switch off, the autoclave waits 60 minutes than enters a Stand By condition, in witch all indicating Leds and the display are switched off.

To reset the machine click one of the pushbuttons.

#### 4.6.- Programming

##### 4.6.1.- Language selection


Press the  and in the meantime the MAINS switch. Keep pressing

until on the display the selected language is visualised.

Press the  to modify the language [forward].

Press the  to modify the language [backwards].

##### 4.6.2. - Hour Selection


Press the  to the hour visualisation.

Press the  to modify the hour [forwards].

Press the  to modify the hour [backwards].


#### 4.6.2.1.- Minutes Selection

Press the  to the minute visualisation.


Press the  to modify the minutes [forwards].

Press the  to modify the minutes [backwards].


#### 4.6.2.2.- Seconds Selection

Press the  to the seconds visualisation. Do not modify.


#### 4.6.3.- Date Selection

Press the  to the day visualisation.

Press the  to modify day [forwards]

Press the  to modify day [backwards]


#### 4.6.3.1.- Month Selection

Press the  to the month visualisation.


Press the  to modify month [forwards]

Press the  to modify month [backwards]

#### 4.6.3.2.- Year Selection


Press the  to the year visualisation.

Press the  to modify year [forwards]

Press the  to modify year [backwards]


#### 4.6.4.- Atmospheric Pressure Visualisation


Press the  to the visualisation of the atmospheric pressure stored.

**DO NOT MODIFY!**

#### 4.6.5.- Water Charge Selection

Press the  to the water charge selection (default: **charging by pump**)

Press the  to modify in **charging by water network**

Press the  and keep pressing until the “beep” sound.

The selected data are storage. Switch off the autoclave by MAINS switch.



## 5 Error Messages

### 5.1.- Overview

Messages of error are emphasized through an Alfa-numerical code, consisting in a letter and three numbers.

The “E” suffix is relative to errors of the operator, which can be corrected by same operator. The “A” suffix is relative to alarms, irregularities of the autoclave, after having carried out the remedy, if the problem persists it is necessary to request telephonic assistance.

**To reset alarms and errors, maintain pressed the buttons here indicated until the intervention of a beep:**

## 5.2.- List of Error Messages

<b>ERROR</b>	<b>CAUSE</b>	<b>REMEDY</b>
E 200	Low water level	Fill tank
E 250	Low level water cleaned up with loading from deionizer	Fill tank
E 300	Maximum level used water	Empty tank
E 401	Door open	Close door with attention
E 451	Door unlocked	Open door and close again
E 461	High temperature in chamber for Vacuum test (above 40°)	Shut off and wait 10 minutes
A 250	Cycle begins without water	Control Waterclean device
A 001	Interrupted cycle	Reset and restart cycle
A 011	Irregular pc display	Shut off and restart autoclave
A 101	Vacuum not reached in 10min.	Reset and repeat the cycle
A 111	Vacuum not maintained for first phase VT	Reset and repeat the cycle
A 121	Vacuum not maintained for second phase VT	Reset and repeat the cycle
A 401	Hatch opened during the cycle or problem closing	To control the closing system
A 501	Pressure not maintained	Reset and repeat the cycle
A 551	Irregular 1° temperature sensor	Reset and repeat the cycle
A 552	Irregular 2° temperature sensor	Reset and repeat the cycle
A 641	Overly high pressure	Reset and repeat the cycle
A 651	Overly high temperature	Reset and repeat the cycle
A 661	Irregular temperature comparison	Reset and repeat the cycle
A 701	Working pressare not reached	Reset and repeat the cycle
A 751	Low temperature	Reset and repeat the cycle
A 761	Irregular pressure sensor	Reset and repeat the cycle
A 801	Irregular unloading pressure	Reset and repeat the cycle
A 901	Main voltage too low	Inform electric energy provider
<b>SUBSTUTUTE RESIN CLEN WATER resin finished</b>		Follow the instruction in Paragraph 14.2 of the manual
<b>D</b>	The printer's door is not closed	Try to close the printer's door
<b>P</b>	Paper is missing	Insert paper
<b>L</b>	There is low voltege on the line	To try to connect the machine to the net electrical worker in an other point; if the problem persists to contact a electrician

## 6 Maintenance and Care

### 6.1.- Waste water discharging

A Led is lighted on front panel when waste water level in tank reaches the allowable maximum.

In this case the autoclave is inhibited.

Connect discharge tube to output connector (see fig. 2) and discharge water.

### 6.2.- Replacing the Sterile Filter

Replace Sterile Filter every three months, or in a shorter time if its colour is dark.

### 6.3.- Cleaning the Chamber


Clean the chamber with a slight abrasive cloth to remove solid particles that may be present. Use only demineralised water to moisten the cloth.

### 6.4.- Door Regulation

To regulate the door lock pressure, it is necessary to adjust the regulator of the same door using the double function wrench (extraction tray + wrench) provided.

Turn anti clockwise (fig.A) to increase lock pressure.


Turn clockwise (fig.B) to decrease the lock pressure.

## 6.5.- Replacing Paper in Printer

Thermographic paper roll 57 mm wide. Thermographic paper can be printed on one side only. Open the printer cover (1), place the paper roll by paying attention to the rotation direction (2), pull the paper until it comes out and close the cover (3); the printer is ready for printing (4).

## 6.6.- Printer Status

If the door of the printer is not perfect closed, on the display appears the letter "**D**". To open and close again the door of the printer.

If it is lacking the paper, on the display appears the letter "**P**". To open the door, to introduce the paper and close again the door.

## 6.7.- Unpackaging and Installing

Autoclave is supported with a wooden pallet and is inclosed in a corrugated board, inserted on the top of the machine.

The packaging is tied with two plastic straps to the pallet.

To unpackage the autoclave first cut the straps, than remove the board with a vertical movement.

Two persons are required to lift the machine, as its weight is 60 Kg.

Take care to the following notes:

- Install the autoclave only in rooms in which only authorized personnel is admitted
- The autoclaves do not have to be installed in places to risk upgrade them of outbreak and/or fire
- Autoclave room needs a good lightening
- Autoclave room needs a good ventilation
- Horizontal support plane load up to 70 Kg
- A 50 mm gap between the autoclave and the wall is needed
- Avoid the presence of water sprays towards the cover of the autoclave
- Do not put anything on the cover
- Do not install the machine near heat sources
- Avoid bending mains cable
- Take care that waste water tank level be at least 100 mm below machine bottom
- Take care that a.c. power supply be in compliance with safety rules.
- Take care that a.c. power supply be 230 V 2000 VA
- Do not insert mains connector through adapters


## 6.8.- Water feeding

**WARNING:** The autoclave can work only with the deionizer supplied as accessory; on this matter, read the user's manual of the deionizer.

# 7 Safety

## 7.1.- Safety labels

## 7.2.- Safety Devices

Safety devices are the following:

- safety pressure valve rated 2.4 bar  $\pm 10\%$
- electromagnetic locking door to avoid opening the door when a cycle is running

## 7.3.- Safety Notes

In quality of constructor Mediline Italia Srl is responsible of the product commercialised, to senses of the enforced laws . The responsibility expires at the moment the device, or parts of the device, is tampered by a non authorized staff or if the parts are not original.

**The location where the autoclave is installed must not be at risk for potential explosion and/or fire.**

The autoclaves must be installed in an atmosphere in compliance with legislative requirement in force on the Country use.

## 7.4.- Waste disposal

This product is subject to Directive 2002/96/EC of the European Parliament and the Council of the European Union on Waste of Electrical and Electronic Equipment (WEEE) and, in jurisdictions adopting that Directive, is marked as being put on the market after August 13, 2005, and should not be disposed of as unsorted public waste. Please utilise your local WEEE collection facilities in the disposition and otherwise observe all applicable requirements.

## 8 Technical Data

### 8.1.- Mechanic

- Working temperature +5 +40°C
- Max altitude 2.000 m
- Max relative humidity at 30°C 80%
- Max relative humidity at 40°C 50%
- Overall dimensions ( wide x high x deep) 510x390x730 mm
- Open door space 370 mm
- Weight ( empty tanks) 58 Kg
- Weight (full tanks) 67 Kg
- Pressure on bearing 2058 N/m<sup>2</sup>
- Volume 0.15 m<sup>3</sup>
- Colour RAL 9016
- Materials Stainless Steel AISI 304 Steel FeP01
- Noise < 70 db A

### 8.2.- Electric

- Supply Voltage 230 V a.c. + / - 10 % single phase
- Power 2.000 W
- Frequency 50 / 60 Hz
- Mains cable 2 + 1 x 1mm<sup>2</sup>
- Fuses 6.3 x 32 - 12 A
- Heat production 3.6 E +6 J \ hour

### 8.3.- Chamber

- Max working pressure 2.2 bar (relative)
- Max vacuum - 0.9 bar (relative)
- Max temperature 138 °C
- Material Stainless Steel AISI 304
- Dimensions diam. 245 x 500 mm

### 8.4.- Clean Water Tank

- Volume 4.5 l
- Allowable cycles read the user's manual of the deionizer
- Material Polietylene

### 8.5.- Waste Water Tank

- Volume 4.5 l
- Allowable cycles 2
- Material Polietylene
- Waste water max temp. 50° C

### 8.6.- Bacteriological Filter

- Diameter 56 mm
- Filtering efficiency 0.3 micrometer

### 8.7.- Trays Holder

- Material Anodized Aluminium
- Dimensions (wide x high x deep ) 192 x 200 x 470 mm
- Provided number 1

### 8.8.- Trays

- Material Anodized Aluminium
- Dimensions (wide x high x deep ) 184 x 17 x 286 mm  
184 x 17 x 140 mm
- Provided number 4+4

#### Trays ( optional )

- Material Anodized Aluminium or Stainless Steel
- Dimensions (wide x high x deep ) 184 x 17 x 460 mm
- Provided number 4

## GUARANTEE

1. The product is guaranteed for a period of 12 months.
2. The guarantee begins from the date of delivery of the machine to client proven by the restitution of the guarantee label attached to the tool, filled in, stamped and signed by the dealer.
3. Repair or substitution during guarantee of a part is effected upon the absolute judgment of TECNO-GAZ company. It does not include the transfer of employees, packaging and shipping costs, which are charged to the customer.
4. Excluded from the guarantee are light bulbs, fuses, and components subject to normal wear.
5. The right to complete substitution of the machine is not recognized
6. The guarantee does not cover damage caused by improper maintenance, negligence, incompetence or other causes not related to the constructor.
7. The guarantee does not include any compensation for direct or indirect damage of any nature towards persons or things, caused by the inefficiency of the machine.
8. TECNO-GAZ is not responsible for breakdowns or damage caused by inappropriate use of the product or in case ordinary maintenance has not been carried out or elementary principles of correct maintenance (negligence) have been ignored.
9. There is no acknowledgment of compensation for machine stoppage.
10. **Mounting, service, overhaul, and control of tools have to be carried out exclusively by service centers authorized by TECNO GAZ to be identified by the card issued after the participation in the certification course.**
11. *Any intervention on the tool carried out by unauthorized staff excludes the right to guarantee as well as the direct liability of the producer. This also applies to interventions performed out of the guarantee period.*
12. Components replaced under guarantee must be returned to TECNO-GAZ postage paid. Failure to return the part within 10 days the cost of the part will be charged to the buyer
13. TECNO-GAZ does not accept return from final users, unless specifically authorized in writing.
14. Return to TECNO-GAZ for repairs, must be controlled by the dealer or by service centers chosen by the final user in compliance with commercial procedures CM-P-003 entitled "CONTROL OF REQUESTS FOR RETURN, REPAIR, SUBSTITUTION OF TECNO-GAZ PRODUCTS AND CLIENT COMPLAINTS".
15. Return to TECNO-GAZ must be documented and authorized as per internal procedures.
16. **Products returned to TECNO-GAZ must have attached documentation of authorized return and a document describing the breakdown, otherwise the diagnostic cost will be charged.**
17. All products in repair must be shipped to TECNO-GAZ postage paid and adequately packaged (we suggest using the original packaging)  
The goods travel at sender's risk.

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## Annexes

**TABLE OF QUALITY LEVEL IN ACCORDING TO DIRECTIVE DIN EN 285**

<b>CEN STANDARD DIN EN 285</b>			
Evaporation residue	≤	10	mg/l
Silicium oxide, SiO <sub>2</sub>	≤	1	mg/l
Iron	≤	0.2	mg/l
Cadmium	≤	0.005	mg/l
Lead	≤	0.05	mg/l
Rest of heavy metals except iron, cadmium, lead	≤	0.1	mg/l
Chloride (Cl')	≤	2	mg/l
Phosphate (P20s)	≤	0.5	mg/l
Conductivity (at 20°C)	≤	15	μs/cm
PH value (degree of acidity)		5 to 7	
Appearance		Colourless clean without sediment	
Hardness (E Ions of alkalin earth)	≤	0.02	mmol/l

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
Made by MEDILINE ITALIA S.r.l. is conform to the provision of law:

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## TECNO-GAZ S.p.A.