Tecno-Gaz S.p.A. Strada Cavalli, 4 • 43038 Sala Baganza (PR) Italy Tel. +39 0521 8380 • Fax +39 0521 833391 • info@tecnogaz.com

Special Edition

TECNO-GGAZANA

Conscious sedation

New delivery System & New Masks









Sedative Analgesie

FOCUSON

Inhalation sedation utilizing nitrous oxide oxygen has been a primary technique in the management of patients' preoccupations and anxieties for more than 150 years and remains so today.





Master Flux Plus

The better the performance, the greater the patient's comfort.







0039 0521.8380



General Information on **Conscious sedation**





A smany as 10 to 30% of adults and children may have some form of fear or anxiety related to dental treatment. There is substantial evidence that these patients will benefit from conscious sedation with nitrous oxide (N_2O) and that this form of sedation is extremely safe and efficient in the trained "dental practitioner's hands."

Conscious Sedation can be defined as: "A technique in which the use of a drug or drugs produces a state of depression of the central nervous system enabling treatment to be carried out, but during which verbal contact with the patient is maintained throughout the period of sedation."

CED writes the 11/05/12 in a CED report:

nhalation sedation utilizing nitrous oxide and oxygen has been a primary technique in the management of patients' preoccupations and anxieties for more than 150 years and remains so today. Nitrous oxide inhalation sedation, when it is supported by behavior management techniques is efficacious for children and adults. A 2008 Cochrane review reported favorable changes in behavior or anxiety when N₂O was used. This has been described as the "standard technique" for pediatric dentistry (NICE 2010).

The technique is safe, valuable and effective for dental procedures and will allow the majority of apprehensive dental patients to be successfully sedated and treated in a much more comfortable and stress-free environment.

Administered properly, by accredited dental practitioners with well-maintained equipment and appropriately trained assistants, the technique has an extremely high success rate and must be part of dental equipment as a fundamental instrument for the management of pain and anxiety.



2

Tel: +39 0521 8380

Special Edition

3

Tecno-Gaz Presentation







(+)

- Tecno-Gaz established in 1979
- Based in Parma, Italy
- Masterflux produced **since 1995**
- Sole manufacturer of conscious sedation devices in Europe
- 12,000 units installed all over the world

Universities: Master Flux Plus used and improved in the University of

Parma, LKH - Univ, Klinikum Graz, Austria, Charité Universitätsmedizin Berlin, Universidad Internacional de Cataluña (UIC), Universidades de Barcelona, Facultad de Odontologia de Murcia, Complutense de Madrid, Università di Torino, Ancona, Bologna, Bolzano, Roma, Napoli, Catania, Sassari, Poznan.



single product, but include a range of services, skills and expertise.



Sedative Analgaesia Masterflux Plus

Conscious sedation

he controlled administration of oxygen and nitrous oxide by inhalation induces a relaxed feeling and reduces emotional stress. Anxiety and the perception of pain are reduced, increasing the amenability of the patient and facilitating collaboration with the practitioner, while maintaining unaltered the vital reflexes.

Why conscious sedation?

omplex operations and increased invasiveness have significantly lengthened the sessions and require a new approach to the patient. With conscious sedation, it is possible to decrease surgical stress, increase patient compliance and, consequently, productivity.



4

~

Master Flux Plus, the conscious sedation machine

aster Flux PLUS is the unique conscious sedation machine entirely researched and developed in Italy.



All parts are manufactured by Tecno-Gaz and are subject to scrupulous individual inspection. To give the machine a positive visual and psychological impact on the patient, the components are enclosed in an elegant unit which contains oxygen and nitrous oxide bottles with relative reducers and attachments.

Tel: +39 0521 8380

Safety devices

Minimum oxygen concentration

There is always a percentage of 30% oxygen present, thereby safeguarding against the risk of HYPOXIA and ANOXIA.

N₂O cut-off

Cuts off the N₂O flow **when the O₂ pressure is zero** (no oxygen)

2

O₂ Emergency

When we press the O₂ device, we cut off the flow of N₂O, **delivering 100% oxygen both to the control balloon and to the patient undergoing treatment.**

Anti-expired gas valve

A special no return valve prevents expired gases from being recycled in the circuit. Additionally, a special device enables the patient to breathe ambient air in the event of a gas failure.



Automatic flow control

The system has been designed to make flow administration quick and easy. Turning the O₂ knob, 1 we set the flow required by the patient in litres per minute, depending on his/her respiratory volume. The special adjustment knob 3 enables the percentage of nitrous oxide delivered to be automatically varied from 0 to 70%. **The columns indicating the flows** have been designed for fast and easy reading.

The couplings have been fitted on the back of the flow control box and their attachments are physically different to safeguard against the risk of inverting the tubes for the

5

NEWS

two gases. Yet another precaution to guarantee safety in the use of the machine.



Master Flux Plus Electronic





Safety devices

he MF ELECTRONIC flowmeter box is a sophisticated electronic engineering project equipped with a proportioner-driven delivery system. An electronic method has been chosen, in which the flow can be precisely adjusted with the touch of a button. This has allowed us to obtain a high precision product which is very reliable and simple to use. Master Flux Plus Electronic is easy to use. Its technical features make it extremely adaptable. Its effective safety devices avoid every possibility of error, and therefore any negative consequences for the patient. The displayed graphic contents are easy to read.





Tel: +39 0521 8380

FILIX PILIS



7



N₂O cut-off

Cuts off the N₂O flow when the O₂ pressure is zero (no oxygen).

Anti-expired gas

A special valve prevents expired gases from being recycled in the circuit. Additionally, a special device enables the patient to breathe ambient air in the event of a gas failure.

Automatic flow control

- The system has been designed to make the flow administration quick and easy. Turning the O₂ knob, we set flow required by the patient in litres per minute, depending on his/her respiratory volume. The special adjustment knob enables the percentage of nitrous oxide delivered to be automatically varied from 0 to 70%.
- The display interface has been designed to ensure a fast read of the gases flow rates.

Tel: +39 0521 8380

Mixed Gas Outlet

Master Flux Plus Electronic





The pressure reducers connecting the flowmeter box to the cylinders are mounted inside the cabinet, out of the sight of the patient; the mounting design is safe with no risk of exchanging the two gas reducers to both the cylinder, and to the flowmeter box.



8

197

Tel: +39 0521 8380

TECNO-GAZ NEWS Fhree totale 8.0 Portable, with 1.PM Preset adult patient 4-hour batteries Party NO 5.0 7.2 - 0.8 45 . 90 settings MENU \$70P hunse totals ±PM Preset paediatric 15.0 22 H MI The 18 40 patient settings 50 . 70 3.5 - 1.5 Totally safe 📿 START . MENU Fhunso totale R.0 Oxygen depletion warning LPM N₂O depletion warning Firsts 900 Therewice: 50 -14/15 Per la Completely configurable MENU START • Flusso.tolale Possibility to record new patients (Aidt 8.0 apportant with progressive numbers and to ±.PM Padalitics retrieve performance data (time 11110 Page. Time (+ 00) These 400 101 1.10% 7.2 ... 0.8 and amount) subsequently Ovoita START MENU ø 1





Web: www.tecnogaz.com Email: marketing@tecnogaz.com

9

Intellivery System & Masks

masks children and adolescents

Workspace

New

The IntelliFlux masks can be adjusted to 180° to work on the lateral part of the mouth and guarantee the ideal access to the work space of the dentist.



Patient Comfort



forehead of the patient

Support area for the





Silicone profile suitable for any type of face. Leaves no marks, gives great comfort, and reduces N₂O emissions in the environment Gas outlet with backflow valves

> Gas inlet with backflow valves

Ease of tube replacement





Faster conscious sedation time

New easy-breathe membranes



Mounting kit for dental chair

10



The New Circuit

Ergonomics

- Perfectly suitable for 4-hand work.
- Circuits tubes located above the patient's head avoiding obstacles to the treatment.
- Can be moved in all positions.

Pleasant treatment

- The patient is sedated, enjoying a peaceful experience, breathing through a mask that will leave no marks on his face.
- The frontal holder and the rubber elastic band ensure perfect stability.



Fast and efficient conscious sedation

• Perfectly sealing mask adapting to all face shapes and curves for a faster and more efficient sedation treatment.

No gas dispersion in dental clinic environment

- Many particular specifications contribute to prevent dispersion such as:
- Backflow valve.
- Adaptable connector to dental unit aspirators.

11

- New easy-breathe membrane.
- Face-adaptable, single use silicone shape.





Prevent Cross-Infections Best practice modalities

Tecno-Gaz, a leader in sterilization, gives information on the best practice in the usage of sedation masks to prevent cross contamination as recommended by European regulations: EN 13060, 17665-1, EN556-1, and the Spaulding classification. Spaulding defines the criticity (critical, semi-critical or not critical) of instruments and objects used during treatment, depending on the use and the possible contact with patient skin or gum and tissues.

What is disposable, and what can be sterilized?



Tubes are defined as non critical. However, Tecno-Gaz proposes 3 types of tubes:

Single use (optional)
Sterilizable at 121°C (Standard)
Sterilizable at 134°C (optional)

Masks

Masks are defined as semi-critical. They can be disinfected or sterilized at 134° C

Rubber parts

Sterilizable at 134° C. Also exists in single use lots.



Tel: +39 0521 8380

TECNO-GAZ

Accessories for Master Flux Plus



Mask

L (XXS) Art. No. MF835ZMF



Mask 2 (XS) Art. No. MF836ZMF



Mask **3** (S) Art. No. MF837ZMF



Mask

4 (M) Art. No. MF838ZMF

One-size silicone head cushion and silicone mask pads in various sizes

Autoclavable at 134°C. Minimum: 30 cycles

1 (XXS) Art. No. 2ZMFA0050 2 (XS) Art. No. 2ZMFA0051



4 (M) Art. No. 2ZMFA0054

Holding strap with **Plastic hook**

Art. No. 2ZMFA0046



Mounting kit for dental chair

Art. No. 3MEDT0007+1ZMFA0064



Breathing balloon

Art. No. 3MEDN0001

Circuit tube

Tube length: 1.8 mt. Not autoclavable

Art. No. 3MEDT0006



Resuscitator balloon Master Flux Plus

for adults, autoclavable with double balloon

Art. No. 2000-R2



Circuit tube

Tube length: 1.8 mt. Autoclavable 134°C. Minimum: 30 cycles Art. No. 3MEDT0018

Circuit tube

Tube length: 50 m Single use

Art. No. 3MEDT0019



Resuscitator balloon Master Flux Plus

for children, autoclavable with double balloon

Art. No. 2100-R2



13



Service/Guarantee/Normative

Guarantee duration: 2 years legal guarantee

CE regulation and normatives: 93/42 EEC Service proposal

Service: spare parts and repairs ensured by the local authorized dealer



14

Tel: +39 0521 8380

In Italy or at the local

distributor workshop

Get more details on our service policy

sending an email to: cme2@tecnogaz.com

TECNO-GAZ

First Aid - the life-saving device

Oxygenation

VITAL TEST is a professional pulse oximeter designed for use in a medical environment. Simple and multi-functional, it can be supplied with a special support, designed to be fixed to a unit, wall or stand. This option allows for constant control of the instrument and permits it to be convenient for continuous use.

Main features

- wavelength or decimal reading of the values on the display
- adult-child reading selection
- on/off sound and adjustable volume
- SpO2 min / max alarms
- Pr min / max alarms
- patient identification
- data reading

On request Vital Test can be supplied programmed with an interface to be connected to a computer, for the registration and printing of all parameters measured.

TECTOHEARTON

Vital Test

Pulse oximeter Art. No. 11250

Defibrillation

Tecnoheart Plus enables safe defibrillation of heart muscle walls since it is equipped with sensors which are able to recognize heart failure due to arrhythmia, ventricular fibrillation and ventricular tachycardia.

- **1 Cover:** the cover is designed to protect the action icons, patient mode button and charge button.
- 2 Status indicator: indicates the status of the unit, the temperature and the battery level.
- **3** Paddle connector: the paddle connector is used to connect the electrodes.
- 4 Action lcons: the LED indicator flashes red under the corresponding action icon.

5 Patient mode switch: once the user has identified the patient based on type, select patient mode (either adult or paediatric patient) using the mode switch.

- 6 Charge button: once preparation for the electric discharge is completed, the current discharge button will flash. Press the button to deliver the discharge.
- 7 Sliding button: the sliding button opens the cover, to turn on push the sliding button to the right.
- 8 SD card slot: the SD card is required to save data and update AED firmware.
- 9 Infrared communication port: the infrared communication port is used to communicate with the PC.

Tecnoheart Plus

Semi-automatic defibrillator Art. No. 11250

100

THE UNIQUE DEFIBRILLATOR WITH EXCLUSIVE FEATURES!



3





Anticipate, advise and guide the change.

WE MANUFACTURE THE FUTURE

BUILDING THE FUTURE, DAY AFTER DAY, IDEA AFTER IDEA.







 Tecno-Gaz S.p.A.

 Strada Cavalli, 4 - 43038 - Sala Baganza - Parma - Italia

 Tel. +39 0521 83.80

 Fax +39 0521 83.33.91

 www.tecnogaz.com

 Cap. Soc. € 280.000 i.v. C.F. e P.IVA/VAT IT00570950345

 R.E.A. PR 138927 Iscr. Reg. Impr. PR 10061

www.tecnogaz.com

