# HELIX

## TEST SYSTEM

PROCESS CHALLENGE DEVICE

### ISTRUZIONI D'USO

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

PROCESS CHALLENGE DEVICE

### OPERATING INSTRUCTIONS

- 1 Open the Helix test system PCD device, separating the cap from the capsule which houses the test strips of the Helix test system PCD device. (fig.1)
- 2
  - a. Fold the strip along the score line. (fig.2a)
  - b. The indicator must face inwards. (fig.2b)
- 3 Insert the folded strip into the chamber. The folded side must be facing towards the outside of the chamber (incoming steam flow). (fig.3)
- 4 Close the Helix test system PCD device hermetically, making sure that the seal ring is positioned properly.(fig.4)
- 5 Run a sterilisation cycle. The sterilisation parameters required for the indicator to change colour are:
  - 134°C of saturated steam for 3.5 minutes
  - 121°C of saturated steam for 15 minutes
- 6 At the end of the test open the capsule of the Helix test system PCD device and remove the indicator.
- 7 If all the air present inside the Helix test system PCD device has been effectively removed, the saturated steam has penetrated inside the cannula until reaching the chamber and a temperature of 134°C has been maintained for at least 3.5 minutes (alternatively: 121°C for 15 minutes), the indicator will have turned from yellow to black to indicate that all the sterilisation parameters have been met.

The HELIX TEST SYSTEM complies with the following Standards: EN 867-5 / ISO 11140-1 Class 2

1) Sufficient removal of air, penetration of steam and temperature reached.

2) Insufficient penetration of steam and insufficient removal of air.

3) Insufficient removal of air and penetration of steam, temperature reached.

4) Insufficient removal of air and penetration of steam, temperature not reached.

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

**VALIDOLINE**

**TECNO-GAZ**  
www.tecnogaz.com

# HELIX TEST SYSTEM PROCESS CHALLENGE DEVICE

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