Circuiti e maschere per  
Circuit and mask for  
Schlauch mit maske für  
Circuitos y mascarillas para  
Circuitos e máscaras para

*Master Flux*

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1. INTENDED USE
2. CONNECTION
3. CLEANING

**1 INTENDED USE**

The Tecno-gaz circuits and masks are accessories of device Master Flux and should only be used with the device Master Flux for the supply of medical gases to the patient.



TECNO-GAZ is not responsible if the circuits and masks are used in a manner different from that stated above.

TECNO-GAZ is not liable for misuse of circuits and masks.

The masks are sold separately from the MASTER FLUX ONLY AS PARTS OF THE DEVICE.

## 2 CIRCUIT AND MASK


### Components CIGUIT

|                 |   |  |
|-----------------|---|--|
| <b>MF820ZMF</b> | CIRCUIT COMPLETE  |    |
| <b>MF835ZMF</b> | MASK size 1   |   |
| <b>MF836ZMF</b> | MASK size 2   |  |
| <b>MF837ZMF</b> | MASK size 3   |  |
| <b>MF838ZMF</b> | MASK size 4   |  |
| <b>SMFA197</b>  | CIRCUIT EXHAUST PASSIVE                                 |  |
| <b>MF824ZMF</b> | CIRCUIT CONNECTION OF<br>ASPIRATOR<br>MASTER FLUX PLUS  |  |
| <b>MF800IMF</b> | CIRCUIT CONNECTION OF<br>ASPIRATOR<br>MASTER FLUX SMART |  |


The circuit **MF824ZMF** / **MF800IMF** is composed of:

|  |   |
|--|---|
| TERMINAL ASPIRATION  |  |
| REDUCTION FOR SUCTION SALIVA   |  |
| ONE-WAY VALVE  |  |
| MASCHIO PORTAGOMMA<br>Only in MF800IMF for the MASTER FLUX SMART version |  |

the circuit **MF820ZMF** is composed of:

|   |   |
|---|---|
| Fitting connection to flowmeter box (*)<br>► Connection 22 M/F 90°  |   |
| Two tubes length 1.8 m (material PVC/PE)<br>► Transparent and flexible tube.<br>► Corrugated tube (rigid spiral) externally (anti-crushing) smooth inside for low compliance and noise.<br>► Flexible connectors for easy insertion.  |  |
| Exhaust fitting connection (material: stirolox) (*)<br>► Reduction 22-15  |  |
| Reducing exhaust connection (material: TPE, PP) (*)<br>► Elastomeric connector 9-11 mm  |  |
| (*) DO NOT USE IN THE MASTER FLUX SMART VERSION   |   |
| PIPE CLIP SUPPORT AND VELCRO BAND<br>► Allows correct positioning of the tubes, the adhesive clip can be secured behind the headboard or via the Velcro strip can be positioned on the headboard without a secure fixing.<br>► It allows a positioning with fixed or rotary movement, according to the needs. |  |

The mask is composed of:


|  |   |
|--|---|
| <p><b>Kit Silhouette nasal and Support front</b><br/> <b>Code 2ZMFA0050 size 1</b><br/> <b>Code 2ZMFA0051 size 2</b><br/> <b>Code 2ZMFA0052 size 3</b><br/> <b>Code 2ZMFA0054 size 4</b><br/> (material: silicon)<br/> ► It reduces N2O environment emissions<br/> ► Create a mask placement comfortable</p> |    |
| <p>Shell (material: PC)<br/> ► Transparent</p>   |    |
| <p>Small support front (material: PC)<br/> ► Transparent<br/> <b>For mask size 1 and size 2</b></p>  |    |
| <p>Large support front (material: PC)<br/> ► Transparent<br/> <b>For mask size 3 and size 4</b></p>  |  |
| <p>Fitting Y (PC) with valves (PC/silicone)<br/> ► Allows the correct flow of the gas mixture by preventing the exhaled gas can flow back into the supply circuit.<br/> ► The small room eliminates the risk of permanence CO2</p>   |  |
| <p><b>HEAD RIBBON WITH FASTENER FOR HEAD RIBBON</b><br/> ► The mask stabilizing and regulating via a velcro the length to adapt it to the patient.</p>   |  |

# CONNECTION INSTRUCTIONS

Check the integrity of the components.

Choose the most suitable mask for the patient.

Place the mask and adjust the head ribbon. The mask must be stable and comfortable for the patient. Connect the pipes, paying attention to the arrow on the Y fitting. The arrow indicates the flow direction. The non-return valves, placed in the fitting, push the gas flow in the direction indicated by the arrow.

**Gas inlet pipe:** Connect a hose to the special corrugated tube (Item G) to the flow control box by inserting the connector connection flow control box ".

Connect the free end to the Y connector placed on the mask, WITHOUT SIDE ARROW.

**Gas inlet pipe:** Connect the other corrugated pipe to the appropriate hose (Pos P) through the elastomeric connector. Connect the free end to the Y connector placed on the mask, SIDE WITH ARROW.

Fix the clip to stabilize the pipes behind the headboard of the chair.

And 'possible fixing with adhesive (fixed) or through the appropriate supplied Velcro (detachable).

Mounting clip  
- Fastening with velcro

Delivery knob  
Flush button  
Dosing knob

**A**  
**B**  
**C**

**G**

Gas inlet hose  
fitting

**F**

Anaesthesia bag

**P**

Gas outlet hose  
fitting

Delivery knob

**A**

**C**

Dosing knob

**B**

Flush button

**P**

Gas outlet hose fitting

**G**

Gas inlet hose fitting

**F**

Anaesthesia bag


EVACUATION EXHALED GASES

The operator can choose the best solution for processing the exhaled gases, using the following connections:

| PASSIVE EXHAUST  |  |
|--|--|
| UNIT-MOUNTED MASTER FLUX   | WALL-MOUNTED MASTER FLUX   |
| <p>Connect the hose, exhaust gas circuit mask to the hose in the center at the center of the unit.</p> <p>Connect the circuit SMA197, to the hose connection at the bottom of the cabinet.</p> <p>Connect the end of the spiral tube to the outside and make sure there is a final filter.</p> <p>The discharge must take place in an open, well-ventilated.</p> | <p>Connect the exhaust gas pipe of the appropriate circuit hose (SMFA257) after the spiral tube (SMFA197).</p> <p>Connect the end of the spiral tube to the outside and make sure there is a final filter.</p> <p>The discharge must take place in an open, well-ventilated.</p> |
|    |   |

| ACTIVE EXHAUST  |  |
|---|--|
| UNIT-MOUNTED MASTER FLUX  | WALL-MOUNTED MASTER FLUX   |
| <p>Connect the hose, exhaust gas circuit mask to the hose in the center of the cabinet.</p> <p>Connect the circuit MF824ZMF/MF800IMF to the hose on mobile.</p> <p>Connect the terminal to the suction.</p>   | <p>MF824ZMF/MF800IMF connect the circuit to the discharge circuit.</p> <p>Connect the terminal to the suction.</p> |
|   |                                 |
| <p><b>SUCTION</b><br/>MF824ZMF/MF800IMF</p>   |  |
| <p><b>SETTING THE VALVE</b></p> <p>Turn the ring perforated valve (SMFA236) until closing holes.</p> <p>Start dispensing OXYGEN setting ~ 10 NI / min.</p> <p>Start the surgical aspirator, set the minimum intake.</p> <p>Observe THE BAG ANESTHESIA.</p> <p>Turn the ring perforated valve gradually opening the holes until you get the ideal state of the BAG ANESTHESIA "HALF FULL HALF EMPTY".</p> <p><b>CAUTION:</b></p> <p><b>By varying FLOW, it may be necessary, adjust the correct intake flows of exhaled gases.</b></p> |                               |

## 4 CLEANING

**Use detergents that are not aggressive to plastic**

Clean the surfaces of the device with a dry cloth.

Before and after use, clean, clean or sterilize the parts in contact with the patient or with his breath.

### Cleaning table for MASTER FLUX parts

Attention:

Use a disinfectant that is compatible with plastics (avoid cream products).

Do not spray the disinfectant directly onto the component to be cleaned, but onto a cloth.

Use disinfectants in the dilution indicated by the manufacturer.

| ITEMS:                                       | DESCRIPTION                                      | IMAGES   |                         |            |                          |               |
|--|--|--|-------------------------|------------|--------------------------|---------------|
| MF835ZMF<br>MF836ZMF<br>MF837ZMF<br>MF838ZMF | SHELL COMPLET FITTING<br>VALVE ENF FRONT SUPPORT |  |                         |            |                          |               |
|  |  | MATERIAL   | THERMAL<br>DISINFECTION | N ° CYCLES | STEAM<br>STERILIZZAZIONE | N °<br>CYCLES |
|  |  | POLICARBONATO -<br>SILICONE  | SI                      | 10*        | 134°                     | 10*           |

|  |   |  |                         |            |                          |               |
|--|---|--|-------------------------|------------|--------------------------|---------------|
| MF835ZMF<br>MF836ZMF<br>MF837ZMF<br>MF838ZMF<br><br>MF824ZMF<br>MF800IMF | SILHOUETTE NASAL AND<br>SUPPORT FRONT<br>- RINGS -<br>EXHAUST CIRCUIT PASSIVE |  |                         |            |                          |               |
|  |   | MATERIAL   | THERMAL<br>DISINFECTION | N ° CYCLES | STEAM<br>STERILIZZAZIONE | N °<br>CYCLES |
|  |   | SILICONE - ALUMINUM<br>- PLASTIC   | SI                      | 10*        | 134°                     | 10*           |

|  |   |  |                         |            |                          |               |
|--|---|--|-------------------------|------------|--------------------------|---------------|
| MF835ZMF<br>MF836ZMF<br>MF837ZMF<br>MF838ZMF<br><br>MF820ZMF<br>3MEDN0002<br><br>SMFA197 | HEAD RIBBON WITH<br>FASTENER -<br>CIRCUIT COMPLETE WITH<br>CONNECTORS -<br>RESERVOIR BAG 3LT -<br>EXHAUST CIRCUIT PASSIVE |  |                         |            |                          |               |
|  |   | MATERIAL   | THERMAL<br>DISINFECTION | N ° CYCLES | STEAM<br>STERILIZZAZIONE | N °<br>CYCLES |
|  |   | NEOPRENE - PE - PVC<br>SBS - TPE - VELCRO  | NO                      | -          | NO                       | -             |

| OPTIONAL  |                                 | MATERIAL               | THERMAL<br>DISINFECTION | N ° CYCLES | STEAM<br>STERILIZZAZIONE | N °<br>CYCLES |
|-----------|---------------------------------|------------------------|-------------------------|------------|--------------------------|---------------|
| 3MEDN0003 | RESERVOIR BAG 3LT -             | SILICONE               | 134°                    | 10*        | 134°                     | 10*           |
| 3MEDT0018 | TUBES LENGTH 1.5 MT             | HYTREL -<br>POLIESTERE | 134°                    | 10*        | 134°                     | 10*           |
| 3MEDT0019 | TUBES LENGTH SINGLE USO<br>50MT | LDPE                   | NO                      | -          | NO                       | -             |

(\*) Maximum recommended number of cycles.

Manual cleaning and manual disinfection should only be considered where no facilities exist for automated washing and disinfection, or automated cleaning is contraindicated.

#### CLEANING – MANUAL

Hand wash using a mild detergent such as a neutral pH liquid soap or tri-sodium phosphate. Rinse thoroughly with clean water, and wipe with an appropriate soft absorbent material that will not leave fibers on the device. Hard brushes or other materials that might damage the products should not be used. Do not use a washing machine. Dry at room temperature. Reservoir Bags and Rubber Masks may be dried by hanging in a warm cabinet at a maximum temperature of 50° centigrade. Do not dry on a hot surface.

#### DISINFECTION – MANUAL

Use a proprietary disinfecting agent. Concentration and immersion times stated by the manufacturer must be strictly adhered to. Compatibility with polychloroprene, stainless steel, polysulfone, polycarbonate and polypropylene materials used in the manufacture of the product must be established prior to conducting any disinfection cycle or routine. Rinse thoroughly with water of a suitable quality, and wipe with an appropriate soft absorbent material that will not leave fibers on the device. Dry at room temperature.


#### STERILISATION

Products can be sterilized by moist heat (autoclaving) up to temperatures of 137° centigrade in accordance with methods specified in ISO17665. Do not sterilize with dry oven sterilizers.

**MASK ASSEMBLY**

After cleaning the device, assemble the mask according to the diagram below.

**ATTENTION TO THE DIRECTION OF FLOW  
ARE THESE VALVES  
BINDING THE FLOW IN THE DIRECTION SHOWN**


Circuit INLET

INLET

Circuit OUTLET


OUTLET

# Exhaled Gas Circuit identification

For easy identification of discharge circuit, during are provided template-circuit, assembly, two blue rings (Cod. 1ZMFA0082).

Figures A B C D show how to use the rings. The indicator is autoclavable.

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**TECNO-GAZ S.p.A**