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0.INDEX

PAR	THEME	PAGE
0	Index	18
1	Instructions for package and storing - disposal	18
2	Technical features	19
3	Supplied accessories	20
4	Accessories on request	20
5	Control panel	20
6	Replacement of the printer roll	21
7	Installation	22
8	Preparation for sterilisation	23
9	Sterilisation table	23
10	Filling water	24
11	Emptying the tank	24
12	Discharge the chamber	24
13	Autoclave Start-Up	25
14	Night cycle	26
15	Maintenance 14.1 Daily routine maintenance – 14.2 Weekly routine maintenance 14.3 Quaterly maintenance	27
16	Autoclave Set-Up	28
17	Trouble Shooting	29
18	Guarantee	30

TECNO-GAZ S.P.A. RESERVE THE RIGHT TO CARRY OUT IMPROVING TECHNICAL CHANGES WITHOUT ANY PREVIOUS NOTICE.

THIS MANUAL IS A SOLE PROPERTY OF TECNO-GAZ S.P.A. AND CANNOT BE COPIED, REPRODUCED OR TRANSFERRED TO THIRD PARTIES WITHOUT THE WRITTEN AUTHORISATION OF TECNO-GAZ S.P.A.

WARNING!

THE VERTICALE AUTOCLAVE IS NOT CERTIFIED FOR THE STERILIZATION OF MEDICAL DEVICES

1. ISTRUCTIONS FOR PACKAGE AND STORING

The autoclave is packaged in the following way:

It is wrapped in a POLYETHYLENE bag with corrugated cardboard internal and contained in a CORRUGATED BOARD BOX whose base has two wooden strips.

The package is not sterile.

The package and the autoclave are fragile, therefore handle with care, transport without shakes and hits and do not turn upside down. Do not lift the autoclave with sudden strokes; **the package handles are used only for the vertical lifting.**

Preserve in dry and protected areas at a temperature (+5;+40)°C.

Don't stack the packages.

PACKAGE OVERALL

(580x1100x590) mm

WEIGHT OF THE PACKAGED MACHINE: 67 Kg

DIMENSIONS:

THE PACKAGE MUST BE PRESERVED DURING THE WHOLE GUARANTEE PERIOD; TECNO-GAZ S.P.A. DO NOT ACCEPT RETURNS WITHOUT ORIGINAL PACKAGE.

TO REMOVE THE PACKAGE, PROCEED AS FOLLOW:

PUT THE DEVICE IN HORIZONTALLY POSITION, OPENS THE LOWER SIDE, BACK TO DEVICE IN VERTICAL POSITION AND SLIP OFF THE PACKAGE.

DISPOSAL

This product is subject to Directive 2002/96/EC of the European Parliament and the Council of the European Union on Waste of Electrical and Electronic Equipment (WEEE) and, in jurisdictions adopting that Directive, is marked as being put on the market after August 13, 2005, and should not be disposed of as unsorted public waste. Please utilise your local WEEE collection facilities in the disposition and otherwise observe all applicable requirements.

2. TECHNICAL FEATURES

DEVICE MECHANICAL FEATURES	
Use room temperature	+5°C ÷ +40°C
Device dimensions (width x height x depth)	500mm,1000mm,450mm
Opened door overall dimensions	1230 mm
Weight	60 Kg
Device furniture colour	RAL 9016
Material making up the device furniture	Steel Fe370
Furniture temperature after 1 working hour	(50÷55)°C
Sound intensity level	Lower than 70 dB
DEVICE ELECTRICAL FEATURES	
Supplying voltage	230 V±10% AC
Number of phases	1
Supplying frequency	50/60 Hz
Max. absorbed power	2500 W
Type of protection	PE
Insulation class	1
Length of the supplying cable L=2200mm	2x1.5mm ² + EARTH
Fuses (6.3X32)	12A delayed
Electrical connector (pin)	SCHUKO 16A+EARTH
DEVICE PNEUMATIC FEATURES	
MAX. using pressure	2.4 bar
FEATURES OF THE EMPTYING TANK	
Capacity of the filling tank	1.2 litres
Material making up the filling tank	High density polyethylene
FEATURES OF THE STERILIZATION CHAMBER	
Dimensions of the sterilisation chamber	Ø 243 x 500 mm
Material making up the sterilisation chamber	STAINLESS STEEL AISI 304

3. SUPPLIED ACCESSORIES		5. DISPLAY BOARD	
The autoclave is sold with the accessories included in the article and contained in the same package; here is a list of these accessories.		The positions indicated on the table are referred on the pictures at the page 3.	
PIECE No.	DESCRIPTION	1A	ON/OFF main switch
1	Door adjusting key	2A	DAY – MONTH – YEAR Display
4	Silicon pipe	3A	HOUR - MINUTE Display
1	Sponge	4A	TEMPERATURE Display
4. ACCESSORIES ON REQUEST		5A	PRESSURE Display
PIECE No.	DESCRIPTION	6A	TIME/MESS. Display
1	Basket with hole Stainless Steel AISI304 dim. Ø 237x430mm	7A	START/STOP Push-button
		8A	SELECT Push-button (Cycle selection)
		9A	+ Push-button (Parameters setup)
		10A	- Push-button (Parameters setup)
THE EQUIPMENT LISTED ABOVE DO NOT REQUIRE ADVICES		11A	121°C Cycle Led normal discharge 15min.
		12A	121°C Cycle Led rapid discharge 15min.
		13A	121°C Cycle Led normal discharge 20min.
		14A	121°C Cycle Led rapid discharge 20min.
		15A	134°C Cycle Led Free
		16A	Used water Led maximum level
		17A	Printer
		18A	Door opening handle
		19A	Rotating wheels
		20A	Autoclave door
21A	Chamber discharge tap		
22A	Tank discharge tap		
23A	Tank condensate discharge		
24A	Cover		

NOTE: FOR ALL REFERENCE AT THE OPERATIONS SPECIFIED ON THIS MANUAL, TO KEEP AT THE PRESENT TABLE.

6. REPLACEMENT OF THE PRINTER ROLL


To replace the paper roll, open the cover of the special compartment and insert a new thermal paper roll (maximum width = 57 mm) and let paper slip out while closing the cover back in place.

At the end of each sterilisation cycle, the autoclave prints a validation report. To cut the report off the autoclave, pull it slightly upwards against the in-built cutter. **After removing a report, make sure that the cover is properly closed.**

Thermal paper can be printed on one side only. So, always check on paper orientation.

To assure long lasting and due filing of sterilisation reports, make sure to store them in a place far from direct light and heat sources.

If the report is hardly readable, perform cleaning of printer head (Caption A on the picture). To do so, first open the paper cover and then remove the roll to perform head cleaning using a clean cloth damped in alcohol.

7. INSTALLATION

The installation is a fundamental operation for the following use and good operation of the device. Here is a list of the points to be followed in order to achieve a good installation.

1. The device must be installed inside a laboratory accessible only to authorised personnel.
2. The room in which VERTICALE AUTOCLAVE is installed must be lighted in according to the laws in force for the lighting standard in the work environment.
3. Position the autoclave at such a height as to make it possible for the operator to check the whole sterilisation chamber and carry out the normal cleaning operations.
4. The room where the device is installed must be enough ventilated.
5. Do not install the autoclave near washing basins, taps, because the device furniture is not waterproof.
6. Do not obstruct the ventilation grids of the device.
7. Do not lean on the door when it is opened.
8. Do not install the device near sources (for instance: Bunsen burners); keep a minimum distance of 40 cm.
9. Install the device in such a way as not excessively bend the supplying cable.
10. Make sure that the electrical system to which the device is connected is produced according to the laws in force and has the suitable dimension for the device features; see PAR.2 of this manual.
11. Connect the device pin to a socket with the characteristics as explained on PAR.2 of this manual or the rating date on the plate of the device.
12. Never connect the device pin to reductions of any type.
13. After the installation, carry out a LOADED sterilisation cycle as explained on PAR.12 of this manual.

ADVERTISING FOR THE PRESSURE EQUIPMENT

The pressure equipment (sterilisation chamber) is a integral part of the autoclave. Moreover the procedures for installation as Par.7, it is necessary to observe as follows:

- 1) Possible operation on the safety valve must be done only by authorised personnel or under authorisation of TECNO-GAZ S.P.A..

To keep the pressure equipment for the better use, moreover the standard maintenance operation described on Par.15, it is necessary to observe as follows:

- 1) Not damage the sterilisation chamber with any tools
- 2) Don't use any corrosive substances or liquids into the sterilisation chamber.

8. PREPARATION TO STERILISATION

In order to obtain a good sterilisation and guarantee a long life at all the autoclave components, make notice of the following points:

1. The instruments must be immersed in a disinfectant liquid or a cleaning solution for the cleaning and the disinfecting, immediately after their use. Respect the instructions of the producer for the solution quantity and exposure time.
2. Position the instruments in a ultrasound device.
3. Rinse the instruments in flowing water. Brush off eventual residual disinfectant liquid from the instruments. **WARNING:** Cleaning and disinfecting chemicals may effect and corrode parts of the autoclave.
4. Carefully dry the instruments to avoid spots formation due to the use of not desalted water .
5. Wash, rinse and dry the supplied basket as well.
6. For a better sterilisation compound instruments such as forceps, shears or alike should be kept opened.
7. Do not stack the instruments. The instruments must be sterilised separately. An overload could compromise the sterilisation.

The information at the previous point, show the importance of instruments preparation.

Also with new instruments, it is necessary to follow these instructions, because they have working residual as oil and grease.

The formation of corrosion agents on the autoclave components depend only on external agents. The introduction of just only one instrument with present trace of corrosion, can be cause of chemical and physical contamination to other instruments and to the autoclave components.

9. STERILIZATION TABLE

CYCLES	TIME OF STERILIZATION	TIPE OF DISCHARGE	WORKING PRESSURE
C1 121°C	15 minutes	NORMAL (liquids)	(1.1÷1.2) Bars
C2 121°C	15 minutes	RAPID	(1.1÷1.2) Bars
C3 121°C	20 minutes	NORMAL (liquids)	(1.1÷1.2) Bars
C4 121°C	20 minutes	RAPID	(1.1÷1.2) Bars
C5 134°C MAX. From 100°C to 140°C	From 1 to 99 minutes	RAPID OR NORMAL (solids or liquids)	(1.1÷1.2) Bar (2.1÷2.2) Bar

The data showed on the table are indicatives; the choice of the sterilisation cycle must be done in according with the specification of the products to be sterilised.

10. FILLING WATER

The water filling of the autoclave is realized from the top of the device; open the door and fill the water inside the chamber. The maximum water level is down of the support of the basket.

NOTE: USE DEMINERALIZE OR DISTILLED WATER ONLY. The use of feedwater with contaminants at levels exceeding those given in the table above can greatly shorten the working life of the device, with damage at the components of the autoclave and then invalidate the manufacturer's warranty.

TABLE OF QUALITY LEVEL IN ACCORDING TO DIRECTIVE DIN EN 285

CEN STANDARD DIN EN 285			
Evaporation residue	≤	10	mg/l
Silicium oxide, SiO ₂	≤	1	mg/l
Iron	≤	0.2	mg/l
Cadmium	≤	0.005	mg/l
Lead	≤	0.05	mg/l
Rest of heavy metals except iron, cadmium, lead	≤	0.1	mg/l
Chloride (Cl ⁻)	≤	2	mg/l
Phosphate (P ₂ O ₅)	≤	0.5	mg/l
Conductivity (at 20°C)	≤	15	µs/cm
PH value (degree of acidity)		5 to 7	
Appearance		Colourless clean without sediment	
Hardness (E Ions of alkan earth)	≤	0.02	mmol/l

11. EMPTYING THE TANK

The device is fitted with a tank for collecting the waste waters used during the sterilisation cycle.

When the discharge tank is full, the max level pilot light EXHAUST TANK, is turned on.

To discharge the tank, proceed as follows:

- Switch off the autoclave
- Remove the cover (A); this is magnetically fixed.
- Open the discharge tap.
- Flow the waters into a container.
- In case of sterilization failure, the used water in discharge tank, can have contaminated residual. This water must be disposed in the black water, use protective latex gloves to carry out the operation.

The device does not work in case the discharge tank is full. In case the operator tries to start a cycle even in this state, the device locks visualising ER3 on the TIME/MESS. display.

The tank is supplied with a breather pipe for condensed steam; position into a collection container.

12. DISCHARGE THE CHAMBER

The device is equipped with an appropriate valve for the discharge chamber from the water used during the sterilization cycle.

To discharge the chamber, follow the directions described below:

- Turn off the autoclave and wait for the chamber to cool down
- Remove the cover (24A), applied magnetically
- Unscrew the cap from the valve (21A)
- Screw the hose connection to the tap (21A)
- Insert the pipe to the hose connection
- Open the drain valve.
- Drain the water into a container
- The water used for sterilization, may contain contaminated waste, so use latex gloves to do this.

After unloading

- Close the valve (21A)
- Remove the pipe and unscrew the hose connection
- Screw the cap on the valve (21A)
- Close the cover (24A)

13. AUTOCLAVE START-UP

WARNING BEFORE USE

A SLAGGING PERIOD OF THE AUTOCLAVE WITH THE DOOR CLOSED, CAN BE DEFORM THE GASKET SEAL AND TO CREATE LACK OF PRESSURE DURING THE FIRST USE. WE ADVISE TO CHECK THE LOCKED DOOR. IF IS TOO SLOW, ADJUSTING IT AS DESCRIBED ON PAR.14.

1. Carry out the installation operations as described in PAR.7, the autoclave is ready to be used. Pay attention to the instructions below:
2. Turn on the autoclave pressing the ON/OFF main switch.
3. Open the door using the lever.
4. Check the level of the distilled water, as described in PAR.10.
5. Close the door using the handle.
WARNING: WHEN CLOSE THE DOOR, TO BE CARE AT THE HANDLE, TO AVOID THE SQUASHING RISK (see at the label positioning on the autoclave – photo at pag.2 and pag.4).
6. Select the cycle using the button SELECT.
7. On the TIME/MESS. display, the initials relating to the set cycle **C1** (if the cycle 121°C 15 min. normal discharge), **C2** (if the cycle 121°C 15 min. rapid discharge), **C3** (if the cycle 121°C 20 min. normal discharge), **C4** (if the cycle 121°C 20 min. rapid discharge) or **C5** (if the cycle 134°C free).
8. To set the parameter of the **C5** cycle, proceed as follow: press and keeping pressing the START/STOP button for 10 seconds; on the TIME/MESS. display the selection page –5- is showed. On the TEMPERATURE display is showed the the cycle time in minutes. Press the button + to increase the value or press the button – to decrease. Press the START/STOP button to pass at the selection page – 6-; On the TEMPERATURE display is showed the the cycle temperature. Press the + button to increase the value or press the – button to decrease. Press the START/STOP button and keeping pressing for 10 seconds, to store the programming; on the TIME/MESS. display the digit **C5** that confirm the data storing. The values of pressure and temperature are automatically calibrated, in according to thermodynamic physical laws. To select the type of discharge at the end of the cycle (normal or rapid), press the – button; the cycle led has a fixed light with the normal discharge and the blinking light with the rapid discharge. This type of selection is possible until to reaching of the sterilisation phase. Now is possible to start the cycle.
9. Press START/STOP button. The heating band element starts to warm, until reached the cycle values. To eliminate the air inside the chamber to obtain a good sterilisation, during the first warm up phase, a solenoid valve is open to permit the outlet.
10. At the end of the sterilisation cycle, the solenoid valve is opens to permit the discharge pressure. In the rapid discharge, the solenoid valve opens immediately after the end of cycle; in the normal discharge, the solenoid valve opens when the internal temperature of the chamber decrease down the 97°C.
11. If at the end of sterilisation cycle, the chamber remain in pressure or in depressure (is not possible open the door), press the button + that open the solenoid valve, to permit at the internal pressure the alignment with the atmospheric pressure.
12. In case after starting the sterilisation cycle you wish to stop it, proceed as follows: Press the START/STOP button, wait that the pressure indicated by the PRESSURE display, is at 0 bar. On the TIME/MESS. display, the writing "Int" (interrupted cycle) is displayed, press the SELECT button, to cancel the alarm. Open the door with the handle.
13. If at the end of the sterilisation cycle on the TIME/MESS. display, an alarm code is displayed, for instance: (AL6, AL7, etc.... see PAR. 16), this means that the sterilisation cycle has not taken place in the right way, so it is necessary to repeat the cycle.
14. At the end of the cycle on the TIME/MESS. display the writing END is displayed followed by an acoustic sign: Now it is possible to open the door by following the warnings below.
15. Before opening the door, make sure that there is no pressure inside the sterilisation chamber; the "PRESSURE" display (5A), must indicate 0 bar. IF THE PRESSURE VALUE IS HIGHER THAN 0, IS NOT POSSIBLE OPENING THE DOOR.
16. USE PROTECTIVE GLOVES TO REMOVE THE LOAD FROM THE STERILIZATION CHAMBER.
17. WHEN OPENING THE DOOR, DO NOT STAY ABOVE OR IN FRONT OF IT TO AVOID BURN DUE TO THE VAPOUR COMING OUTS.
18. BEFORE YOU REMOVE THE COVER (24A), TURN OFF THE AUTOCLAVE

14. NIGHT CYCLE

VERTICALE AUTOCLAVE is programmed to interrupt heating after 60 inactivity minutes (STAND BY). Then, is possible to run sterilisation cycles without operator, since when the cycle end, after 60 minutes, the autoclave is automatically turned off.

All the LEDS are turned off and only the ON/OFF main switch, remains on.

To check that sterilisation has taken place, the operator must press START/STOP button or SELECT button, to reset the autoclave and to verify the result of the sterilisation..

If the TIME/MESS. display indicates END, this means that the sterilisation has taken place.

If during the sterilisation cycle some anomalies are revealed, an alarm or error message is displayed on the TIME/MESS. Display and by the printer report.

On the PAR.16, are explained some possible anomalies during the sterilisation cycle.

15. MAINTENANCE

CAUTION : DISCONNECT ALWAYS THE MAIN SUPPLY BEFORE EACH TECHNICAL INTERVENTION

14.1. DAILY ROUTINE MAINTENANCE

CLEANING THE DOOR SEAL

Clean the seal on the door and on the outer edge of the sterilisation chamber sealing the gasket using a wet cloth or the soft part of the supplied sponge.

This cleaning operation must be carried out to remove possible dirt which can cause a pressure loss in the sterilisation chamber and a possible seal cut.

Check the water level in the tank(see PAR.10)

14.2. WEEKLY ROUTINE MAINTENANCE

CLEANING THE STERILIZATION CHAMBER

To replace the water into the chamber. Removing the panel (A), open the tap and discharge the water into a container.

To remove the deposit on the chamber floor, clean with the abrasive part of the supplied sponge.


To moisten the sponge, only use distilled or demineralized water.

14.3. QUARTERLY MAINTENANCE

- Lubricate the hinges and the closing pin with silicon oil.

- Adjusting the door:

If the autoclave handle is too loose while closing, adjust as shown in the following pictures. Rotate in clockwise direction (fig. 11) to increase the seal pressure on the chamber. In case it has loosened too much and you are not able to close the handle, rotate in counterclockwise direction.

16. AUTOCLAVE SET-UP

CAREFULLY READ THIS PARAGRAPH BEFORE PROGRAMMING. ALWAYS AND STRICTLY PERFORM PROGRAMMING WITH AUTOCLAVE DOOR WIDE OPENED.

To set the parameters which are visualised on the display, proceed as follows:

Turn on the autoclave with the ON/OFF main switch, placed on the the control panel keeping pressed the START/STOP push-button; after some seconds on the TIME/MESS. display, the programming page –1- (HOUR-DATE) is visualised; now you can release the START/STOP push-button.

On the TEMPERATURE display the blinking seconds will appear.

With the + and - buttons is possible to modify the parameters.

Press the SELECT push-button to pass at minute selection. On the HOUR/MINUTE display the blinking digit indicate the data to modify.

With the + and - buttons is possible to modify the parameters.

Press the SELECT push-button to pass at hour selection. On the HOUR/MINUTE display the blinking digit indicate the data to modify.

With the + and - buttons is possible to modify the parameters.

Press the SELECT push-button to pass at day selection. On the DAY-MONTH-YEAR display the blinking digit indicate the data to modify.

With the + and - buttons is possible to modify the parameters.

Press the SELECT push-button to pass at months selection. On the DAY-MONTH-YEAR display the blinking digit indicate the data to modify.

With the + and - buttons is possible to modify the parameters.

Press the SELECT push-button to pass at year selection. On the DAY-MONTH-YEAR display the blinking digit indicate the data to modify.

With the + and - buttons is possible to modify the parameters.

Press the START/STOP push-button in order to pass at the second programming page –2- (atmospheric pressure alignment). The display will show 0.00; keep the door open until pressure in the autoclave adjusts to the atmospheric pressure which varies according to the altitude. **CARRY OUT THE ATMOSPHERIC PRESSURE ALIGNMENT ON THE FIRST INSTALLATION OF THE AUTOCLAVE.**

Push the - button to confirm the input data. An acoustic signal confirm the input data.

Press the START/STOP push-button in order to pass at the third programming page –3- (peripheric selection).

On the TEMPERATURE display the blinking digit means the selection.

With the + and - buttons is possible to modify the parameters. 0=None – 1=Internal Printer – 2= External Printer – 3=Computer – 4=Internal Printer+Computer.

Press the START/STOP push-button in order to pass at the four programming page (autoclave serial number). This parameter is setted by TECNO-GAZ S.P.A.. The modify of this parameter can be done only by authorisation of TECNO-GAZ S.P.A. to store the settings.

Press and keeping pressed the START/STOP button until you hear an acoustic signal confirm that all data are successfully stored and to exit from programming.

On the TIME/MESS. display is visualised the last selected sterilisation cycle.

The autoclave is ready to start.

17. TROUBLE SHOOTING

The TIME/MESS. display indicates possible failures and errors which can occur during the autoclave operation.

MESSAGE	FAILURE DESCRIPTION	WHAT TO DO
Int	Sterilisation cycle discontinued. It happens when a sterilisation cycle is stopped by pressing on the START/STOP button.	Press the SELECT button to reset the alarm.
ER3	Water level error. The used water tank is full.	Empty the tank by means of the outlet tap located on the front side of the autoclave.
ER4	Open door signal upon cycle start.	Check for the door to be properly shut and repeat closing.
AL4	Open door signal during cycle.	Press SELECT to reset the alarm. Make sure that the door provides for a tight sealing against the sterilisation chamber. If required, adjust the door according to PAR.14.3 Check the door safety switch. If the display keeps showing alarm message, call for technical service.
AL5	Pressure fall. It takes place if during the sterilisation, the pressure on the chamber falls sensibly.	Press SELECT to reset the alarm. If the display keeps showing alarm message, call for technical service.
AL6	TIME OUT, waiting for pre-heating. The autoclave is not capable of reaching the pre-heating temperature.	Press SELECT to reset the alarm. If the display keeps showing alarm message, call for technical service.
AL8	ALARM FOR LACKED PRESSURE EXHAUSTION OF THE STERILIZATION CHAMBER.	Press SELECT to reset the alarm. If the display keeps showing alarm message, call for technical service.
AL9	General alarm. Cause: 1. Power supply lacking. 2. Voltage reduction higher than 10%.	Press SELECT to reset the alarm. If the display keeps showing alarm message, call for technical service.

The errors (E....) are indicated on the display for 4 seconds and they are accompanied by a sound signal. They are often caused by the user lack of attention.

The alarms (A....) are accompanied by sound signals. They remain visualised on the displayed and they are stored by the microprocessor up to when the SELECT button is pressed.

TECNO-GAZ S.P.A. REFUSES ALL RESPONSIBILITY FOR TECHNICAL REPAIRS CARRIED OUT BY NON AUTHORIZED PERSONNEL.

18. GUARANTEE

1. The product is guaranteed for a period of 12 months.
2. The guarantee begins from the date of delivery of the machine to client proven by the restitution of the guarantee label attached to the tool, filled in, stamped and signed by the dealer.
3. Repair or substitution during guarantee of a part is effected upon the absolute judgment of TECNO-GAZ S.P.A. company. It does not include the transfer of employees, packaging and shipping costs, which are charged to the customer.
4. Excluded from the guarantee are light bulbs, fuses, and components subject to normal wear.
5. The right to complete substitution of the machine is not recognized
6. The guarantee does not cover damage caused by improper maintenance, negligence, incompetence or other causes not related to the constructor.
7. The guarantee does not include any compensation for direct or indirect damage of any nature towards persons or things, caused by the inefficiency of the machine.
8. TECNO-GAZ S.P.A. is not responsible for breakdowns or damage caused by inappropriate use of the product or in case ordinary maintenance has not been carried out or elementary principles of correct maintenance (negligence) have been ignored.
9. There is no acknowledgment of compensation for machine stoppage.
10. **Mounting, service, overhaul, and control of tools have to be carried out exclusively by service centers authorized by TECNO-GAZ S.P.A. to be identified by the card issued after the participation in the certification course.**
11. **Any intervention on the tool carried out by unauthorized staff excludes the right to guarantee as well as the direct liability of the producer. This also applies to interventions performed out of the guarantee period.**
12. Components replaced under guarantee must be returned to TECNO-GAZ S.P.A. postage paid. Failure to return the part within 10 days the cost of the part will be charged to the buyer
13. TECNO-GAZ S.P.A. does not accept return from final users, unless specifically authorized in writing.
14. Return to TECNO-GAZ S.P.A. for repairs, must be controlled by the dealer or by service centers chosen by the final user in compliance with commercial procedures CM-P-003 entitled "CONTROL OF REQUESTS FOR RETURN, REPAIR, SUBSTITUTION OF TECNO-GAZ S.P.A. PRODUCTS AND CLIENT COMPLAINTS".
15. Return to TECNO-GAZ S.P.A. must be documented and authorized as per internal procedures.
16. **Products returned to TECNO-GAZ S.P.A. must have attached documentation of authorized return and a document describing the breakdown, otherwise the diagnostic cost will be charged.**
17. All products in repair must be shipped to TECNO-GAZ S.P.A. postage paid and adequately packaged (we suggest using the original packaging)
The goods travel at sender's risk.

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GB

This manual must always be kept with the product, in complying with the Directives of European Community. TECNO-GAZ reserves the right to modify the enclosed document without notice. TECNO-GAZ reserves the property of the document and forbids others to use it or spread it without its approval.

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TECNO-GAZ S.p.A.