Master Flux

Plus

LABEL:

<i>Istruzioni per l'uso</i>	ITALIANO	
<i>Instructions for use</i>	ENGLISH	<i>Translation</i>
<i>Mode d'emploi</i>	FRANCAIS	<i>Traduction</i>
<i>Bedienungsanleitung</i>	DEUTSCH	<i>Übersetzung</i>
<i>Instrucciones de uso</i>	ESPAÑOL	<i>Traducción</i>
<i>Instruções de uso</i>	PORTUGUÊS	<i>Tradução</i>

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TABLE OF CONTENTS

- 1. UTILISATION AND INTENDED USE**
- 2. SAFETY**
 - 2.1 SAFETY MARKING
 - 2.2 SAFETY DEVICES
 - 2.3 HANDLING
 - 2.3 DISPOSAL
- 3. TECHNICAL DATA**
- 4. UNPACKING**
- 5. COMPONENTS**
 - 5.1 COMPONENTS SUPPLIED LIST
 - 5.2 UNIT-MOUNTED COMPONENTS
 - 5.3 WALL-MOUNTED COMPONENTS
 - 5.4 USER INTERFACE
 - 5.5 MASKS
 - 5.6 MASKS MEASURES
 - 5.7 MASKS CIRCUIT
- 6. INSTALLATION**
 - 6.1 INSTALLATION OF UNIT-MOUNTED MASTER FLUX PLUS
 - 6.2 INSTALLATION OF WALL-MOUNTED MASTER FLUX PLUS
 - 6.3 CONNECTING CIRCUITS AND ACCESSORIES
 - 6.4 EVACUATION OF EXHALED GASES
 - 6.4.1 DIRECT CONNECTION TO THE OUTSIDE OF THE OFFICE
 - 6.4.2 CONNECTION TO A SURGICAL ASPIRATOR
- 7. OPERATION OF THE APPLIANCE**
 - 7.1 USER INTERFACE COMMANDS AND DEVICES
 - 7.2 INSTRUCTIONS FOR USE
 - 7.3 RESUSCITATION WITH MASTER FLUX PLUS
- 8. MAINTENANCE**
 - 8.1 CYLINDERS CHANGING
 - 8.2 CLEANING
 - 8.3 SCHEDULED MAINTENANCE
- 9. ALARMS**
- 10. PROCEDURES FOR AFTER SALE SERVICE**
- 11. SUPPLIED ACCESSORIES**

1

UTILISATION AND INTENDED USE**INTENDED USE**

MASTER FLUX PLUS is a device designed for conscious sedation or sedative analgesia. Conscious sedation consists in administering a mixture of nitrous oxide and oxygen to the patient.

USE

The concentrations of the administered mixture are:

O₂ ≥ 30%

N₂O ≤ 70%

It is administered through a nasal mask. The patient must inhales the mixture by breathing actively through his/her nose.

Conscious sedation is suitable for all patients of an orthodontic practice, both adults and children.

Conscious sedation allows the patients to undergo orthodontic procedures limiting the fear and anxiety that all patients experience when they enter an orthodontic office.

**IT IS NOT SUITABLE:**

- DURING PREGNANCY,
- FOR PATIENTS WITH DRUG ADDICTIONS
- FOR PATIENTS SERIOUS PULMONARY INFECTIONS
- INDIVIDUALS WITH SERIOUS MENTAL ILLNESS.

MASTER FLUX PLUS must only be administered by specifically trained medical staff.

Do not use this device to administer general anaesthesia or as a part of or in combination with a general anaesthesia system.

The sale of this device is limited by law to doctors and orthodontists only.

ATTENTION

Do not attempt to repair, modify or calibrate this device. Unauthorised repairs, modification or abuse of this device negatively affects its performance and invalidates the warranty.

Before each use make sure the gas supplies are connected correctly

The device administers drugs solely for conscious sedation

The contra-indications for the use of the drug are stated in the safety sheet and documents attached to the drug supply. If they are not provided in attachment, request it from the gas supplier.

Do not remove the label

Do not replace the supplied gearboxes with other ones not provided by the manufacture.

2

SAFETY

**Use of the PULSE OXIMETER must be routine during conscious sedation.
The meaning of the detected values:**

Saturation values O ₂	Hypoxia level
≥ 95%	No hypoxia
90-94%	Low-grade hypoxia
85-89%	Moderate hypoxia
84%	Serious hypoxia

2.1 Safety marking

The device complies with Directive 93/42 EEC as amended.
The conformity declaration is provided in attachment to the device.

2.2 Safety devices

The device is fitted with safety devices:

- If the oxygen (O₂) pressure drops, the delivery of nitrous oxide (N₂O) is automatically shut off.
- An incorporated mixing system regulates the N₂O so that there is always a minimum quantity of O₂ in the mixture of at least 30%.
- If the gas mixture flow is shut off, the patient will breathe air from the outside through a specific safety valve.
- If the flow changes direction, returning to the device, a one-way valve protects the patient from re-inhaling exhaled gases and from accumulations of carbon dioxide (CO₂).
- The flow setting knob controls the overall flow; it automatically regulates the flow of oxygen and nitrous oxide according to the set percentage.
- The 'FLUSH' button allows high flows of pure oxygen to be delivered.
- The maximum delivery of nitrous oxide (N₂O) is 70% of the total flow. The amount of oxygen (O₂) can never drop below 30%.

WARNINGS & USEFUL NOTIONS

DURING GAS ADMINISTRATION IT IS ESSENTIAL TO BE PRESENT AT ALL TIMES, TO KEEP THE OXYGEN AND NITROUS OXIDE GAS FLOWS UNDER CONTROL BY CHECKING THE FLOW VALUES DISPLAYED IN THE FLOW METERS AND THE FLOW CONTROL BOX.

2.3 Handling

Please only use the supplied handles to move the Master Flux Plus furniture unit.


Figure 1

	<p>Do not push the furniture unit without using the appropriate handles.</p>
	<p>Do not sit on the furniture item or place any weights on it.</p>
	<p>Before moving the furniture unit please carefully put the panel and tubes back into the fixing system, trying not to stumble or break anything.</p>

2.4 Disposal

The disposal of the packaging, device and single components should take place in compliance with the regulations applicable to the treatment and elimination of potentially bio-contaminating waste in the country of use.

Do not dispose of the packaging materials as household waste.

The device does not contain electrical / electronic components and contains no battery.

3

TECHNICAL DATA

POWER SUPPLY

Power supply pressure	3.5 bar \cong 50 PSI (max 6 bar \cong 87 PSI)
Maximum flow rate	10 NL/min. (litres per minute)

WEIGHT UNIT-MOUNTED MASTER FLUX PLUS:

Weight without cylinders	42 Kg
Sustainable maximum weight, tanks compartment, wheels level	180 Kg
Sustainable maximum weight, upper tanks compartment	45 Kg

CYLINDERS FOR UNIT-MOUNTED MASTER FLUX PLUS:

Maximum height	950 mm
Maximum diameter	140 mm


WEIGHT WALL-MOUNTED MASTER FLUX PLUS:

Weight of the flowmeter box	6,3Kg
-----------------------------	-------

FLOW METER:

Repeatability	+/-10% flow indicated
---------------	-----------------------

Flow rate is then measured against the flat top edge of the float, reducing parallax errors.

CONSUMPTION:

Below is an example of consumption, to be used as a theoretical reference

10 Lt MEDICAL OXYGEN CYLINDER
CHARGED AT A PRESSURE OF 200 Bar

AVAILABLE VOLUME	2000 gaseous Litres
------------------	---------------------

10 Lt NITROUS OXIDE CYLINDER
CHARGED WITH 7 Kg of LIQUID GAS

AVAILABLE VOLUME	4600 gaseous Litres
------------------	---------------------

BY SETTING:

DELIVERY	10 NLt/min
MIXTURE	50 %

AUTONOMY:

OXYGEN	400 min
NITROUS OXIDE	950 min

Taking into account the average gas consumption parameters, TECNO-GAZ advises replacing a nitrous oxide tank every 2 oxygen tanks.


Figure A


Figure B

04

UNPACKING

IT IS FORBIDDEN TO STACK MASTER FLUX PLUS UNITS

The packing shall not undergo any shock and shall be handled with care avoiding letting it roll or fall.

Do not use appliances with clear damages due to the transport.

Remove the packaging and check the appliance condition.

The packaging shall be stored for the whole warranty period.
The manufacturer does not accept returns deprived of their original packing.

The devices must be kept safe in a dry location and at a temperature between +5 / +40 °C.

05

DESCRIPTION OF COMPONENTS

5.1

Components supplied list

ARTICLE	16010ZMF wall-mounted Master Flux Plus with Reducers with AFNOR Oxygen and UNI Nitrous Oxide	1601ASD wall-mounted Master Flux Plus with Reducers with Oxygen and Nitrous Oxide DIN attachments	160IAS wall-mounted Master Flux Plus with Reducers with Oxygen and Nitrous Oxide UNI assemblies	MF009ZMF Master Flux Plus unit-mounted without gear reducer kit	MF006ZMF Master Flux Plus unit-mounted with Reducers with AFNOR Oxygen and UNI Nitrous Oxide	1600ASD Master Flux Plus unit-mounted with Reducers with Oxygen and Nitrous Oxide DIN	1600ASF Master Flux Plus unit-mounted with Reducers with Oxygen and Nitrous Oxide AFNOR	1600AS Master Flux Plus unit-mounted with Reducers with Oxygen and Nitrous Oxide UNI
UNIT-MOUNTED	SMFA301	SMFA301	SMFA301	2ZMFA0020	2ZMFA0020	2ZMFA0020	2ZMFA0020	2ZFMA0020
MASTER KIT	22ZMFA0065	22ZMFA0063	22ZMFA0062	2ZMFA0057	2ZMFA0059	22ZMFA0058	22ZMFA0061	22MFA0060
REDUCER N ₂ O	SMFA360	SMFA360	SMFA360		MF011ZMF	MF011ZMF	MK005ZMF	MF011ZMF
REDUCER O ₂	SMFA361	SMFA358	MF003ZMF		1513A4F	1513A4D	1513A4F	1513A4
MASKS								

MF820ZMF Master Flux circuit

MF824ZMF Master Flux Plus suction circuit (new)

MF835ZMF Intelliflux 1 complete mask Master Flux plus 0-4 years

MF836ZMF Intelliflux 2 complete mask Master Flux plus 4-10 years

MF837ZMF Intelliflux 3 complete mask Master Flux plus > 10 years

5.2

Unit-mounted COMPONENTS


Figure 2

FIGURE	POSITION	Q.TY	CODE	DESCRIPTION
3	1	1	SMFA144	PLASTIC COVERING
	2	1	SMFA506	FLOW MEASUREMENT BOX
	3	3	SMFA223	SPACER
	4	3	CM25021	M5 FLAP SCREW
	5	3	3MECQ0010	SPACER
	6	3	CM20015	WASHER
	7	3	CM25019	M6 FLAP NUT
	9	2	SMFA237	GAS ELIMINATION HOSE FITTING

5.3

Wall-mounted COMPONENTS


Figure 3-4

WALL MOUNTED (Figures 3-4)

FIGURE	POSITION	Q.TY	CODE	DESCRIPTION
4	10	1	SMFA313	PIN
	11	1	SMFA364	O ₂ TUBE
	12	1	SMFA365	N ₂ O TUBE
	13	1	SMFA301	WALL-MOUNTED SUPPORT
	14	1	SMFA501	WALL PLATE
	15	1	CM30011	PIN SECURING KNOB

5.4

USER INTERFACE


Figure 5

A	FLOW knob
B	FLUSH button
C	% N ₂ O knob
D	N ₂ O flowmeter
E	O ₂ flowmeter
F	Rubber bag
G	Mask connection nozzle

5.5

MASKS


Figure 6

①	Support area for the patient's forehead
②	Gas inlet with non-return valve
③	Gas outlet with non-return valve and PC printed arrow to indicate the flow direction
④	Velcro elastic band with plastic hooks

5.6

MASKS MEASURES

MF835ZMF MASK 1 0-4 years	MF836ZMF MASK 2 4-10 years	MF837ZMF MASK 3 >10 years	MF838ZMF MASK 4

Figura 7

5.7

MASKS CIRCUIT

Figure 8

COMPLETE CIRCUIT OF MASKS


Figure 9

For parts of individual components, refer to the "Circuits and Masks" user manual.

06

INSTALLATION

Remove the packaging and check the conditions of the equipment.
Do not use equipment that presents evident damage due to transportation.

6.1 Installation of unit-mounted MASTER FLUX PLUS**FLOWMETER BOX INSTALLATION** (Figure 2)

see enclosed Quick Guide 0ZMFI0007

Screw the 3 steel spacers (3) onto the flowmeter box (2).

Remove the protections from the flowmeter box supply nipples.

Install the flowmeter box on the unit by inserting the steel spacers (3) into the relative holes on the unit, securing them to the unit using 3 plastic spacers (5), 3 washers (6) and the three thumb nuts (7) supplied.

FRONT COVER INSTALLATION (Figura 2 Pos 1)

see enclosed Quick Guide 0ZMFI0007

Position the plastic cover near the three holes on the unit.

Attach the plastic cover by screwing on the three thumb screws (4) supplied.

EXHAUST NOZZLE INSTALLATION (Figura 2 Pos 9)

see enclosed Quick Guide 0ZMFI0007

If you are using the exhaled gas exhaust tube, screw the supplied nozzles onto the plastic fittings.

CYLINDER INSTALLATION (NOT SUPPLIED) (Figure 10)

Attach the reducers to the cylinders (see instruction manual for reducers).

For 5 lit cylinders install the shelf and attach it using the relative shelf stops.

For 10 and 14 lit cylinders position the shelf on the bottom of the unit.

The reducer pressure gauges (M) must be placed in such a way that they are visible to the operator.

Lock the cylinders to the piece of furniture using the special Velcro strap supplied.

Secure the cylinders to the unit using the velcro strap supplied.

Screw the O₂ reducer (R1) tube end nut onto the flowmeter box supply nipple (I) (RIGHT THREADING).

Screw the N₂O reducer (R2) tube end nut onto the flowmeter box supply nipple (H) (LEFT THREADING).

Close the rear door with the two magnets located on the top of it.


Figura 10

FIGURE	POSITION	Q.TY	CODE	DESCRIPTION
6	R1	1	1513A4	O ₂ GEARBOX (ITA)
	R2	1	MF011ZMF	N ₂ O GEARBOX (ITA)
	L	1		VALVE
	M	1		PRESSURE GAUGE

6.2 Installation of wall-mounted MASTER FLUX PLUS


Figure 11

Figure 12

(FIGURE 5, 10)

Attach the MASTER FLUX PLUS SUPPORT to the wall (13) using 4 Ø14 mm M8 wall dowels.

Attach the MASTER FLUX PLUS PLATE to the wall (14) using 4 Ø 8(mm) wall dowels.

Screw the pin (10) onto the base of the flowmeter box.

Slot the pin into the end of the support (15).

Lock the pin to the support using the relative threaded knob.

Screw the nut on the N2O tube onto the N2O flowmeter box (H)supply nipple (RIGHT threading).

Connect the nut on the O2 tube to the O2 flowmeter box supply nipple (LEFT threading).

Connect the free ends of the tubes to the wall outlets (press the couplings and turn to lock the coupling to the outlet).

NOTES ON THE CENTRALISED SYSTEM

Installation must be carried out by specialised technicians. The centralised system must be set up in compliance with regulations in force European Directive 93/42/EC

Connect the wall outlet to the equipotential point on the electrical system, with the supplied loop connectors.

The supplied reducers must be connected to the cylinders of the centralised system.

The outlet of the **N₂O** reducer includes a flexible tube that is connected on one end to a male fitting with **1/4" G** threading for connection to the system.

The outlet of the **O₂** reducer includes a flexible tube that is connected on one end to a male fitting with **1/8" G** threading for connection to the system.

The fastening plate is not supplied.

6.3 Connecting circuits and accessories:


Figure 13

FIGURE	POSITION	Q.TY	CODE	DESCRIPTION
8	A	1		'FLOW' KNOB
	B	1		'FLUSH' BUTTON
	C	1		'% N2O' KNOB
	D	1		N2O FLOW RATE GAUGE
	E	1		O2 FLOW RATE GAUGE
	F	1	3MEDN0002	RESERVOIR BAG
	G	1		MASK CONNECTION NOZZLE
	H	1		NIPPLE N2O
	I	1		NIPPLE O2

Connect the Rubber bag (F) to the flowmeter box by pressing it onto the relative nozzle.

The mask kits are connected to the device by the nozzle (G).

The gas exhaust tube will be connected in different ways, based on the type of gas evacuation system that is being used. (see 6.4 - EVACUATION OF EXHALED GASES).

6.4 Evacuation of exhaled gases:

To safeguard the work environment and the staff, the exhaled gases must be conveyed outside of the work environment.

Two systems can be used to discharge the gases externally:

- Direct connection to the outside of the office
- Connection to a surgical aspirator.

It is important to evaluate:

- Whether the aspirator is suitable for the evacuation of the gases (**N₂O and O₂**)
- The aspirator exhaust must be outdoors and ventilated.

Make sure it is compliant with regulations in force in the country of use.

6.4.1 Direct connection to the outside of the office

UNIT-MOUNTED MASTER FLUX PLUS

Screw the hoses (**Figure 2 Pos 8/9**) to the unit.

Remove the supplied valve.

Connect the mask circuit gas exhaust tube to the nozzle (**Figure 15**).

Connect the gas evacuation spiral tube (**Figure 14**) to the nozzle on the unit.

Connect the end of the spiral tube to the outside and make sure there is a final filter (Figure 14 Pos 34).

WALL-MOUNTED MASTER FLUX PLUS

Connect the gas exhaust tube to the mask circuit to the nozzle (**Figure 14 Pos 36**).

Connect the end of the spiral tube to the outside and make sure there is a final filter (**Figure 14 Pos 34**).


Figure 14 (Art SMFA197)

FIGURE	POSITION	Q.TY	CODE	DESCRIPTION
9	34	1	SMFA197	FILTER - SMFA260
	35	1		TUBE – SMFA199
	36	1		FITTING – SMFA257

Figure 15

6.4.2 CONNECTION TO A SURGICAL ASPIRATOR

UNIT-MOUNTED MASTER FLUX PLUS

Screw the hoses (Figure 2 Pos 8/9) to the unit-mounted.

Connect the mask circuit gas exhaust tube to the supplied fitting (**Figura 17**).

Connect the suction circuit (Figure 16) to the hose fitting.

Connect the fitting to the surgical aspirator tube (Ø11 mm).

WALL-MOUNTED MASTER FLUX PLUS

Remove the tube (Figure 16 Pos 43).

Connect the tube, gas exhaust of the mask circuit, to the valve (Figure 16 Pos 44).

Connect the fitting to the surgical aspirator tube (Ø11).

VALVE ADJUSTMENT

Turn the perforated ring nut on the valve (**Figura 16 Pos 44**) until the holes close.

Start gas delivery by setting ~10 l/min.

Start the surgical aspirator by setting the minimum aspiration value.

Observe the bar (**Figura 13 Pos F**).

Turn the perforated ring nut on the valve gradually opening the holes until the bag reaches its ideal condition of **“HALF FULL HALF EMPTY”**

ATTENTION

By changing the flow delivery it may be necessary to adjust the correct aspiration flows of the exhaled gases.


Figure 16 (Art MF824ZMF)

FIGURE	POSITION	Q.TY	DESCRIPTION
9A	41	1	REDUCTION
	42	1	TERMINAL PART
	43	2	PIPES
	44	1	VALVE

Figure 17

7

OPERATION OF THE APPLIANCE

7.1 USER INTERFACE COMMANDS AND DEVICES


Figura 18

CONTROL KEYS (FIGURE 10)

A	FLOW KNOB – DELIVERED FLOW REGULATION
B	EMERGENCY FLUSH BUTTON – ONLY OXYGEN DELIVERY
C	NITROUS OXIDE % REGULATION KNOB

CONTROL ELEMENTS

E	PURE OXYGEN FLOW
D	PURE NITROUS OXIDE FLOW
F	BALLOON (CHECK – ADEQUATE FLOW FOR THE PATIENT)

7.1 Instructions for use

Before using the unit-mounted MASTER FLUX PLUS version, lock the wheels that have brakes on them by pressing the tabs downwards.

If it becomes necessary to move the appliance unit (e.g. from one room to another), please unblock the wheels and push the device gently using the handles (see Paragraph 2.3)

	Avoid shocks and sudden movements that could cause accidental falls.
	Move and position the MASTER FLUX PLUS only on non-inclined and smooth surfaces.
	Make sure the cylinders are stable and secured correctly.
	Close the cylinders before moving the device.
	Do not use the handles to lift the device.
	Do not move the device using the handles.
	Do not expose the controls to involuntary knocks.

When using MASTER FLUX PLUS in its wall-mounted version, it is possible to move and position the flowmeter box as required by gripping the relative handle (Fig 4 Pos. 10).

	DO NOT LEAN ON AND DO NOT REST OBJECTS ON THE FLOWMETER BOX (BOTH WALL AND UNIT-MOUNTED VERSIONS). DO NOT FORCE THE ARM OF THE WALL-MOUNTED VERSION.
	SLOWLY OPEN THE VALVES ON THE O₂ AND N₂O CYLINDERS TURN THE VALVES COUNTER-CLOCKWISE (FIGURA 10 POS L).

Read the pressure on the pressure reducer pressure gauges, (see TECHNICAL DATA) to estimate the autonomy of the oxygen cylinder.

Please make sure that the balloon and tubes are inserted well into their fittings.

Choose the most suitable mask for the patient.

Place the mask on the patient fitting it on well around his/her nose.

Secure the mask to the loop on the back of the head rest of the dental chair.

Make sure that the tubes for the mask circuit do not get bent or choked.

The operator must be stationed in such a way that he/she can visually control:

- **THE RUBBER BAG** (Fig. 18 Pos. F)
- **THE FLOW METERS** (Fig. 18 Pos. D, E)

	During the procedure constantly check the flows displayed on the flow meters (Figure 13 Pos D, E).
	Before leaving the patient, even briefly, shut off sedation.

Make sure that the %N₂O knob (Fig.18 Pos C) index displays 0.

Using the FLOW knob (Fig.18 Pos A) administer a flow of 3 Lt/min of oxygen for 30 seconds.

INSTRUCT THE PATIENT TO BREATHE IN AND OUT FROM THE NOSE.

Use the FLOW knob (Pos A) to adjust the flow.

Gradually adjust the flow, check the status of the bag (Fig. 18 Pos F) which must take on a shape that we will define as **half full half empty**.

Make the patient breathe for at least one minute before **evaluating the flow**.

Once the ideal flow and the pulmonary capacity of the patient have been determined, adjust the gas elimination valve as explained in paragraph 6.4 'Evacuation of exhaled gases'.

It is possible to change the percentage of nitrous oxide using the % knob (Figure 18 Pos C).

The flow meters display the amount of Nitrous oxide and Oxygen that is delivered in litres/min.

The sum of the two flows is the amount initially set based on the patient's respiratory capacity.

The instructions below are examples of administration. Every user should set the flow and percent according to the patient subject to medicine administration.

PROCEDURE SPECIFICATIONS FOR CHILDREN:	PROCEDURE SPECIFICATIONS FOR ADULTS:
<p>The flow for CHILDREN is approximately 4-5 Lt/min.</p> <ul style="list-style-type: none"> - Turn the 'FLOW' knob clockwise to set the flow at 4/5 l/min and instruct the patient to breathe for ~1 minute. - Check the flow setting on the O₂ flow column (Figure 18, Pos E). - The bag must be half full half empty. 	<p>The flow for ADULTS is approximately 6-8 Lt/min.</p> <ul style="list-style-type: none"> - Turn the 'FLOW' knob clockwise to set the flow at 6-8 l/min and instruct the patient to breathe for ~1 minute. - Check the flow setting on the O₂ flow column (Figura 18, Pos E). - The bag must be half full half empty.

Using the % N₂O knob (Figure 18 Pos C) administer nitrous oxide, set it at **20%**, wait a few minutes and see how the patient reacts.

If the patient is not sedated gradually increase the percentage (ex 30%).

THE PERCENTAGE OF NITROUS OXIDE ADMINISTERED TO THE PATIENT MUST BE: THE LOWEST POSSIBLE LEVEL TO OBTAIN GOOD SEDATION (BASE LINE).

Once the patient's BASE LINE has been established, note down the FLOW value and the %N₂O on his/her medical record, this data will be useful for future procedures.

THE GAS MIXTURE THAT IS ADMINISTERED TO THE PATIENT:
MUST NOT EXCEED 70% NITROUS OXIDE
MUST NOT BE LOWER THAN 30% OXYGEN

If the patient reports symptoms of discomfort (nausea, dizziness, head ache) immediately stop administration.
The “FLUSH” button (Fig. 13 Pos. B) is located on the device.
By pressing this button the bag will fill with oxygen.
By pressing the bag we will administer the patient with elevated quantities of pure oxygen.

END OF THERAPY

The treatment will end when the orthodontist has finished the procedure on the patient.

To end treatment:

- Shut off the NITROUS OXIDE supply by turning (counter-clockwise) the %N₂O knob (C – figure 18) to the end stop (0%).
- Administer (as at the beginning of treatment) a flow of 3 Lt/min of oxygen for 2 minutes
- Shut off the OXYGEN supply by turning (counter-clockwise) the FLOW knob (A – figure 18) to the end stop.
- Loosen the loop that secures the tube and take the mask off of the patient.

Before signing the patient out keep him/her under observation (waiting room) for ~5÷10 minutes.

The doctor may deem it necessary to increase the observation period.

TANKS CLOSURE

Once the doctor has finished using the device it is **IMPORTANT** to close the cylinders.

Close the valve on the Nitrous oxide cylinder (N₂O), (turning it clockwise).

Close the valve on the Oxygen cylinder (O₂), (turning it clockwise).

TANKS STORAGE


At the end of the day, store the cylinders in a deposit in compliance with the laws in force.

7.3 RESUSCITATION WITH MASTER FLUX PLUS

*(KIT Art. 1504/A Not included – see figure 19)

RESUSCITATION procedures WITH MASTER FLUX PLUS must be carried out with the %N₂O knob that regulates the percentage of nitrous oxide in position 0%.

SLOWLY open the oxygen cylinder valve (turn the valve counter-clockwise); make sure that the nitrous oxide cylinder is closed (turn the valve clockwise).

Take the MASTER FLUX PLUS mask kit off by sliding it off of the spiral tube.

Insert the relative fitting (30) into the spiral tube.

Remove the bag (**F**) and put the cap (**T**).

Connect the end of the oxygen treatment tube (31) to the fitting (30).

Connect the other end of the oxygen treatment tube to the bottom stud of the resuscitation bag (32) through the relative coupling.

Connect the mask (33), to the valve on the resuscitation bag.

Turn the **FLOW** knob to deliver Oxygen at the required flow rate, CHECK THE FLOWS DISPLAYED ON THE FLOW COLUMNS.

The nitrous oxide flow meter =0

The oxygen flow meter = SET FLOW

Place the mask on the patient.

Act on the resuscitation bag in accordance with the resuscitation instructions.

Once use of the device is finished it is **IMPORTANT** to close the valve on the Oxygen cylinder (by turning it clockwise).


Figure 19

FIGURE	POSITION	Q.TY	CODE	DESCRIPTION
11	29	1	1504/A OPTIONAL	SPIRAL TUBE - CM88051
	30	1		FITTING – SMFA278
	31	1		OXYGEN TUBE – CM84022
	32	1		RESUSCITATION BAG – 1025-2
	33	1		RESUSCITATION MASK – CM84014

8 MAINTENANCE

8.1 CYLINDERS CHANGING

If cylinder pressure is very low ~10÷20 bar, proceed with changing or re-charging the cylinders.

 **We recommend replacing the nitrous oxide bottle every 2 oxygen tanks.**

Slowly close the valves (turn counter-clockwise).

Discharge the system by opening it from the “FLOW” and ‘%’ knob, then close it back up.

Take the pressure reducers out, keeping them in a covered and clean location. (See installation)

Re-charging, checking and periodic inspection of the cylinders must be carried out by specialised companies in compliance with regulations in force in the country of use.

8.2 CLEANING

 **Use detergents that are not aggressive to plastic**

Clean the surfaces of the device with a dry cloth.

Before and after use clean, wipe or sterilise the parts that come into contact with the patient or his/her breath.

 Clean the protection Plaxiglass® of the flow meters with an antistatic detergent

MASTER FLUX PLUS components cleaning table

FIG/POS	MATERIAL	CODE	DESCRIPTION	STEAM STERILISATION	NUMBER OF SERIALIZATION CYCLES
3 / F	NEOPRENE	3MEDN0002	RESERVOIR BAG 3LT	NO	-

 For cleaning the mask (silicone parts and shell) refer to the "Circuits and Masks" user manual.

 **THE SUPPLIED ITEMS VARY ON THE BASIS OF THE DEVICE CONFIGURATION**

8.3 SCHEDULED MAINTENANCE

Check the integrity of the components before each use

If MASTER FLUX PLUS has not been used for over a month make before an operational check.

MASTER FLUX PLUS is complete with security parts that are regulating the correct check of the device. Anyhow we recommend to check the functioning of the device doing the periodical standard tests.

8.2.1 - Function test: check that the device is working properly accordingly the instructions listed in par. 7.1 'Instruction for use'

8.2.2 - Check the floating balls in the flow meters; they must be able to rotate on themselves and move free in the measuring cone when regulating the flow. Be also sure that the set values are stable.

8.2.3 - Check the button FLUSH: press it and check that the N₂O floating ball is going to 0 (zero), oxygen on the other hand must reach max providing level. Leave o₂ FLUSH button: the O₂ supply must be interrupted.

8.2.4 - Check the hold of N₂O supply. Choose gas supply at 50% concentration. Close the bottle as if the bottle would be empty. Check that after some seconds the N₂O supply is completely stopped.

8.2.5 - Check for leaks: Close the "FLOW" knob . Open the valves on the Oxygen and Nitrous oxide cylinders and read the values provided by the pressure gauges. Close the cylinder valves. Wait ~ 5 minutes and read the new values displayed on the pressure gauges. If the values are the same, the device does not have any leaks.

Use only original spare parts and accessories

PRESSURE REDUCERS

FOLLOW THE INSTRUCTIONS CONTAINED IN THE DEVICE MANUALS.

FULL-CIRCUIT PANELS

FOLLOW THE INSTRUCTIONS CONTAINED IN THE CIRCUITS AND DEVICE PANELS MANUALS.

SITUATIONS THAT MAY ARISE DURING OPERATION:

Problems	Possible causes	Solution
It does not deliver oxygen	Empty cylinder	Re-charge the cylinder
	Pressure reducer not working	Change the reducer
	Problems inside the flowmeter box	Have the device checked by Tecno-gaz
It does not deliver nitrous oxide	No oxygen	Charge the cylinder
	Pressure reducer not working	Change the reducer
Nitrous oxide gauge pointing upwards	Incorrect manoeuvre	Clear the 'FLOW' knob to zero Restore correct operation
	Internal problem	Have the device checked by Tecno-gaz

10 PROCEDURES FOR AFTER SALE SERVICE

In case of damage or for overhaul contact directly the phone service of:
 TECNO –GAZ S.p.A.
 Strada Cavalli 4, 43038, Sala Baganza, Parma, ITALIA.
www.tecnogaz.com

PHONE	+39 0521 83.80
FAX	+39 0521 83.33.91
@	cme2@tecnogaz.com

The service staff shall decide whether the return to the manufacturer or the intervention of a technician is necessary; after checking the device the service staff will draw up a quote to be submitted to the distributor which, on its turn, will submit it to the final customer for subscription.

Upon reception of the subscribed quote the device will be repaired and then shipped back within the times indicated on the quote form.

In case of shipment of the device to the manufacturer follow the compulsory indications below:

- Use the original packaging, in case you do not have it anymore, use a suitable packaging. The goods shipment risk is charged to the sender.
- Only send the flow meter box and the pressure reducers, placed individually in clean cellophane.
The reducers must not be contaminated (detergents, alcohol, oils ...). If assistance does not consider the packaging to be suitable, the reducers will be scrapped.
- Write down and include into the package a document describing accurately the anomaly or the service required.
- Free port shipment is required, otherwise the supported transport costs will be charged.

All no original packages we will receive shall be disposed of.

The machine will be shipped back in its original packaging (the cost for the packaging will be charged to you) through the customer's forwarder.

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TECNO-GAZ S.p.A.