Istruzioni per l'uso  
Instructions for use  
Mode d'emploi  
Bedienungsanleitung  
Instrucciones para el uso

# R400

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The “R400” oxygen reducer has been designed and constructed under the 93/42/EC directive in full compliance with the UNI EN ISO 10524-1 technical standard requirements.

The test performance is evidenced by the application of the CE marking plate, and by the declaration of conformity enclosed to this Instructions Manual.

On the reducer body the production parcel is punched, see PICT. 2.

The adhesive tape wrapping the reducer reads: the manufacturer's name, the outlet rated Pressure ( $p_2$ ) the inlet rated Pressure ( $p_1$ ), the Reference flow ( $q_1$ ), the Safety valve opening pressure ( $p_{RV}$ ), the CE marking followed by the number of the Body certified to the release of the product declaration of conformity under the 93/42/EC directive, see PICT. 2.

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Technical service and repair interventions are carried out by Tecno-Gaz at its own test and repair laboratory.

TECNO-GAZ S.p.A.  
Strada Cavalli n°4 – 43038 Sala Baganza – Parma – Italia  
Tel. +39 0521-833926 – Telefax +39 0521-833391  
e-mail: [info@tecnogaz.com](mailto:info@tecnogaz.com)  
<http://www.tecnogaz.com>

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**WARRANTY**

The oxygen pressure reducer R400 is covered by a 12-month warranty, except for wear parts (fittings).

This warranty is valid starting from the date the reducer is delivered to the customer, evidenced by the return of the warranty coupon duly filled-in, stamped and signed by the retailer.

In case of controversy, the date indicated on the delivery note will be considered that valid.

Repairing and replacement under warranty conditions of any components will be carried out exclusively at the manufacturer discretion; the warranty does not cover travel expenses of personnel, neither packaging nor transportation costs.

The warranty does not cover damages caused by wrong maintenance, negligence, unskilfulness or any other cause the constructor is not to be held responsible for, and consequently no kind of compensation is acknowledged to people or property.

The right to have the whole reducer replaced is not acknowledged.

No compensation for the reducer stop is acknowledged.

The warranty will void automatically in case the reducer is tampered with, repaired or modified by the customer or by unauthorised third parties.

For interventions, the customer has to apply exclusively to the retailer or to the service centres of the manufacturer.

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See PICT. 6

**WEIGHT:**

With female 1/4" G outlet: ..... 0.885 kg  
 With diaphragm: ..... 0.910 kg


$p_1$ : ..... 200 bar (20.000 kPa)  
 $p_2$ : ..... 4 bar (400 kPa)  
 $p_3$ : ..... 9 bar (900 kPa)  
 $p_4$ : ..... 4,6 bar (460 kPa)  
 $p_{4MAX}$ : ..... 5 bar (500 kPa)  
 $p_5$ : ..... 4,6 bar (460 kPa)  
 $q_1$ : ..... 40 NI/min  
 $p_{RV}$ : ..... 6 bar (600 kPa)

LENGTH - DIAMETER: ..... See PICT.(3)

CONSTRUCTION MATERIALS: Brass - Viton - Tecaflon

- $p_1$ : Infeed rated pressure
- $p_2$ : Outfeed rated pressure referring to  $p_3$  and  $q_1$
- $p_3$ : Minimum infeed pressure
- $p_4$ : Outfeed maximum pressure referring to  $p_3$  and  $q=0$
- $p_{4MAX}$ : Outfeed maximum pressure referring to  $p_1$  and  $q=0$
- $p_5$ : Outfeed maximum pressure referring to  $p_1$  and  $q_1$
- $q_1$ : Rated flow
- $q_d$ : Flow delivered by the diaphragm reducer
- $p_{RV}$ : Safety valve opening pressure

OUTFEED PRESSURE ACCORDING TO INFEED PRESSURE IF  $q_1 = 40$  NI/min

**PAR. 6****IMPORTANT**

-   **Explosion risk !!!!! Do not lubricate any part or accessory of the reducer.**
-  The R400 reducer must be used by duly trained staff.
-  Carefully read this Manual instructions for a correct use of the PRESSURE REDUCER.
-  Do not attach the reducer to a bottle containing gases other than Oxygen.
-  Check for the device integrity after unpacking.
-  Use, packaging, transport, and storage temperature ranges between (- 20°C) and (+60°C).
-  Before using the reducer and after every mounting or intervention on the reducer check its correct operation following control instructions described in PAR.11.
-  Every 2 months check for the efficiency and the correct operation of the REDUCER as per PAR.11.
-  Every 6 months check the outlet pressure, as per PAR.13.
-  Do not stress mechanically the reducer when applied to the bottle (ex.: do not lean on it, do not hang objects, do not seize to lift the bottle etc..) to avoid breaking the device.
-  Be sure the bottle is secured correctly and steady.
-  Use the reducer under strict visual control.
-  Position the manometers so that the user can easily see them.
-  Fittings for the connection to bottles with O-rings seals can be clamped without wrench.
-  Nuts connecting the reducer to the bottle with rigid seals have to be clamped with suitable wrench not to damage the hexagon.
-  TECNO-GAZ is not responsible for a wrong use of the device.
-  Always close the bottle valve when the device is not used.
-  Vent the gas inside the reducer in safe areas keeping at a distance from flames, flammable materials or heat sources.
-  Do not unscrew the reducer nut from the bottle valve when the reducer is under pressure.
-  Vent the reducer before unscrewing from bottle.
-  Slowly open the bottle valve.
-  Do not fit foreign bodies into the vent valve.
-  Do not unscrew or loosen parts of the reducer that are not mentioned in the instructions.
-  After use, store the bottle and the reducer in a safe and steady place, far from heat sources and electric circuits.
-  Before connecting the reducer to the bottle visually check (without disassembling) for the absence of impurities held by the filter (Pict.5). If residues are present apply to the Tecno-Gaz service centre.
-  Every three years the reducer must undergo a thorough overhaul.
-  Do not touch with greasy hands and do not pour liquids into the connection inlet hole

**PAR. 7****PACKAGING, TRANSPORT AND STORAGE**

**The use, packaging, transport, and storage temperature ranges between (-20°C) and (+60°C).**

The reducer is bagged into a cellophane fitted into shaped protection housed in a corrugated cardboard box.

The box also contains the Use Manual. If the reducer is equipped with flowmeter and humidifier see the respective user manuals enclosed.

In case of long stop put the reducer in its original packaging and store it in a dry place, protected against bad weather at a temperature of (-20 to +60)°C, relative humidity 40 - 70%

**PAR. 8****USE OF THE REDUCER**

**The R400 reducer must be used by duly trained staff**

The reducer has to be used for delivering medical oxygen.

The R400 reducer is designed and manufactured for direct connection to the bottle.

The nut and the inlet connection vary according to the use country.

1. Connect the reducer by screwing the nut to the bottle valve, See PICT.4. Tighten only with hands, without any spanner.
2. Open the bottle valve very slowly by rotating it counterclockwise. Visually check the reducer.
3. Carry out a leak test as per instructions of PAR.11.
4. Start delivering.
5. After use close the bottle.
6. Safely vent the gas inside the reducer under PAR.6.
7. After use store the bottle with the reducer in a safe and steady place, far from heat sources and electric circuits.

**PAR. 9****MAINTENANCE**

**For an effective and safe operation of the REDUCER the user has to carry out some periodical checks.**

**Should these checks indicate the need of repair interventions apply to the closer TECNO-GAZ service centre.**

1. Carefully clean the reducer surfaces following instructions of PAR.10 "CLEANING".
2. Visually check the O-ring wear, if the seal is worn or damaged replace it as follows:
  - a) Remove the worn or damaged O-ring using a tool and be careful not to scratch the O-ring housing in the inlet connection, see PICT.5.
  - b) Carefully fit the new O-ring into its housing.
  - c) Be sure the O-ring is correctly fitted.
  - d) You can now use the REDUCER.

**PAR. 10**

**CLEANING**

- 1 - Clean the device 1 per week using a soft cloth slightly soaked in alcohol.
- 2 - Do not touch with greasy hands and do not pour liquids into the connection inlet hole.

**PAR. 11**

**LEAK TEST**

 Every 2 months check for the REDUCER efficiency and correct operation.

**1. LEAK TEST:** (See PICT.4)

- a) Close gas delivery (close the valve of the flow meter or of any other downstream device).
- b) Open the bottle valve rotating the knob counterclockwise.
- c) Write down the pressure value on the manometer.
- d) Close the bottle valve rotating the knob clockwise.
- e) Wait for 5 minutes.
- f) Check pressure on the manometer.
- g) If pressure at point f) equals pressure at point c) you can use the reducer.
- h) If pressure at point f) is lower than pressure at point c), follow instructions below.
  - Slowly open the bottle valve rotating the knob (A) counterclockwise.
  - Check for gas leak using a bubble promoter.
  - Bubbles form on the inlet connection/nut
  - Close the bottle valve rotating the knob clockwise.
  - Carefully dry the reducer
  - Vent the reducer
  - Dismount the reducer from the bottle valve unscrewing the nut.
  - Be sure the O-ring is not worn or damaged, in case of worn or damaged seal replace it following instructions in PAR.9.
  - Mount the reducer on the bottle valve as described in PAR.8.
  - Restart from point a).

 **IF BUBBLES FORM ON THE REDUCER BODY DO NOT USE THE REDUCER.**

Apply to the TECNO-GAZ service centre

➤  **IF THE LEAK IS DETECTED ON THE OVERPRESSURE VALVE, PICT 4, DO NOT USE THE REDUCER.**

Apply to the TECNO-GAZ service centre

**PAR. 12**

**CONNECTION TO OTHER MEDICAL DEVICES**

The PRESSURE REDUCER can be connected to other medical devices to check oxygen delivery, such as:

- Flow meter, Hose connection, Hose connection with diaphragm 6NI/min, Hose connection with diaphragm 10NI/min, Hose connection with diaphragm 15NI/min, Flow meter with humidifier, Flow meter, Couplers, Terminal units, etc...

For the connection of the reducer to other medical devices the user has to comply with the following instructions:

- a) Connect the pressure reducer only to medical devices with the CE marking (93/42/EC).
- b) Use a dynamometric wrench for clamping.
  - The torque value to clamp a coupler, hose connection or diaphragm to the reducer body must be higher than or equal to 20 Nm.
  - The torque value to connect manometers or flow meters to the reducer body must be higher than or equal to 20 Nm.
- c) Carry out checks suggested by the manufacturer on every device before connecting to the reducer.
- d) Check correct operation after connection.

**PAR. 13**

**OUTLET PRESSURE CHECK**

Every **6** months check the outlet pressure using a sample manometer with a 12 bar (1.200 kPa) full scale to be connected to the reducer outlet.

The manometer has to show a steady closing pressure of (4 to 5)bar =(400 to 500)kPa for 5 minutes.


**IF PRESSURE INDICATED ON THE MANOMETER IS HIGHER, DO NOT USE THE REDUCER**

Apply to the TECNO-GAZ service centre

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




**GB**

This manual must always be kept with the product, in complying with the Directives of European Community.

TECNO-GAZ reserves the right to modify the enclosed document without notice.

TECNO-GAZ reserves the property of the document and forbids others to use it or spread it without its approval.

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**TECNO-GAZ S.p.A. a socio unico**