**MANUALE UTENTE  
USER MANUAL  
ANWENDERHANDBUCH**

# SKYLAB

**Art. 4500/S**

**ITALIANO**

**ENGLISH**

**DEUTCH**

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


- Sensore

**R:**

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Declaration Of Conformity	3
Warranty condition	4
Loss of warranty	5
Delivery terms	6
General warnings and safety rules	7
Useful contacts and addresses	7
Symbols	7
Technical Features	9
Instructions For Packaging, Transportation, Storage And Lifting	10
Machine Use	10
Operator's Safety	11
Operating Principle	12
Machine Description	13
List Of Mechanic Spare Parts	14
Installation Instructions	17
List Of Supplied Accessories	17
Control List	17
Materials Used For Sandblasting	18
List Of Tool Spare Parts	19
Tool Description And Operation	20
Instructions For The Machine Use	22
Emissions Of Toxic Substances During Operation	24

Routine Maintenance	25
Extraordinary Maintenance And Repair	27
Instructions For Machine Demolition	27
List Of Electric System Components	28
List Of Components Of The Pneumatic System	28
Summary Of Instructions To Comply With	29
Drawings And Diagrams	30

**DECLARATION OF CONFORMITY**

( according to ISO/IEC 17050-1)

Complies with the following pertinent regulations in force: Machine Directive 2006/42/CEE, Electromagnetic Compatibility Directive 2004/108/CEE, Low Tension Directive 2006/95/CEE and their amendments and updating.

The original declaration of confirmity is provided in attached to the manual.

**WARRANTY CONDITION**

- 1) **DURATION:** the product is covered by a 12-month (twelve) guarantee.
- 2) **EXCLUSION:** the guarantee excludes:
  - a) faults originating from non-performance of scheduled maintenance and due to negligence or improper use of the product by the User;
  - b) periodic checks and maintenance;
  - c) repair or replacement of parts subject to wear that are fragile or have an unpredictable life, unless they prove faulty at the time of their delivery;
  - d) faults caused by personnel or originating from technical personnel transfer and transportation;
  - e) faults or damage originating from improper use or operating errors;
  - f) faults or damage originating from contamination in the water and air infeed systems, extraordinary chemical agents or electrical events;
  - g) faults or damage originating from the use of detergents, disinfecting agents, sterilizing fluids or processes that are not described explicitly in the instruction manual for operation and maintenance;
  - h) natural color change of plastic components.
- 3) **INSTALLATION AND TESTING REPORT:** an essential prerequisite for guarantee validity is the return of the installation and testing report after all its sections have been filled in by the User and the Fitter. The report must be returned within 15 days after installation, otherwise the guarantee becomes null and void.
- 4) **LIMITS:** the guarantee entitles its holder to repairs or replacement of faulty components free of charge. The right to replace the entire machine is excluded. As for components either applied or added to TECNO-GAZ products or by third parties, having their own guarantee certificate, the conditions, limits and exclusions indicated in the relevant certificates apply.

Under no circumstances shall carrying out of one or more than one repair during the guarantee term change the guarantee expiry date.
- 5) **CLAIMS:** in the event of the Purchaser making a claim relating to the application of the guarantee and the quality or conditions of the delivered equipment, the said Purchaser shall in no way suspend and/or delay payments.
- 6) **DISPUTES:** if a dispute arises on the application and interpretation of this Certificate of guarantee, the Court of PARMA (Italy) shall be the court of jurisdiction, regardless of the location in which the purchase agreement was stipulated.
- 7) **DEROGATIONS:** any special derogations from the guarantee terms described in this section shall not imply the granting to the Purchaser of any right and shall be considered as granted for a specific case only.

8) MISCELLANEOUS: for any other issue, which is not covered by this Certificate of Guarantee, the rules described in the Civil Code and existing relevant Laws in the Republic of Italy apply.

## LOSS OF WARRANTY

The guarantee becomes null and void in the event that:

- a) the equipment presents damage caused by a fall, exposure to flames, liquid spillage, natural events, bad weather conditions or by any other cause that is not due to manufacturing defects;
- b) installation has not been performed in line with TECNO-GAZ instructions and by authorized personnel;
- c) the equipment is repaired, modified or tampered with by the Purchaser or by non-authorized third parties;
- d) when asking for an intervention under guarantee, the product serial number has been removed, deleted, counterfeited, etc.;
- e) the installation and testing report is not returned suitably filled in and signed within 15 days from installation;
- f) the Purchaser suspends and/or delays payment of amounts due, for any reason whatsoever, relating to the purchase and/or maintenance of the equipment;
- g) scheduled periodic maintenance or any other instruction envisaged in the manual for operation and maintenance is not carried out.

**DELIVERY TERMS**

1. All goods delivered must be packaged in their original packaging. If the goods are delivered in a non-genuine packaging, a charge will be applied to restore the original packaging.
2. All goods must be delivered CARRIAGE PAID.
3. The client shall ask a written authorization to TECNO-GAZ sales dept. before returning the goods by means of the "Form for Client's Returned Goods Management" (CM-M-002) in compliance with the requirements of TECNO-GAZ Quality Management System (ref. CM-P-003, CM-I-009).
4. The document accompanying the returned goods shall show:
  - the TECNO-GAZ product ode/item no.;
  - the number of the TECNO-GAZ purchase invoice or shipment document;
  - the product serial number and batch number;
  - a description of the goods returned;
  - the number of the Form for Client's Returned Goods Management" (CM-M-002) or a copy of the same attached.
5. Returned goods are considered as accepted only after they have been examined. If they prove non-conforming, they may be sent back to the consignor at TECNO-GAZ full discretion.
6. Charges will be applied to restore and test the returned product.
7. Parts replaced under guarantee must be returned, carriage paid, to TECNO-GAZ.
8. If the parts are not returned, their cost will automatically be debited to the recipient.
9. TECNO-GAZ does not accept goods returns from end users.
10. Goods sent to TECNO-GAZ for repair are managed as returned goods.

## GENERAL WARNINGS AND SAFETY RULES

-  Make sure that the equipment is supplied with the correct voltage as shown in the plate.
-  Do not remove the plate.
-  Make sure that the machine is equipped with a grounding system.
-  Use a dry cloth to clean the machine.
-  Before any operation is performed, disconnect the power cable from the socket.
-  Use original spare parts only.

*Non-compliance with the rules above relieves TECNO-GAZ of all responsibility.*

## USEFUL CONTACTS AND ADDRESSES

For technical service please apply to the dealer that invoiced the device, or directly to:

TECNO-GAZ SpA

Str. Cavalli n°4 – Sala Baganza – Parma – ITALIA.

Tel. +39 521 8380 - Fax +39 521 833391

e-mail : [info@tecnogaz.com](mailto:info@tecnogaz.com)

<http://www.tecnogaz.com>

## SYMBOLS

Pursuant to DIRECTIVE 2002/96/EC, this symbol shows that the product must not be disposed of as urban waste at the end of its operating life.

General warning and information for the user.

Pay attention to the instructions provided alongside this symbol.

Conformity with the Directive 2004/108/CE and Directive 2006/95/CE and further amendments.

**PE**

External protection.

Protection GROUNDING SYSTEM.

(Operating) GROUNDING SYSTEM.

Alternate current.

Warning: high voltage. Risk of electrocution!

Read the instruction manual carefully.

**POS.**

Position (number identifying a component in the figures or diagrams).

**Fg.**

Figure

**$\Omega$**

Ohm (Electric Resistance unit).

**s**

Seconds (time unit).

**W**

Watt (Power unit).

**Hz**

Hertz (Frequency unit).

**mm**

Millimeters (Length unit).

**A**

Ampere (Electric Current unit).

**V**

Volt (Electric Voltage unit).

**Pa**

Pascal (Atmospheric Pressure unit).

**°C**

Degrees Centigrade (Temperature unit).

**kg**

Kilograms (Weight unit).

**TECHNICAL FEATURES**

<b>MECHANIC:</b>		<b>ELECTRIC:</b>	
Height	600mm	Supply Voltage	230 V±10%
Width	600 mm	Frequency	50 Hz
Depth	530 mm	Phases	1 + GND
Weight	36 Kg $\cong$ 352 N	Fuse Current	5 A
Material	Iron	Maximum Power	1 KW
Sand tank capacity	500cc	Protection	PE
Sand tank material	PVC	Class	1
		Power Cord 2 M	2X1.5 mm <sup>2</sup> +GND
		Connector (Plug)	16 A + GND

<b>SOUND POWER LEVEL:</b>		<b>PNEUMATIC:</b>	
Standard Operation Noise (Iso/Dis 3746))	81.4 dB(A)	Feeding	Dry non lubrified compressed air
Maximum Noise (Iso/Dis 3746)	96.0 dB(A)	Feeding Pressure:	8 Bars
		Max Consumption 12mc/H	200 NL/Min

### INSTRUCTIONS FOR PACKAGING, STORAGE, TRANSPORTATION, AND LIFTING

The Sandblasting machine has a weight of 36 kg and a compact size, hence it needs no special packaging.

It is bagged into a cellophane bag and then fitted into a foam polystyrene protection to be housed into a corrugated board box sealed with adhesive tape.

During transportation the package has to be protected from bumps or shocks and humidity.

#### ATTENTION:

Do not disperse packaging materials into the environment; recycle reusable products and give the remaining ones to a company specialising in waste disposal.

Packed machines have to be stored in dry places at a temperature of -10 to + 40°C.

Machines cannot be stacked.

#### CAUTION:

**The only seizing point for the manual lifting and transportation of the machine are the holes in the gloves bushings.**

**Do not lift the machine through the side doors.**

### MACHINE USE

Is a small compressed air-fed sandblasting machine to be used for the surface treatment of objects.

Such objects have to be of correct size so as to be housed in the closed cabins and allow a safe operation (closed door).

During processing, they must not react with materials used for sandblasting, that is to say they do not have to send out gas substances (that cannot be captured by filters).

The correct operation of the machine depends on the absence of humidity inside the cabin, consequently no wet object or fluid has to be introduced to avoid sands wetting.

## OPERATOR'S SAFETY

During the sandblasting machine operation the operator may run two types of risk:

- 1) If handling the tool when the door is open, he could aim the jet to himself and send sand into his eyes.
- 2) If the sandblasting treatment is carried out with the door open, the dust suction does not occur and, consequently, dusts pollute the workplace.

To eliminate such risks, the sandblasting machine is equipped with two safety systems inhibiting the machine use in case one of the door is open or the dust aspirator is not working.

**The first safety device is linked to the door lock pict. 17 and includes a pneumatic sensor and a sensor indicator located on the locking lever.**

**The machine works only if every door is closed.**

**If the sensor's checking doesn't close the little nozzle with continual loss of the sensor, the sensor doesn't allow the pneumatic feeding of tools.**

**To lock the door approach it to the seal, press the lock handle downward and rotate it clockwise by 90° - see pict. 15-16 – then release the handle.**

**Once the operation is completed, be sure the handle cannot rotate in the opposite direction unless it is pressed downward again.**

**The sensor indicator is a wear part, periodical checks are suggested.**

**The second safety device is a check of the aspirator start. Only in these conditions sandblasting is possible.**

**The aspirator conveys the dusts deriving from the sandblasting processing to a set of filters. If the maintenance is correct, the air getting into the environment will be clean.**

**OPERATING PRINCIPLE**

The operation of the mini sandblasting machine is characterised by two actions:

The action of the sand on the object and the action of the aspirator on dusts produced during the processing.

Two independent circuits allow the accomplishment of the above described operations.

Sand contained in the tank **pos 140 pict 2** is sucked by an ejector housed in the tool **pict 4**.

Through the venturi effect, the air conveys sand into the tool nozzle and gives it a speed.

Coming out the nozzle **pos 202 pict 4**, sand is thrown to the object.

The sand speed, weight, and toughness abrade the object and thus remove part of the material and fragment the sand.

Heavy particles fall into the container and are then continuously recycled.

Lighter dusts are removed from the cabin through the suction mouth positioned on the lamp, **pos 69 pict 2**.

The aspirator keeps the cabin atmosphere clean and conveys air and dusts in a filtering bag **pos 253 pict 1**.

Dusts are filtered while the air passes through the bag and is conveyed to a second filter, **pos 230 pict 1**, that further cleans the air before releasing it into the atmosphere.

## MACHINE DESCRIPTION

Is a mini sandblasting machine with three independent cabins and allows the simultaneous work of three operators.

Each cabin is equipped with a tool with built-in control button **pos 215-217 pict2**, a compressed air gun to clean objects **pos 53**, a lighting lamp **pos 6-7**, a rack for objects **pos 118**, a sand tank **pos 140** and a sealed transparent door **pos 227 pict1**, for an easy visual inspection.

Gloves are of the super-resistant type, (any type of glove can be fitted in) and are fixed to the cabin by means of a screw ring nut that allows an easy replacement.

The processing dusts are sucked by the turbine **Pos 5 pict 2** and held by a filtering bag **pos 253 pict 1**.

Before reaching the environment the air is cleaned by the final filter **pos 230 pict 1** this filter is made in special paper.

In front of the cabin there is the control panel **pict 3** housing suction and lighting switches, the pressure reducer used to adjust the tool pressure and the manometer displaying the pressure value.

The three cabins can perform a 120° rotation to the right or to the left with respect to the central cabin marked by the writing SKYLAB.

This manual rotation is obtained by means of several bearings connected to a rest **pos 130 pict 1** supporting the whole machine.

Thanks to this movement an operator can have access to each of the three cabins without moving around the machine.

The cabin cleaning and the contaminated sands removal is carried out by suction.

The machine is equipped with a suction nozzle hose **pos 251-252 pict 1** to be fitted into the gloves safety valve **pos 175 pict 1** to carry out the above described operation.

The handle **pos 3 pict 1** located on the central cabin left hand side is used to adjust the turbine power from zero to the maximum value.

For dusts suction during processing the turbine power is kept at a low value (enough to guarantee a good visibility inside the cabin) while for the cleaning phase the turbine power has to be the maximum one.

**LIST OF MECHANIC SPARE PARTS**

<b>Pos</b>	<b>Q.ty</b>	<b>Description</b>	<b>Cod.</b>
<b>105</b>	3	Side door	SSKA105
<b>112</b>	1	Complete furniture	SSKA112
<b>117</b>	1	Tray	SSKA117
<b>118</b>	1	Central grid -	SSKA118
<b>119</b>	1	LH grid	SSKA119
<b>120</b>	1	RH grid	SSKA120
<b>121</b>	1	Motor adjuster bracket	SSKA121
<b>122</b>	1	Filter box	SSKA122
<b>123</b>	1	Filter lid	SSKA123
<b>124</b>	3	Gloves counter nut	SSKA124
<b>130</b>	1	Rest	SSKA130
<b>131</b>	1	CE plate	SSKA131
<b>133</b>	1	Sand tank plug	SSKA133
<b>134</b>	1	Sand tank sieve	SSKA134
<b>135</b>	1	Turbine protection	SSKA135
<b>136</b>	1	Turbine seal bushing	SSKA136
<b>140</b>	3	Sand tank	SSKA140
<b>141</b>	1	Filter plug	SSKA141
<b>142</b>	1	Induction manifold	SSKA142
<b>143</b>	3	Gloves bushing	SSKA143
<b>144</b>	3	Gloves ring nut	SSKA144
<b>149</b>	3	RH glass hinge plate	SSKA149

<b>150</b>	3	LH glass hinge plate	SSKA150
<b>153</b>	3	Bracket	SSKA153
<b>154</b>	3	Rotation stop	SSKA154
<b>155</b>	3	Positioner	SSKA155
<b>156</b>	3	Lower bearing support	SSKA156
<b>157</b>	6	Bearing spacer	SSKA157
<b>158</b>	1	Lower cable lock	SSKA158
<b>159</b>	1	Upper cable lock	SSKA159
<b>160</b>	2	Stop striker	SSKA160
<b>161</b>	6	Hoses bushing	SSKA161
<b>162</b>	6	Nut for bushing	SSKA162
<b>163</b>	12	Washers for bushing	SSKA163
<b>164</b>	6	Sand plunger	SSKA164
<b>165</b>	3	Turbine spacer	SSKA165
<b>166</b>	18	Inserts for sand tank	SSKA166
<b>167</b>	3	Reflector plate	SSKA167
<b>168</b>	1	Cable pressing washer	SSKA168
<b>169</b>	1	Cable pressing washer	SSKA169
<b>170</b>	2	Suction hose terminal	SSKA170
<b>175</b>	1	Gloves valve	SSKA175
<b>215</b>	1	Big tool diam. 3.5 complete	SSKA215
<b>217</b>	2	Big tool diam. 1.5 complete	SSKA217
<b>226</b>	3	Control panel sticker	SSKA226
<b>227</b>	3	Door glass	SSKA227

<b>228</b>	3	Door seal	SSKA228
<b>229</b>	3	Protective film	SSKA229
<b>230</b>	1	Final filter	SSKA230
<b>233</b>	3	Cabin sticker	SSKA233
<b>234</b>	1	Filter sticker	SSKA234
<b>235</b>	1	Gloves valve spring	SSKA235
<b>236</b>	6	Spring for lower bearing support	SSKA236
<b>237</b>	3	Spring for positioner	SSKA237
<b>238</b>	6	Glass locking washer	SSKA238
<b>239</b>	3	Cabin lock	SSKA239
<b>241</b>	3	Plug handle	CM30005
<b>243</b>	3	RH hinge	SSKA243
<b>244</b>	3	LH hinge	SSKA244
<b>251</b>	1	Suction nozzle	CM61002
<b>252</b>	1	Suction hose	CP01016
<b>253</b>	1	Filtering bag	SSKA253
<b>254</b>	3	Gloves (Pairs)	CM61001
<b>255</b>	1	Turbine adjustment knob	CM30008
<b>256</b>	1	Hose elastollan diam. 8/6	CP01009
<b>257</b>	5	Anti-vibration rests	CM70002
<b>258</b>	3	Cable stops diam. 6	CE02011
<b>259</b>	3	Spiral hose diam. 4/2	CP01017
<b>261</b>	1	Turbine pipe	CP01015
<b>262</b>	1	Handle hood	CM30009
<b>276</b>	8	Electric cable seals diam. 8	CE02015

<b>277</b>	1	Turbine seal	CM50054
<b>278</b>	3	Sand tank seal	CM50002
<b>280</b>	2	Suction hose seal	CM50018
<b>281</b>	2	Suction plug seal	CM50001
<b>283</b>	1	Filtering bag seal	CM50003
<b>284</b>	7	Seal nylon 1/8	CM50007
<b>285</b>	1	Glove valve seal	CM50012
<b>290</b>	9	Bearing 608 2Z	CM60003
<b>291</b>	2	Turbine pipe clamps	CM17005

## INSTALLATION ISTRUCTIONS

The Sandblasting machine has to be positioned on a stand at a suitable height for the operator.

In case of simultaneous use by operators the machine has to be positioned in the middle of the room or in a position allowing the presence of several people around it. A part in case of the above-described conditions, can be located in any position provided the access to one of the cabins is possible.

Cabin rotation, allows the use of any of the cabins, consequently the machine should be positioned so as to allow the rotary movement.

Be sure the stand supporting the machine has a capacity suitable to the machine weight (at least 50 kg) and is steady.

Has to be connected to:

a monophase electric socket with 220 V and 50 Hz frequency. Such socket has to be equipped with ground wire, magnetothermal switch with differential (overload cutout) and has to distribute a maximum power of 1 Kw.

To a DRY and NON LUBRICATED compressed air source with a minimum pressure of 7 BARS and a distribution capacity of at least 200 nl/min.

## LIST OF SUPPLIED ACCESSORIES

The Sandblasting machine comes with the following accessories:

N° 1 Suction hose with nozzle	<b>Pos.251-252 Pict.1</b>
N° 2 Filtering bags	<b>Pos.253 Pict.1</b>
N° 3 Glass protective films	<b>Pos.229 Pict.1</b>
N° 1 Pritt roller (glue)	
N° 1 CE Use and Maintenance manual	
N° 1 Plug for sand tank	<b>Pos.133 Pict.2</b>
N° 1 Sieve for sand tank	<b>Pos.134 Pict.2</b>

## CONTROL LIST

All cabins are equipped with a control panel that includes:

**Pos. 1** Green luminous switch for cabin lighting.

**Pos. 2** Green luminous switch for dusts suction.

**Pos. 51** Pressure reducer to adjust the tool pressure.

**Pos. 52** Manometer showing the tool pressure.

• **Suction force** is adjusted using the knob **pos 3 pict 1** (SUCTION ADJUSTER).

## MATERIALS USED FOR SANDBLASTING

Abrasives used for sandblasting are numerous and each has different features to satisfy all requirements of the dental sector.

The most commonly used are:

CORUNDUMS used to remove investment and for metal surface treatment.

ALUMINIUM OXIDES used for pickling and frames preparation.

CRYSTAL MICROBALLS for glazing.

The above-described abrasives are available in different granulometries, the selection of the abrasive depends on the grade of processing to be obtained.

Coarse grain sands incise remarkably the object and create on the surface a high roughness while fine grain sands, although removing material, smooth surface.

The abrasion grade depends on the features of the used sand with regard to the processed type of metal and on the speed at which the sand strikes the object.

The machine gives sands the maximum speed when tools are fed at a 4-5 Bars (BIOXIDE and BALLS tools) and 6-6.5 Bars (CORUNDUM tool) pressure.

The sandblasting machine can be used with all abrasives available on trade provided their size does not exceed 500 my; for very fine sands, lower than 25 my, their compatibility with the venturi suction system has to be checked for.

**LIST OF TOOL SPARE PARTS**

Pos.	Q.ty	Description	Cod.
<b>176</b>	3	Big handle	SSKA176
<b>178</b>	3	Big tool body	SSKA178
<b>180</b>	1	Big tool tip	SSKA180
<b>181</b>	2	Big tool tip	SSKA181
<b>183</b>	1	Air nozzle diam. 2	SSKA183
<b>184</b>	1	Air nozzle diam.1.2	SSKA184
<b>185</b>	1	Air nozzle diam. 0.8	SSKA185
<b>188</b>	3	Control lever	SSKA188
<b>189</b>	3	Button rod	SSKA189
<b>190</b>	3	Button plug	SSKA190
<b>191</b>	3	Hose coupler M 6	SSKA191
<b>192</b>	3	Hose coupler M 5	SSKA192
<b>193</b>	3	Washer diam.4x9	CM20002
<b>194</b>	3	O-ring 2012	CM50036
<b>195</b>	6	O-ring 2010	CM50050
<b>196</b>	3	O-ring 2056	CM50038
<b>198</b>	3	O-ring D.	CM50068
<b>199</b>	3	O-ring 2018	CM50048
<b>200</b>	3	Spring for big tool	SSKA200
<b>202</b>	1	Mix nozzle diam.7x3.5	SSKA202
<b>203</b>	1	Mix nozzle diam.6x2.5	SSKA203
<b>204</b>	1	Mix nozzle diam 4x1.5	SSKA204
<b>206</b>	3	TCEI screw M3X8	CM15037
<b>207</b>	3	Hose elastollan neutral	CP01008
<b>209</b>	3	Hose elastollan light blue	CP01005
<b>211</b>	1	Big tool tip	SSKA211
<b>WHOLE TOOLS</b>			
<b>215</b>	1	Big tool diam. 3.5	SSKA215
<b>216</b>	1	Big tool diam. 2.5 (ONLY ON DEMAND)	SSKA216
<b>217</b>	1	Big tool diam. 1.5	SSKA217

## TOOLS DESCRIPTION AND OPERATION

The tool is comprised of:

Handle **pos 176**.

Tool body **pos 178**.

Control lever **pos 188**.

Air nozzle **pos 183-4-5**.

A tip with mix nozzle **pos 180-1**.

The tool body houses an ejector that sucks the sand from the container, conveys it to the tip and ejects it through the hard metal mix nozzle, **pos 202-3-4**.

Each mix nozzle corresponds to an air nozzle **pos 183-4-5**; if the diam. 3.5 **Pos 180** tool tip is replaced with a diam. 1.5 **pos 181** tip, the diam. 2 **pos 183** nozzle too has to be replaced with a diam. 0.8 **pos 185** nozzle.

### ATTENTION:

When replacing the air nozzles **pos 183-4-5**, be sure the O-rings **pos198** are correctly housed with a perfect tightness.

The tool tip is screwed to the body; before screwing, remove dust, if any, from threads and seal **pos 196**, using a compressed air gun, then cover with talc to reduce the seal – screw friction.

If the stick is pressed, **pos 188**, the air flows out from the nozzles thus creating the vacuum inside the tip and, consequently, sucking the sand through the white tube **pos 207**.

The maximum vacuum value is 0.7 Bar corresponding to a tool feeding of 4-5 bars (diam. 0.8 nozzles) **pos 185** and 6-6.5 Bars (diam. 2 nozzles) **pos 183**.

**An increase or a decrease in the feed pressure to the above mentioned values causes a drop in the ejector yield and consequently vacuum and dust suction drop too.**

**A further element affecting the tool yield is the position of the tip with regard to the air nozzle; indeed, if the tip is completely screwed to the body it prevents the sand from reaching the venturi.**

**In the above-described condition only air flows from the tip.**

**After unscrewing the tip, the sand flows through the venturi and is ejected.**

**The highest yield point is reached by unscrewing the tip of one or two turns to the “completely screwed” position pict 5.**

**When the sand jet is non constant but intermittent, it means the tip position is not correct.**

### IN CASE OF OBSTRUCTION OF THE TOOL DUCTS:

Before unloading the container and ducts we suggest to try to eliminate the obstruction by capping the tool tip by a finger toe.

In this situation, by pressing the control lever, ducts are pressurised and sand is forced to go back to the tank **pict 5**.

In case the above described intervention does not solve the problem, remove the tool tip and check for the presence of foreign bodies.

### ATTENTION:

These operations have to be performed as if a processing phase, that is the door closed and the hands fitted into the gloves.

If there is no foreign body, the obstruction must be located near the sand tank.

In this case, the intervention with open door is necessary and it is hence necessary to stop the pneumatic feeding, unload the container and identify the obstruction cause.

## INSTRUCTIONS FOR THE MACHINE USE

After machine installation, remove the plug **pos 141 pict 1** (rotate it to the right or to the left while pulling it upward) and be sure that the filtering bag **pos 253** is correctly housed **pos117**, otherwise, fit the bag routine maintenance.

The first cabin, c, marked on the glass by the "**CORUNDUM**" indication, is equipped with a big tool with a 3.5 mm nozzle. This cabin can process course sands up to **500 my**.

The second cabin is positioned to the right of the first one and is marked by "**BIOXIDE**", and is equipped with a big tool with a 1.5 mm nozzle. This cabin can process sands up to **125 my**.

The third cabin is positioned to the left of the first one and is marked by "**BALLS**", and is equipped with a big tool with a 1.5 mm nozzle. This cabin can process sands up to **125 my**. To load sands open the transparent door, lift the grid **pos 118 pict 2**, pour sand into the tank **pos 140** (never exceed the first notch of the tank **pos 140 SAND LEVEL.**); re-place the grid.

Among accessories you can find: plug **pos 133 pict 2** and sieve **pos 134** to be fitted into the sand tank.

The plug has to be fitted to prevent sand recycling; for instance in case of roughing operations when uncontaminated sand is used.

The sieve holds the bigger sand particles to avoid ducts obstruction; for instance, during the investment removal, big sized fragments could be produced.

Put an object into the cabin, close the door, turn on the internal lamp and the aspirator. Adjust suction through the knob **pos 3 pict1**. The correct level of suction is indicated by gloves: when they are pushed inside the cabin without deformation of their volume (they do not inflate) they indicate the suction is correct.

Adjust the tool pressure through the reducer **pos 51 pict3** to the following values:

**CORUNDUM tool 6-6.5 Bars**

**BIOXIDE and BALLS tool, 3.5-4.5 Bars.**

Fit your hands into the gloves, hold the tool by the right hand and support the object by the left hand, press **pict5** and aim the jet to the surface to be processed.

The strongest abrasive force of the sand is near the nozzle, the more you move away from it the slower the speed and the bigger the jet width. see **pict6**.

When the processing is over, before turning off the aspirator, unlock the door and slightly move the glass away from the seal.

In few seconds the aspirator can remove dust completely from the cabin and at door opening no dust will disperse in the environment.

Before every sandblasting intervention, it is important to check the filtering bag efficiency.

It is impossible to decide an amount of working hours for filter replacement since it

depends on the different working conditions.

Indeed, some types of processing generate a lot of dust, for instance the investment removal. The use of very fine sands, or recycled sands, wear the filter.

Also the simultaneous use of cabins reduces the filtering bag efficiency more rapidly. The Operator has to decide on the filter replacement according to the carried out processes.

In any case, we suggest to replace the bag quite often since the machine operation and the operators' health depend on its correct action.

The cabin cleaning is carried out by means of a sucking tube **pos 252 pict 1**, as follows:

Remove the glove protection valve **pos 175 pict 1** and fit the tube, turn on the turbine and take it to the highest speed by means of the regulator **pos 3**, and with the nozzle **pos 251** suck sand from the cabin.

**AFTER CLEANING THE CABIN, THE BAG HAS TO BE REPLACED AND THE ACCESSORIES CLEANED TOO.**

The door is covered by a transparent film (replaceable) to protect it against abrasion.

## EMISSIONS OF TOXIC SUBSTANCES DURING OPERATION

Materials used for sandblasting, are not harmful for the skin but can be very dangerous in case of inhalation

Sandblasting generates very fine dusts, the aspirator has to convey the cabin air, saturated with such particles, to a container (filter) made of porous material (paper). The paper filter pores have a smaller size than that of dust particles, the latter are hence held while the air passes through to reach a second filter with characteristics specially conceived for health protection; the second filter cleans the air from the filtering bag once more before air can reach the environment.

The above described measures to protect people's health are effective if instructions are complied with, **routine maintenance and extraordinary maintenance**.

If filters are performing correctly, dust emissions during the machine use are negligible.

**TECNO GAZ is not responsible for the bad functioning of the cleaning circuit in case it is due to the non respect of the relevant instructions or to the unskilfulness of the machine user.**

The wrong functioning of the air cleaning system due to defects of the original components has to be notified to the local retailer in due time.

### ATTENTION:

DO NOT USE THE MACHINE WITH DOORS OPEN.

DO NOT USE THE MACHINE WITHOUT GLOVES.

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Routine maintenance concerns the following components (**PICT1**):

- 1) FILTERING BAG. **pos 253**
- 2) FINAL FILTER PROTECTION **pos 114**
- 3) FINAL FILTER **pos 230**
- 4) GLOVES **pos 254**
- 5) GLASS PROTECTION **pos 229**

**1) FILTERING BAG**

The filtering bag has to be replaced frequently; if bags are not replaced frequently you might cause the obstruction of the machine ducts and put at risk the operators' health.

It is not possible to quantify the exact duration of a filter since it depends on many factors.

The operator has to decide on the filter replacement according to the type of processing: consider that when the filter is close to saturation visibility inside the cabin is reduced. In this case, **do not increase the turbine revolutions but replace the filter.**

**ATTENTION:**

**After every cabin cleaning the filtering bag has to be replaced.**

The saturated filter (obstructed) does not stop dust anymore and that means a very heavy work for the final filter.

**The "ONE WAY" filtering bag cannot be cleaned nor reused.**

To replace the filtering bag remove the plug **pos 141 pict 1** (rotate to the right or to the left while pulling it upward); hanging from the plug is the filter holder.

Take the assembly in a dry place where dust cannot be harmful and wear a protective nose- and-mouth mask.

Take the plug by one hand and the holder by the other **pict 7**.

Push the holder to the plug until reaching the end-of-stroke, then rotate the holder counterclockwise to the end-of-stroke, and separate slowly the holder from the plug.

When the filtering bag neck has come out of the holder, seize it with the fingers of your left hand and extract it completely (be sure the sand contained in the bag does not tear it).

Dispose of the full bag in a suitable dump and brush the holder and the plug.

To re-assemble the filter, open a new bag, fit the plug lower side and wrap the bag around the plug **pict8**.

Fold the exceeding part of the bag clockwise on the plug.

Fit the bag into the holder without curling it up; pins along the plug diameter have to pierce the filter paper and be fitted into the slots on the upper edge of the

holder.

Fix the holder on the plug, rotate it clockwise until reaching the end-of-stroke and then lock it by pulling it to the plug.

Before fitting the filter into the machine, clean the induction manifold with an aspirator **pos 142 pict 1**.

To facilitate the filter introduction cover seals **pos 282** with talc.

## 2) FINAL FILTER PROTECTION

During the extraction of the filtering bag some dust located between the induction manifold **pos 142 pict1** and the plug **pos 141 pict1**, falls into the suction chamber.

To protect both the turbine and the final filter from possible obstructions has been equipped with a fine texture fabric bag holding the dust.

Such component is not a filter however it requires the same maintenance interventions.

**The fabric bag has to be cleaned each time the filtering bag is replaced: remove the fabric bag contents, shake it open air or vacuum clean the fabric.**

This component has to be replaced in case of tear or wear.

The fabric bag inside the suction chamber is not secured but simply pressure fixed.

In the upper edge the of the fabric a round spring (open ring) is positioned inside a tape **pict 9**.

To extract the fabric bag introduce your hand into the bag and by your fingers extract the bag upper edge by pulling it upward **pict 9**.

When the bag is mounted be sure the edge and hence the spring leans on the three projecting screws inside the chamber.

## 3) FINAL FILTER (pict 10)

The final filter holds the dust escaped from previous treatments and protects people's health.

The final filter has to be serviced periodically; approximately every ten replacements of the filtering bag the final filter condition has to be checked.

The final filter can be cleaned and reused; through the gun blow the air from the outside to the inside.

**ATTENTION: THE GUN AIR PRESSURE HAS TO BE LOW (LOWER THAN 1 BAR) AND THE JET HAS TO BE KEPT AWAY FROM THE PAPER NOT TO TEAR IT.**

The performance of the above-mentioned filter is limited and depends on the performance of previous filters, the filter can be re-used only for a reasonable amount of times.

**It is not possible to assess the filter condition visually, hence it is good rule to replace it to avoid environment contamination.**

When intervening on this filter, the filter box **pos 122-123**, the turbine duct and the protection sponge **pos 113** have to be cleaned; if possible, use an aspirator.

To dismount the FINAL filter, UNPLUG THE MACHINE FROM POWER SUPPLY, remove the side door **pos 105** with the 'filter' marking, unscrew the screws located next to the filter box **pos 320** and extract the cover.

**When the filter is fitted into the box, respect the seal position and do not forget to fit the protection sponge into the filter frame. See pict 10.**

#### 4) GLOVES (PICT 12)

To replace gloves, unscrew by 1/4 of turn ccw the glove locking ring nut **pos 144** positioned inside the cabin, then remove and extract the glove.

To mount gloves, fit them into the bushing **pos 124** and fit the ring nut into the joint and rotate cw by 1/4 of turn. See **pict 11**.

Thanks to the special system, any type of gloves can be fitted.

#### 5) GLASS PROTECTION (PICT 12)

The transparent film beneath the glass is fixed with glue (supplied with the machine) on the external perimeter.

To remove the film lift an edge and pull it.

Clean the glass with alcohol, apply the glue to the new film, apply on the hinges holes too and be sure the film perfectly adheres to the glass.

### EXTRAORDINARY MAINTENANCE AND REPAIR

Extraordinary maintenance: periodical interventions to check and clean the machine internal parts.

Such interventions are performed on-site by authorised staff or at TECNO GAZ premises.

Repair interventions too can be performed on-site by authorised staff.

Apply to the nearest retailer to ask for technical service.

### INSTRUCTIONS FOR MACHINE DEMOLITION

The sandblasting machine is made of iron material and inert electromechanic components.

At the machine demolition no special instruction has to be complied with since all materials are reusable.

Do not dispose of the machine in unattended areas, apply to authorised dumps.

**LIST OF ELECTRIC SYSTEM COMPONENTS**

Pos.	Q.ty	Description	Cod.
1	3	Bipolar switch, luminous, green, 250 V 10 A	CEQ0001
2	3	Bipolar switch, luminous, green, 250 V 10 A	CEQ0001
4	1	Turbine regulation control unit 220 V 50 HZ with potentiometer	SSKA004
5	1	Motor with turbine 220V 50HZ with brushes 0.8 KW	CPG0001
6	3	Lamp E14 15Watt 230 V	CEE0005
7	3	Colourless lamp holder	CEE0004
8	1	Double fuses holder for min. 5x20 fuses	CEF0001
9	2	Fuses min. 5x20 5A	CEF0009
30	1	Passband filter	CEZ0002
99	1	Electropilot	CPY0027

**LIST OF COMPONENTS OF THE PNEUMATIC SYSTEM**


Pos	Q.ty	Description	Cod.
51	3	Pressure reducer 1/8" with connections	SSKA071
52	3	Manometer diam. 53 1/8" flange and bracket	CPP0001
53	3	Air gun	CP02001
85	1	Safety kit	SSKA085
99	1	Electropilot	CPY0027
100	1	3-way valve	CPY0028
145	1	Base	SSKA145

**SUMMARY OF INSTRUCTIONS TO COMPLY WITH**












- 1) To lift the machine, fit your hands in the gloves holes.
- 2) Has to be fed with DRY AND NON LUBRIFIED compressed air.
- 3) When sand is loaded in the cabin, comply with indications concerning the type and the amount of sand.
- 4) Do not introduce wet objects into the cabin.
- 5) DO NOT USE THE MACHINE IF DOORS ARE OPEN
- 6) Always use the aspirator at the set speed.
- 7) Be sure the filtering bag is always present.
- 8) Replace the filtering bag at frequent intervals.**
- 9) Always replace the filtering bag at each cabin cleaning.**
- 10) At each replacing of the filtering bag: clean the plug, the holder, and the fabric bag with seals.
- 11) Cover seals with talc.
- 12) Be careful not to tear the filtering bag.
- 13) Periodically check the FINAL filter, clean and replace it.
- 14) Do not forget to apply the transparent film to the glass.
- 15) The feeding pressure of the big tool 3.5 diam. is of 6-6.5 bars
- 16) The feeding pressure of the big tool 1.5 diam. is of 4-5 bars
- 17) DO NOT USE THE MACHINE WITHOUT GLOVES.
- 18) DO NOT START THE TOOL WHEN THE DOOR IS OPEN, THE SAND JET CAN INJURE THE OPERATOR OR NEARBY PEOPLE.
- 19) WHEN THE FINAL FILTER CONDITION IS CHECKED, OR IN CASE OF ANY OTHER INTERVENTION INSIDE THE MACHINE, UNPLUG THE MACHINE FROM POWER SUPPLY.
- 20) In case of interventions on tools with doors open, disconnect the pneumatic feeding and discharge the machine system (DETACH THE TUBE).**

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**TECNO-GAZ S.p.A.**