Istruzioni per l'uso
Instructions for use
Mode d'emploi
Bedienungsanleitung
Istrucciones para el uso

ASTRA S

- Cod. 2190-U (3 Lt)
- Cod. 2192-U (6 Lt)
- Cod. 2182-U (9.5 Lt)
- Cod. 2183-U (14 Lt)
- Cod. 2184-U (18 Lt)
- Cod. 2185-U (21 Lt)

ASTRA D

- Cod. 2191-U (3 Lt)
- Cod. 2193-U (6 Lt)
- Cod. 2186-U (9.5 Lt)
- Cod. 2187-U (14 Lt)
- Cod. 2188-U (18 Lt)
- Cod. 2189-U (21 Lt)

I

GB

F

D

E

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Abstract:

The efficiency of a decontamination procedure by sonication for different dental instruments after experimental microbial and viral contamination was tested. Both germicidal and virucidal activity of sonication in presence or absence of a cationic bio-biguanid disinfectant was assessed following three different disinfection/sterilisation protocols. Dental instruments were contaminated with a mixed culture of Streptococcus faecium, Staphylococcus sp., Pseudomonas aeruginosa, Mycobacterium sp., Escherichia coli and Bacillus subtilis, or with Polio virus type 1 and Herpesvirus simplex type 1 (HHV1, following the new herpesvirus denomination), exposed to ultrasonic treatment in an ultrasonic bath and the surviving microorganisms titered. The results showed that an effective disinfection of dental instruments, expressed by an equal or higher than 4 logs microbial and viral reduction, can be obtained after 15 min or 10 min sonication in the ultrasonic cleaner equipped with a Sweep System Technology. Conversely, by the combined action of chemical disinfection and ultrasonic treatment in the same device, a sterilising effect was obtained after only 5 min for microbial and after 10-15 min for virally contaminated instruments. The synergistic effect of chemical and physical means, as already accepted as an effective cleaning procedure of medical instruments, can therefore be applied as preamble to the final autoclave sterilization process of dental instruments (surgical) potentially contaminated by organic fluids and dental material harbouring pathogenic microbes and viruses.

AUTHORS: *Patrizia Di Gennaro, _Giuseppina Bestetti, °Antonia Radaelli, §Manuela Paganini, §Carlo De Giuli Morghen, §Margherita Neri

*Department of Environmental Sciences, University of Milano-Bicocca, 20126 Milano, Italy. §Departments of Medical Pharmacology and °Pharmacological Sciences, University of Milano, 20129 Milano, Italy.

_Corresponding author: Prof. Giuseppina Bestetti, Department of Environmental Sciences, University of Milano-Bicocca, P.zza della Scienza, 1-20126 Milano Italy, Phone: 0039-2-64482925; Fax: 0039-2-64482996; E-mail: Giuseppina.Bestetti@unimib.it

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The ASTRA ultrasonic units allow various types of objects to be cleaned safely and easily. The small size of these units means that they can be placed directly on the work surface and fulfil various requirements.

General features of the units

- Simple and easy-to-read controls.
- High working frequency.
- Low noise levels.
- High power ultrasound generators.
- Heating temperature and cleaning time control by microprocessor (only for D series)


Fig. 1


Fig. 2


Fig. 3


Fig. 4

Warranty

The information contained in this manual may be changed without prior notice.

TECNO GAZ S.p.A. declines responsibility for errors in the manual, for accidental damage or that caused by the supply, performances or use of this equipment. This product is guaranteed against material and manufacturing defects for 12 (twelve) months from the date of purchase of the product according to the clauses of the warranty certificate provided.

In the event of a fault occurring in the unit during the warranty period, the firm TECNO GAZ S.p.A. may decide whether to repair or replace the product found to be defective.

Work during the warranty period

For work or repairs under guarantee, this unit must be sent to TECNO GAZ S.p.A. via the retailer.

The warranty does not cover transfer costs of technical staff, nor costs of despatch or transport risks, which are to be borne by the purchaser. The purchaser shall in all cases bear the despatch costs, including duties, of products sent to TECNO GAZ S.p.A. from another country.

Warranty restrictions

The warranty covers replacement or repair of the components recognised as being unsuitable due to manufacturing defects, inclusive of the necessary labour.

The above is not applicable to faults caused by improper or inadequate maintenance or use by the customer, by unauthorised changes made to the unit, by the use of the product in environmental conditions other than those specified in this manual for said product or by inadequate preparation of the installation area. No compensation shall be

payable during the period of time required for replacement or repair of the unit.

Replacement shall in any case be decided solely by the manufacturer and only in the event of the unit being ascertained as wholly unsuitable and impossible to repair.

Compensation of direct or indirect damages of whatsoever kind to persons or property due to the use or non-use of the unit shall not be due.

Safety

In order to retain the original safety features of the product, the customer must not replace any part thereof nor make unauthorised changes.

WARNING

This symbol draws attention to one or more procedures whose partial or complete non-observance may cause partial or total damage to the product or injure the operator.

Before carrying out the procedures indicated after this symbol, make sure the conditions specified have been fully understood and observed.

Safety symbols:

- Symbol of the instructions manual. Consult the manual to avoid any damage to the product or injury to the operator.

- High voltage symbol. Do not open or tamper with the device. Live electric parts inside the unit.

DECLARATION OF CONFORMITY

We, the firm
 TECNO-GAZ S.p.A.
 Strada Cavalli, 4
 43038 Sala Baganza (Parma)

declare under our responsibility that the ultrasonic units of the **ASTRA 3S - ASTRA 6S - ASTRA 9.5S - ASTRA14S - ASTRA 18S - ASTRA 21S - ASTRA 3D - ASTRA 6D - ASTRA 9.5D - ASTRA14D - ASTRA 18D - ASTRA 21D** where to this declaration pertains, conform with the standards and directives indicated below:

EN 55011B;

EN 61010

according to the requirements of the directives 89/336/EEC, 73/23/EEC and 93/42/EEC.

Data 15-11-2005

Managing Director
 Tecno Gaz S.p.A.

TECHNICAL FEATURES ASTRA ULTRASONIC UNITS

Model	ASTRA 3S	ASTRA 6S	ASTRA 9.5S	ASTRA 14S	ASTRA 18S	ASTRA 21S
Supply Voltage	230/240V - 50/60Hz					
Input Power	130 W	180 W	200 W	300 W	300 W	400 W
Input Power with heating	305 W	355 W	500 W	800 W	800 W	1400 W
Weight in Kg	2,8	3,5	6,3	7,1	7,5	14,5
EXTERNAL Dimensions (in mm)	270 170 210	325 270 210	400 270 370	440 340 425	440 340 425	600 330 425
TANK Dimensions (in mm)	240 140 100	300 240 100	300 240 100	330 300 150	330 300 200	500 300 150
Number of transducers	2	4	3	4	4	6
Ambient conditions	Temperature from 5 to 40°C, rel. hum. 80% up to 31°C with linear decrease up to 50% at 40°C					
Installation category	Class II according to EN 61010					

Model	ASTRA 3D	ASTRA 6D	ASTRA 9.5D	ASTRA14D	ASTRA 18D	ASTRA 21D
Supply Voltage	230/240V - 50/60Hz					
Input Power	130 W	180 W	200 W	300 W	300 W	400 W
Input Power with heating	305 W	355 W	500 W	800 W	800 W	1400 W
Weight in Kg	2,8	3,5	6,3	7,1	7,5	14,5
EXTERNAL Dimensions (in mm)	270 170 210	325 270 210	400 270 370	440 340 425	440 340 425	600 330 425
TANK Dimensions (in mm)	240 140 100	300 240 100	300 240 150	330 300 150	330 300 200	500 300 150
Number of transducers	2	4	3	4	4	6
Ambient conditions	Temperature from 5 to 40°C, rel. hum. 80% up to 31°C with linear decrease up to 50% at 40°C					
Installation category	Class II according to EN 61010					

INDEX:

Chapter 1

Introduction	15
Electrical connection	15
Before starting up the unit	15

Chapter 2

Control panel and description of symbols	16
--	----

Chapter 3

Working instructions for ASTRA models in the S and D series	17
---	----

Chapter 4

Accessories for the ASTRA units	18
---------------------------------	----

Chapter 5

Cleaning with ultrasounds	18
Useful advice and suggestions	19
SONICA® liquid detergents	19
Maintenance	20

Chapter 6

After-sales service	20
---------------------	----

Chapter 1

INTRODUCTION

Dear customer,

We wish first of all to thank you for having chosen a ASTRA ultrasonic unit. All ASTRA ultrasonic cleaning units are totally reliable and every component has been designed and produced to ensure top performances at all times.

Please complete the Warranty Certificate, enclosed with the instructions manual, carefully and throughout and send it off straightaway.

In this way you can benefit from a twelve-month warranty as from the purchase date, according to the terms stated in the certificate. This instructions manual is an essential part of ASTRA equipment. It contains important information about installation and operation security, maintenance and use. This instructions manual should always accompany the ASTRA equipment, an erroneous installation or use may cause damages to people, animals or things. TECNO GAZ S.p.A. may not be considered responsible for damages caused by improper use or installation.

ELECTRICAL CONNECTION

Before plugging in the unit to the power socket, check that the voltage indicated on the rating plate, attached to the rear of the unit (⑫), corresponds to that of your work place. (See Fig. 4).

Earthing of the unit is compulsory by law. Any deliberate cutting of the earth wire inside or outside the unit or removal of the earth terminal from the plug is forbidden in that it makes use of the unit itself hazardous.

The manufacturer declines all responsibility in respect of persons or things arising from failure to observe this rule. The electrical safety of this unit is only guaranteed when it is correctly connected to an effective earthing system, as provided by current electrical safety laws.

If you do not have an electrical system with proper earthing, do not connect the unit to the socket and consult a specialist electrician as soon as possible.

WARNING

This unit requires earthing.

BEFORE STARTING UP THE UNIT

WARNING

Pour the washing liquid into the tank as far as the horizontal reference notch printed inside the tank. Make sure that the level of the liquid is always within this level.

Make sure that the electric socket has been disconnected before filling up the tank. Any accidental leakage of liquid could cause a short circuit or electrocution and would therefore be dangerous for the operator.

Check that the unit has not been damaged. Do not use units which have been damaged during transport; if in doubt consult the retailer or the manufacturer directly.

Install the unit on a flat and stable surface, capable of bearing the weight of the unit, the tools or objects which are placed inside for washing and the relevant liquid. Handle it with care.

Install the ultrasonic cleaning unit away from sources of heat. Also ensure that it is not installed near sources of moisture, on wet work surfaces or near sources of dust.

Make sure that the feet of the unit are correctly positioned to ensure circulation of air.

Do not start up the unit if: the power cable or plug has been damaged, the unit does not operate properly, it has been damaged or if it has fallen.

In these circumstances there is the risk of electric shock, fire or other accidents. Do not attempt to adjust the unit personally, but instead contact only technical specialists or the manufacturer.

Do not immerse the power supply cable or plug in water. Keep the power supply cable away from hot surfaces.

Do not allow the supply cable to hang from the edges of tables or furniture units.

Do not leave or use this unit outdoors.

Do not under any circumstances tamper with the electronic system of the ultrasonic cleaning unit as there is the risk of electric shock. For any repairs always contact the after-sales service of the manufacturer.

Do not replace the power supply cable. If the cable is damaged through wear or other reasons, switch off the unit immediately and contact the assistance service.

Do not lift and carry the unit when it is full of liquid. The handles have been designed exclusively for lifting and carrying the unit only when it is empty and disconnected from the supply socket.

WARNING

This ultrasonic cleaning unit only operates with water or detergents for ultrasonic treatments recommended by the manufacturer.

Only pour the quantity of liquid required for washing. **Do not use acid or highly alkaline** solutions such as sodium hypochlorite as these substances, coming directly into contact with the steel tank, cause a series of microscopic holes with irreparable and above all hazardous consequences for the operation of your unit. All acid substances or those alkaline substances which may release corrosive compounds such as chloride or other chemical compounds, or disinfecting substances with a glutaraldehyde base, used together with strong ultrasonic cavitation, cause highly accelerated corrosion even in stainless steel.

If you must necessarily use substances which could damage the steel tank, use glass or plastic beakers. Remember that hydrofluoric acid breaks glass beakers and that its use is highly dangerous for your health.

Remember that any acid particles, micronised by intense ultrasonic cavitation, are dispersed in the work area, **corroding your equipment and damaging your health too.**

Do not use substances such as petrol, benzol or benzene, or other harmful or explosive or flammable solvents. Only use solutions suitable for the type of work to be performed.

To prevent damage to the ultrasonic cleaner, change the solution regularly, don't operate the cleaner dry, don't place parts or containers directly on the bottom of the cleaning tank; use a tray or wire to suspend items. Failure to comply may cause transducers damage and will void your warranty.

Don't allow the solution to drop more than 1 cm below the operating level line with heat or ultrasonics on. Failure to comply may cause transducers and/or heater damage and will void your warranty.

Important instructions for use of the ultrasonic cleaning unit

Before switching on the unit make sure the steel tank has been filled sufficiently with water or liquid detergent and, if the version of your unit is equipped with a liquids drain device, check that the tap is fully turned off. In this way accidental leakage of liquid substances will be avoided.

Always check the level of liquid so that it does not drop below 4 centimetres from the upper edge of the tank. If this should happen, the heating element and ultrasonic transducers and the electronic circuit could suffer serious damage.

Do not place your hands in the tank during operation.

Aggressive liquids which could damage the tank should be placed in a special beaker which is then placed in the tank containing water which transmits the ultrasonic waves to the beaker and hence to the solution it contains. Remember not to rest the glass beaker on the bottom of the tank; a special perforated support for these glass containers is available.

Only use this unit for the professional use foreseen in this manual and for the purpose for which it has been designed. This unit has been designed for washing and detaching surplus materials from instruments, prostheses, drills, probes, forceps, metal objects and metal parts in general, objects and test-tubes in glass, crystals, stones and archaeological objects, rings, spectacles, parts which are difficult to reach by hand etc.

The ASTRA cleaning units also provide ideal solutions for emulsion, accelerating chemical reactions, mixing solutions, degassing liquids, decontamination and for dissolving sediments.

Chapter 2

CONTROL PANEL AND DESCRIPTION OF SYMBOLS

Fig. 1 Front panel of ASTRA models in the S series

- ① Heating on/off switch
- ② Timer/Switch knob
- ③ Cleaning time in minutes
- ④ Symbol  (infinite cleaning time)

Fig. 2 Front panel of the ASTRA models in the D series

- ⑤ Main Switch
- ⑥ Heating resistor "On/Off" red pilot light .
- ⑦ Dip temperature programming buttons, 40°-50°-60°C
- ⑧ Heating  button
- ⑨ Ultrasounds  button.
- Cleaning times programming button, 5-10-15' min.
- Unit "On" green pilot light

 Resistance symbol

 Cleaning time symbol (in minutes)

 Cleaning temperature symbol (degrees centigrade)

Fig. 4 Rear side view of the ASTRA units

-  Unit rating plate.
-  Socket.

Chapter 3

WORKING INSTRUCTIONS FOR MODELS IN THE S SERIES (see Fig.1)

- 1) To switch on the unit rotate the Timer/Switch knob clockwise (2) and set the required cleaning time.
- 2) At the end of the time set the machine will switch off automatically. If operation of the machine is to be interrupted at any time, return the Timer/Switch knob to the "0" setting.
- 3) In order to set a cleaning time longer than 15 minutes, rotate the Timer/Switch knob anticlockwise. In this way the cleaning time can be set by hand, and the unit switched on or off according to needs (see Fig. 3).
- 4) In order to actuate heating with thermostat control at 60°C, press the illuminated button marked Heater (1)  60°C
- 5) After the washing cycle, in order to de-actuate heating, press the Heater switch (1)  60°C

WORKING INSTRUCTIONS FOR MODELS IN THE D SERIES (see Fig. 2)

- 1) To switch the unit on or off press the main switch (5).
- 2) After the main switch has been pressed, the green ON pilot light (11) and that of the 40°C temperature button will light up. The 40°C temperature button pilot light will light up each time the main (5) switch is pressed as the microprocessor always sets the lowest temperature and the shortest cleaning time (in any case, in order to start the cleaning cycle or switch on heating, the (8) "  " and (9) "  " buttons must be pressed.
- 3) Press one of the buttons for 5-10-15 minutes according to needs (10).
- 4) In order to start or stop the cleaning cycle press the "  " button (9). In this way the cleaning cycle will be started for the time you have set. At the end of this time the machine will stop automatically.
- 5) To double washing times and to obtain washing timing of 10-20-30 minutes, press the minutes key (10) twice corresponding to the time you require (within a maximum interval of 5 seconds). The machine will give a double Beep to confirm that double pressing of key (10) has been correctly set. Press the " Start/Stop " (9) "  " key to start or stop the washing cycle. Note: If you have set double washing times, the machine will give a double Beep after pressing the key (9)  before the washing cycle starts.

HOW TO REGULATE THE CLEANING TEMPERATURE (see Fig. 2)

- 1) Before switching on the heating check that there is liquid inside the tank and that the level is not below 4 cm from the upper edge. If the heating is switched on without liquid, there is the risk of burning the resistor and damaging the entire unit.
- 2) Set the required temperature of the washing dip by pressing one of the three buttons: 40°-50°-60°C (7). On pressing one of the three buttons, the pilot light corresponding to the temperature button set will light up.
- 3) In order to switch on the resistor press the heating (6) "  " button. The corresponding pilot light will light up.
- 4) On the control panel the red heating pilot light (6) marked by the symbol  will also light up. If the pilot light is on, this means that the resistor is in operation and that the liquid is heating. When the pilot light goes out this means that the resistor is no longer in operation since the liquid has reached the temperature set previously.

Note:

Whenever a button of the membrane pushbutton panel of the control panel is pressed, a beep is heard which confirms pressing of the button.

WARNING

When the heating function is enabled with temperatures higher than 50°C, do not place your hands in the washing liquid. Risk of scolding and burns. Use gloves and adequate means of protection.

Chapter 4

ACCESSORIES FOR ASTRA ULTRASONIC UNITS

The following accessories are available for the ultrasonic units on request:

- **Stainless steel lid.**

Use the steel lid to close the ultrasonic tank.

- **Take great care** with the condensation produced through the heating steam inside the lid since, during opening of the latter, condensation liquid could drip onto the electrical parts such as the socket. The lid therefore has to be raised vertically allowing the condensation to drip only inside the ultrasounds tank.

- **Rectangular basket in stainless steel mesh.**

The stainless steel mesh basket is a very useful accessory in that it allows washing of the objects to be optimised and the objects themselves to be drained inside the tank after washing. By folding the handles, the basket can be suspended with precision over the tank.

- **Glass beakers.**



Special beakers with an adjustable black rubber support ring are available. They are used to save detergent liquid whenever small quantities of objects are to be cleaned, or if washing is to be performed with a detergent solution in a beaker and rinsing simultaneously in another beaker. They are also used to wash objects with particularly aggressive solutions which must not come into direct contact with the steel tank.

- **Stainless steel beaker-holder.**

The perforated lid is used to support the beakers during washing. Remember to adjust the support ring of the beaker so that the bottom of the beaker is always in contact with the liquid of the tank, but at the same time does not come into contact with the bottom of the steel basket.

- **Small round basket**

In stainless steel mesh to be inserted in the beaker for washing small parts.

Chapter 5

CLEANING WITH ULTRASOUNDS

What are ultrasounds?

Ultrasounds are vibrations of a material medium which are similar to sound waves, but with frequencies which are too high to be audible to the human ear. The frequency limit of human hearing ranges from around 10 KHz to approximately 18 KHz. The audible frequency limit decreases with the increase in age. More particularly ASTRA ultrasonic cleaning units operate at a frequency of as many as 45 KHz, i.e. at a frequency much higher than the audible one. All the ASTRA units consist of an ultrasounds generator and one or more piezoelectric transducers (according to the model) which are attached to the outside base of the tank in stainless steel containing the washing liquid. The electronic generator produces a continuous signal at a frequency of 45 KHz, and pilots the piezoelectric transducers which transform the electrical signal into a mechanical vibration. This vibrational energy is transmitted to the liquid in the tank at the frequency of 45,000 oscillations per second. These pressure and vacuum oscillations create an enormous quantity of microbubbles in the liquid which, by imploding in extremely rapid succession, create enormous impact energies between the detergent liquid and the surface to be cleaned. This is known as "cavitation" and provides an efficient and safe cleaning method in less time.

USEFUL ADVICE AND SUGGESTIONS

The type of detergent to be used. In order to wash any instrument or object, use tap water as the main substance and dissolve in it a quantity of detergent or solution which varies according to the type of material residue to be removed from the object. Only use pure detergent if expressly specified in the instructions for use of the product. If for example an object soiled with greasy or fatty material is to be cleaned, use a small quantity of SONICA® UG solution, diluting it in water as indicated in the instructions for the product and clean in ultrasounds.

The SONICA® UG detergent enables removal of many types of impurities.

When small quantities of plaster have to be removed. Use pure SONICA® RG solution. Pour it into a beaker or directly into the tank. Its strong action will allow total removal of the plaster even from dental prostheses.

When many objects have to be cleaned. Do not overload the ultrasonic cleaning unit. Preferably arrange the objects rationally on the base of the rectangular basket, place the latter in the tank and clean as appropriate. Too many objects washed at the same time reduce the efficiency of cleaning with ultrasounds.

When the objects are very dirty. Use long cleaning times, heating the washing liquid according to the object to be washed.

When to replace the detergent solution. To achieve efficient cleaning change the detergent solution often, above all if it is used to clean surgical instruments.

Shape and size. There are no special rules for the shape of objects to be cleaned, however avoid inserting particularly heavy and large objects in the tank.

The liquid must always cover the entire object.

Cleaning time. The cleaning time always depends on the type of object, and the type and quantity of the material to be removed.

SONICA® LIQUID DETERGENTS

SONICA® UG Concentrated detergent for general use. Suitable for cleaning objects and instruments made of plastic, metal, rubber and silicone, dental instruments, for getting greasy deposits off precision mechanical parts and for cleaning metallographic coupons. Also suitable for cleaning lenses, spectacle frames, electrostatic metal filters and gas-masks. Does not contain dangerous substances. Dilution ratio: 1:1 to 1:10 depending on strength needed. pH 6.5 to 7.5

SONICA® RC Ready-to-use detergent for removing cement deposits from dental instruments. Not suitable for use with aluminium. pH 13

SONICA® AF Concentrated detergent for removing limestone deposit stains and large quantities of oxides and rust. Not suitable for aluminium. Suitable for cleaning very dirty limestone-encrusted laboratory glassware. Dilution ratio: 1:10 to 1:20 depending on thickness of deposits. Acidic product. pH 0

SONICA® AC Concentrated detergent for removing limestone traces and small amounts of oxides and rust. Suitable for removing tartar from dental prostheses. Dilution ratio: 1:5. pH 2.0-2.5

SONICA® RA Ready-to-use detergent for removing alginate deposits from impression trays without damage even after lengthy immersion. The detergent is also suitable for metal impression trays. pH 7.5-8.5

SONICA® RG Special neutral detergent for removing plaster deposits from all types of surfaces without damage e.g. articulator plates, spatulas, miscellaneous instruments, all types of prostheses. Can also be used for gypsum investments. When plaster comes into contact with the cleaning liquid small bubbles form and show that the detergent is acting on the plaster. The product should not be diluted. pH 7.5-8.5

SONICA® CL 4% is a concentrated aqueous solution with decontaminating and cleansing action developed for using with SONICA® Ultrasonic cleaners. The solution is concentrated.

Application Fields

Decontamination and simultaneous detergency of surgical instruments and medical devices. Temporary sterile conservation of surgical instruments.

Action Mechanism

Chlorexidine reacts with the negatively charged groups that are found on the cell surface. For this reason, when it comes into contact with bacterial suspensions, it is absorbed immediately, causing an irreversible loss of cytoplasm components, which damages the cell membrane and inhibits enzymes. Cetrimide is a quaternary ammonium salt that reduces surface tension in the point of contact and has precipitant, complexing and denaturing effects on bacteria proteins, causing enzyme changes in the cytoplasm membrane.

Germicidal Activity

Chlorexidine acts on vegetative bacteria, yeasts, fungi, several protozoa, viruses (HIV) and salmonella. Cetrimide acts as a bactericide above all on gram-positive bacteria; it has varying fungistat properties and acts as a virucide against lipophilic viruses. The action of the product is superior to that of its single components. Its bactericidal actions are reduced by the presence of organic material (serum).

In hospital usage SONICA® CL 4% is mainly recommended for its action against gram-negative and gram-positive bacteria, as well as lipophilic viruses (HIV).

From the tests made, SONICA® CL 4% has shown that it acts against such bacteria as *Citrobacter freundii*, *Serratia marcescens*, *Pseudomonas aeruginosa* ATCC 15442, *Pseudomonas aeruginosa* ATCC 14502, *Proteus rettgeri*, *Staphylococcus aureus* ATCC 114, *Staphylococcus aureus* ATCC 6538, and *Klebsiella pneumoniae*.

Composition

100 g of solution contain: Active Ingredients: Chlorhexidine gluconate g 1,50; Cetrimide g 15,00; Excipients: Isopropyl alcohol (F.U.I.) g 6,00; Essence g 2,00; Colouring E110 g 0,10; Purified water q.b. to g 100,00

Abstract:

The efficiency of a decontamination procedure by sonication for different dental instruments after experimental microbial and viral contamination was tested.

Both germicidal and virucidal activity of sonication in presence or absence of a cationic bio-biguanid disinfectant was assessed following three different disinfection/sterilisation protocols. Dental instruments were contaminated with a mixed culture of *Streptococcus faecium*, *Staphylococcus* sp., *Pseudomonas aeruginosa*, *Mycobacterium* sp., *Escherichia coli* and *Bacillus subtilis*, or with Polio virus type 1 and Herpesvirus simplex type 1 (HHV1, following the new herpesvirus denomination), exposed to ultrasonic treatment in an ultrasonic bath and the surviving microorganisms titered. The results showed that an effective disinfection of dental instruments, expressed by an equal or higher than 4 logs microbial and viral reduction, can be obtained after 15 min or 10 min sonication in the ultrasonic cleaner equipped with a Sweep System Technology. Conversely, by the combined action of chemical disinfection and ultrasonic treatment in the same device, a sterilising effect was obtained after only 5 min for microbial and after 10-15 min for virally contaminated instruments. The synergistic effect of chemical and physical means, as already accepted as an effective cleaning procedure of medical instruments, can therefore be applied to obtain a safe and effective sterilisation of dental instruments potentially contaminated by organic fluids and dental material harbouring pathogenic microbes and viruses.

AUTHORS: *Patrizia Di Gennaro, _Giuseppina Bestetti, °Antonia Radaelli, §Manuela Paganini, §Carlo De Giuli Morghen, §Margherita Neri

*Department of Environmental Sciences, University of Milano-Bicocca, 20126 Milano, Italy. §Departments of Medical Pharmacology and °Pharmacological Sciences, University of Milano, 20129 Milano, Italy.

_Corresponding author: Prof. Giuseppina Bestetti, Department of Environmental Sciences, University of Milano-Bicocca, P.zza della Scienza, 1-20126 Milano Italy, Phone: 0039-2-64482925; Fax: 0039-2-64482996; E-mail: Giuseppina.Bestetti@unimib.it

MAINTENANCE OF THE UNIT

Cleaning is the only maintenance operation normally required. It must be performed with the unit switched off in electrical terms. Only use mild detergents for cleaning the inside of the tank and a soft cloth for cleaning the outside parts such as the control panel and casing. **Do not leave deposits of dirt** inside the tank above all if your model is fitted with a liquids drain device. If the liquid outlet hole is obstructed by dirt residues, clean it by means of a blunt flexible stick, taking care not to perforate the rubber tube connected.

Chapter 6

AFTER-SALES SERVICE

If the unit does not function:

- Check that the plug is inserted properly in the power socket;
- Check the level of liquid in the tank.

WARNING

The after-sales service must only be provided by technicians trained by the manufacturer. It is dangerous for anybody else to carry out repairs. If you require assistance, contact your retailer as soon as possible. He will send the unit directly to the manufacturer for the appropriate repairs. Otherwise the manufacturer can be contacted directly at the following address:

TECNO-GAZ S.p.A.
Strada Cavalli, 4
43038 Sala Baganza (Parma)
Tel. +39 0521 833926
Fax +39 0521 833391
info@tecnogaz.com
www.tecnogaz.com

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