

www.tecnogaz.com



KYRI - TYGI

Imbottiture memory foam | Cartella colori

Memory foam upholstery | Color chart



Tecno-Gaz S.p.A.

Strada Cavalli, 4 - 43038 - Sala Baganza - Parma - Italy

Tel./Ph. +39 0521 83.80 Fax +39 0521 83.33.91 - www.tecnogaz.com

Cap. Soc. € 280.000 i.v. C.F. e P.IVA/VAT IT00570950345 - R.E.A. PR 138927 Iscr. Reg. Impr. PR 10061



Tecno-Gaz Spa.



Scheda tecnica | Technical information sheet




















Composizione / Composition	99 % PVC-1 % PU
Peso / Weight	1,140 gr/m ² (tolleranza/tolerance +/-5%)
Spessore / Thickness	1,3 mm (tolleranza/tolerance +/-2)

Avvertenze per il lavaggio: si sconsiglia di lavare con solventi e detersivi a base di alcool.

Cleaning instructions: we recommend to avoid cleaning with alcohol based solvents and detergents.

Cartella colori | Color chart KYRI - TYGI

 70 Nero opaco Matt black	 75 Giallo Yellow	 92 Verde brillante Glossy green	 96 Blu pacifico Pacific blue
 71 Nero brillante Glossy black	 76 Arancio Orange	 93 Verde acqua Water green	 97 Blu elettrico Electric blue
 72 Grigio antracite Anthracite gray	 77 Mirtillo Blueberry	 94 Azzurro brillante Glossy light blue	 98 Blu Blue
 73 Argento Silver	 78 Lilla Lilac	 95 Azzurro Light blue	 99 Bianco White
 74 Grigio chiaro Light gray	 91 Verde chiaro Light green		