



KyriEvo

**The solution for
multidisciplinary surgery room**



KyriEvo

Multidisciplinary integrated surgery room




Quality & Design
MADE IN ITALY

Solutions for managing areas of the dental practice.


Digital Workflow Area


Requirement

To efficiently manage the exchange of information and facilitate the needs of the diagnostic phase and the manufacture of products.

Solutions

Implementation of integrated technologies and flexible operating environments.





Sterilisation Area

■ **Requirement**

Productivity by reducing downtime and allocating staff to higher-value activities.

▶ **Solutions**

Implementation of FAST sterilisation cycles with support of new-generation autoclaves and multifunctional tanks.



Multifunctional Practice Area


■ **Requirement**

Ensure flexibility in the availability of instruments and in management of the clean/dirty cycle.

▶ **Solutions**

Kyri Evo
Kyri Cart
Smart Trolley Evo
Kyri Bridge

Operational flexibility in the dental practice.




Technology.

Kyri Evo is the flexible solution that meets all the needs of the modern dental practice.



Kyri Evo: The heart of multidisciplinary.

A multi-purpose system ready to adapt
to your operating needs.



Kyri Evo. Positions.

Seat




Fowler



Trendelenburg



Surgical table
Height adjustable
from 50 to 81 cm



Integrated system.



Reception

Surgery

Endodontics

Hygiene

Emergency

Implantology

Prosthetics

Kyri Evo: Unique features.

Bars for armrest and drip tray accessories
(removable on request)

Suction arm



Leg safety lock

Operator-programmable memories allowing automatic chair positioning



Slimmer backrest structure: to facilitate positioning of the dentist's legs

The multifunctional surgical chair

The Kyri Evo chair was designed, developed and produced based on the advice and feedback of the greatest dental surgeons and implantologists, facing daily ergonomic pitfalls, hygienic constraints and patient stress. 4 motors that allow customisation of the position for improved patient comfort.

Paediatric headrest
position with child
booster



**Kyri Evo is available
in 20 new colours**
with memory and
standard upholstery



Enveloping backrest for
arm and lumbar support

Seamless Upholstery for
maximum hygiene

Raised leg section

3D headrest
The orbital headrest
provides comfort
for patients during
lateral surgery



Minimum height 50 cm
allowing easy access
even for patients with
reduced mobility

New suction arm
Ambidextrous use.
Possibility of housing up
to 3 instruments. Easily
accessible filter for
cleaning. Wall-mounted
version available.



**Removable
and rotatable
armrest**



Kyri Evo: Your practice's professional image.


The surgical chair that meets the specific needs
of each dental discipline.



Ideal for all

Kyri Evo is the ideal choice for all dental professionals. Versatile, ergonomic and integrated, it exceeds expectations by offering a working experience far superior to that provided by the classic dental unit. Kyri Evo is the best choice for modern dental practices that want to optimise their performance by changing their modus operandi.


Kyri Evo: Operating ergonomics.



Comfort and ergonomics of dental staff

The Kyri Evo dental chair was designed with a particular focus on operating ergonomics and patient comfort. With adjustable functions and a wide range of positions, it ensures optimal posture for the practitioner and a comfortable experience for the patient. The improved ergonomics also help to reduce fatigue for professionals during long treatments, while maintaining efficiency and precision.

The ultimate in ergonomics for high-level performance.



Everything within sight.

Kyri Evo radically changes the concept of workspaces. The dentist can position himself in operating mode without taking his attention away from the patient's mouth, having all the necessary instruments at his disposal with 1st, 2nd and 3rd class movements.

Kyri Evo: Work with four and six hands.



Maximum freedom of movement without obstruction.

The extreme accessibility around the Kyri Evo chair allows you to organise your work with 4 and 6 hands with the aid of the Kyri Cart, Smart Trolley Evo and Kyri bridge devices, becomes particularly easy.

Kyri Evo solution. The multidisciplinary room.

Kyri DSS

Surgical aspirator



Kyri Cart

Surgical cart



KyriEvo



ErgoSeat

Professional orthopaedic seat





ST Evo
SMART TROLLEY

Smart Trolley Evo

Professional furniture line
for dentistry



Orion Line

LED dental lamps




Orion DS



Orion Tech



Orion Light




Kyri Bridge

Bridge table



Kyri Evo solution. Workspace organisation.



Kyri Evo solution.

Flow management.

1

Reception

The multidisciplinary room is configured so as to accommodate the patient, keeping instruments out of view, thereby reducing emotional impact.



2

Use

The dentist and staff can freely avail of the Smart Trolley Evo, the cart and the MasterFlux Smart sedation device as required.



3

Sterilization

When surgery is over the post-operative Smart Trolley Evo is transferred to the sterilisation room for conditioning of the instruments and disposal of material.



4

It returns with clean instruments and replenishment of materials

The intra-operative Smart Trolley can be used to replenish instruments from the warehouse.



The management of workflows within the dental practice is greatly simplified thanks to the wheeled mobility provided on all the complementary devices: Kyri cart, Smart Trolley Evo, Kyri bridge.

Smart Trolley Evo: an image of order and cleanliness.

A dedicated, tidy, room, equipped with professional equipment and devices puts the patient at ease by creating a positive environment.



Smart Trolley Evo: the modular, rational and dynamic system.



intra.electro.post.



Freedom of movement and flow management

Smart Trolley Evo. The evolution of daily operations thanks to the new line of professional furniture for dentistry. Even more efficient. Even smarter..... Flexibility is the key. Operating in a functional manner, with total freedom of movement. Make your surroundings adapt to you.

Surgical lighting

Orion DS

Led lamp with dual function for surgical and dental lighting. It consists of 30 elliptical reflectors divided into 6 modules of 5 LEDs each for a luminous intensity of 140.000 Lux. With 24 cm circular and 20×14 cm elliptical illuminated field.



Light intensity

from 40.000 to 140.000 Lux

Colour temperature

4.500°K and 5.000°K

Illuminated field

Circular 24 cm
Elliptical 20×14 cm

Orion Light

The lightest and most manageable. Dual function for surgical and dental lighting. Light intensity of 140.000 Lux with a circular illuminated field of 24 cm and an elliptical illuminated field of 20×14 cm. This allows the patient's mouth to be moved without having to manipulate the lamp.



Light intensity

from 40.000 to 140.000 Lux

Colour temperature

4.500°K and 5.000°K

Illuminated field

Circular 24 cm
Elliptical 20×14 cm

Orion Tech

Orion Tech is a surgical lamp developed specifically for the dental field, with characteristics calibrated to the needs of the dentist.



Light intensity

from 40.000 to 130.000 Lux

Colour temperature

4.500°K and 5.000°K

Illuminated field

Circular 27 cm
Elliptical 20×14 cm

Surgical aspiration


Kyri DSS

The mobility of Kyri DSS allows the suction arm to be placed closer and better in the work space and facilitates cleaning in the sterilisation room. When not in use, the suction arm is stored under the worktop.




Suction arm

Chair or wall version



Kyri Cart

Kyri Cart offers a high degree of customisation. The 5 instrument compartments allow a multitude of configurations: the basic configuration (syringe, turbine, motor and ultrasound) can be combined with a surgical micromotor.



Kyri Bridge

Mobile instrument cart for implantology and surgery. This space provides a natural connecting and communicating element between the dentist and the assistant. The assistant positions the instruments according to the predetermined sequence and the doctor uses them progressively.

Ergoseat

Professional orthopedic chair with compensated movement system between the seat and backrest, enabling an active position at all times, preventing postural problems.



Tecno Stool

Saddle type stool

Assistant seat



The solution for multidisciplinary surgery room.



MTGZD0358 Rev. 00/11/2023 GraphoServiceParma



Kyri Evo

Scan the QR to download the PDF brochure



Our company for a SUSTAINABLE FUTURE



Scan the QR and download our brochures



Tecno-Gaz S.p.A.

Strada Cavalli, 4 - 43038 - Sala Baganza - Parma - Italy

Tel. +39 0521 83.80 Fax +39 0521 83.33.91 - www.tecnogaz.com

Cap. Soc. € 280.000 i.v. C.F. e P.IVA/VAT IT00570950345 - R.E.A. PR 138927 Iscr. Reg. Impr. PR 10061



www.tecnogaz.com

All rights are reserved. Any variations on images and subjects can be done without notice. Tecno-Gaz S.p.A is not to be considered responsible for damages caused by lack or wrong information here mentioned. Product images are only for advertising purposes.

